



CERTIFICATION by CUSTODIAN of RECORDS
Tucson Police Department Crime Laboratory

TS4

I hereby certify that these documents are a true and correct copy of the periodic maintenance* for the Intoxilyzer 5000 Serial Number 66-001759 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.

A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.080AC solution are 0.070AC to 0.090AC.

Date: 5/19/06 Cal 1 0.080 AC Cal 2 0.081 AC

Date: 6/15/06 Cal 1 0.079 AC Cal 2 0.080 AC

Date: N/A Cal 1 AC Cal 2 AC


Between the dates of 5/19/06 and 6/15/06 the instrument was operating properly and accurately.

The instrument was taken out of service for repair on N/A.
The instrument was returned to service on N/A.

Jennifer L. Coombs
Jennifer L. Coombs

6/20/06

ADPS QAS Permit # QAS 29886
SUBSCRIBED AND SWORN to before me this 20, day of June, 2006.



Kathy A. Kinney
Notary Public

*Includes documents specified in A.A.C. R13-10-104.A.5

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 5000*

1 of 6 files

STANDARD QUALITY ASSURANCE PROCEDURE *AM*

QA SPECIALIST Jennifer Coombs 461160 AGENCY Tucson Police Department
DATE June 15, 2006 TIME 1248
INTOXILYZER SERIAL # 66-001759 LOCATION DUI HQ

() 1. Display reads "READY TO START" or "PUSH BUTTON TO START TEST."

DIAGNOSTIC TESTS

- () 1. DVM Test check. Setting should be between .010 and .600. Mode selection switch S2 on, S1 and S3 off or keyboard menu selection "H". Reading is .369.
- () 2. Display Test check. Mode selection switch S1 on and S2 and S3 off or keyboard menu selection "V".
- () 3. Printer Test check. Mode selection switch S1, S2, S3 off or keyboard menu selection "P".
- () 4. Clock time check. Mode selection switch S10 on or keyboard menu selection "E".
- () 5. Date check. Mode selection switch S11 on or keyboard menu selection "E".

OPERATIONAL TESTS

- () 1. Alcohol-free subject Test result 0. 000 AC
- () 2. Error recognition logic system functioning
Invalid test printed
- () 3. Proper sample recognition system
Invalid sample printed
Deficient sample printed
- () 4. Standard Calibration Check standard 0. 080 AC
Results: 0. 079AC 40.080 AC

Instrument operating accurately and properly. YES NO

COMMENTS: * moved clock 1246 → 1249

SIGNATURE Jennifer L Coombs

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD

Tucson Police Department Crime Laboratory

3 of 6
file

TUCSON POLICE DEPT.
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-001759
06/15/06 13:46

TEST	BRAC VALUE	TIME
AIR BLANK	.000	13:45
CAL. CHECK	.080	13:45
AIR BLANK	.000	13:45

SUBJECT'S NAME

DUI HQ
INSTRUMENT LOCATION

TIME FIRST OBSERVED

46160

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS



INTOXILYZER™ INSTRUMENT PRINTER CARD

© 1986 by CMI INC

TUCSON POLICE DEPT.
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-001759
06/15/06 13:44

TEST	BRAC VALUE	TIME
AIR BLANK	.000	13:43
CAL. CHECK	.079	13:44
AIR BLANK	.000	13:44

SUBJECT'S NAME

DUI HQ
INSTRUMENT LOCATION

TIME FIRST OBSERVED

46160

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS



INTOXILYZER™ INSTRUMENT PRINTER CARD

© 1986 by CMI INC

4 of 6
by file

SN 66-001759
E767.14
INVALID TEST

06/15/06
13:08

SAMPLE INTRODUCED
AT IMPROPER TIME.

SUBJECT'S NAME	
TIME FIRST OBSERVED	INSTRUMENT LOCATION DUI HQ
OPERATOR 461160	
ADDITIONAL INFORMATION AND/OR REMARKS	

INTOXILYZER[®] INSTRUMENT PRINTER CARD



© 1986 by CMI INC

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R13-10-104(A)

5 of 6
file

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 5000*

STANDARD CALIBRATION CHECK PROCEDURE

64

QA SPECIALIST Jennifer Coombs 46160 AGENCY Tucson Police Department
DATE June 15, 2006 TIME 1248
INTOXILYZER SERIAL # 66-001759 LOCATION D111 HQ

- (✓) 1. Pour a standard alcohol concentration solution into a clean dry simulator jar and assemble the simulator. Ensure that a tight seal has been made. Standard value: 0. 080 AC.
- (✓) 2. Turn on the simulator and allow the temperature to reach 34°C ± 0.2°C.
- (✓) 3. Set Intoxilyzer mode selection in the ACA mode by switching mode selection switch #9 on or selecting "C" on keyboard menu.
- (✓) 4. Attach simulator to the simulator entrance port on the Intoxilyzer.
- ✓ ✓ 5. Intoxilyzer 5000 display reads "READY TO START" or "PUSH BUTTON".
- ✓ ✓ 6. Push Start Test button or press enter on keyboard.
- ✓ ✓ 7. Insert card in response to display.
- ✓ ✓ 8. Air blank completed.
- ✓ ✓ 9. Standard Calibration Check completed. Test results 0. 01990 #0.080 AC.
- ✓ ✓ 10. Air blank completed.
- ✓ ✓ 11. When test complete, remove printed record. Attach the record to the completed checklist.
- (✓) 12. Return mode selection switch #9 to off position after all calibration checks are complete or type Q and enter on keyboard.

SIGNATURE Jennifer L. Coombs

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD

DPS Form Exh E-2 (Rev. 05-1)

TUCSON POLICE DEPT.
 INTOXILYZER - ALCOHOL ANALYZER
 AZ MODEL 5000 SN 66-001759
 06/15/06 13:13

TEST	BRAC VALUE	TIME
AIR BLANK	.000	13:10
*SUBJECT TEST	.042	13:13
AIR BLANK	.000	13:13

* DEFICIENT SAMPLE - VALUE
 PRINTED HAS HIGHEST OBTAINED.

SUBJECT'S NAME _____
 TIME FIRST OBSERVED _____
 INSTRUMENT LOCATION DUI HQ

OPERATOR Holloo

ADDITIONAL INFORMATION AND/OR REMARKS _____

EMI INC.

INTOXILYZER INSTRUMENT PRINTER CARD
 © 1986 by CMI INC.

by
bot to file

TUCSON POLICE DEPT.
 INTOXILYZER - ALCOHOL ANALYZER
 AZ MODEL 5000 SN 66-001759
 06/15/06 13:09

TEST	BRAC VALUE	TIME
AIR BLANK	.000	13:09
INVALID SAMPLE	.MDA	13:09
AIR BLANK	.000	13:09

BEGIN NEW DEPRIVATION PERIOD

SUBJECT'S NAME _____
 TIME FIRST OBSERVED _____
 INSTRUMENT LOCATION DUI HQ

OPERATOR Holloo

ADDITIONAL INFORMATION AND/OR REMARKS Sweet Breath

EMI INC.

INTOXILYZER INSTRUMENT PRINTER CARD
 © 1986 by CMI INC.