



CERTIFICATION by CUSTODIAN of RECORDS  
 Tucson Police Department Crime Laboratory  
 270 S. Stone Avenue Tucson, AZ 85701-1917

I hereby certify that these documents are a true and correct copy of the periodic maintenance\* for the Intoxilyzer 8000 Serial Number 80-002201 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.

TS4  
1-22-2010

A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.

Date: 1/19/10 Cal 1 0.100 AC Cal 2 0.098 AC

Between the dates of - and - the instrument was operating properly and accurately.


The instrument was taken out of service for repair on -.  
 The instrument was returned to service on 1/19/10.

Sarah Cione  
 Sarah L. Cione

1/20/10  
 ADPS QAS Permit # 33355

SUBSCRIBED AND SWORN to before me this 20<sup>th</sup> day of January, 2010.

A.M. Hernandez  
 Notary Public



\*Includes documents specified in A.A.C. R13-10-104.A.5

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY  
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

TS4  
1.22.2010

QA SPECIALIST Sarah Cione AGENCY Tucson Police Department

DATE 1/19/10 TIME 1435

INTOXILYZER SERIAL # 80-002201 LOCATION Crime Lab

(  ) 1. Display reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- (  ) 1. Clock time check.
- (  ) 2. Date check.

OPERATIONAL TESTS

- (  ) 1. Alcohol-free subject test result 0. 000 AC.
- (  ) 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed
- (  ) 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.
- (  ) 4. Standard Calibration Check standard 0. 100 AC  
Result: 0. 098 AC

Instrument is operating properly and accurately. YES  NO

COMMENTS  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SIGNATURE Sarah Cione

2 of 12 sc

INTOXILYZER 8000  
Location: ARIZONA  
Serial Number : 80-002201  
01/19/2010

14:37:09

BY

Standard Lot#: 826001I  
Last Changed By: S. CIONE #50616

QAS: S. CIONE #50616  
TUCSON PD

Operator: S. CIONE #50616  
TUCSON PD

Subject: SQAP, SQAP  
DOB: 11/11/1111  
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	14:37:57
Diagnostic Test	Pass	14:38:24
Air Blank	0.000	14:38:52
0.100 Cal Check	0.099	14:39:11
Air Blank	0.000	14:39:40
Subject Test	0.000	14:39:58
Air Blank	0.000	14:40:25
Five Minute Wait	ABT*	
Air Blank	0.000	14:41:02

\*Sequence Aborted

**Not a Successfully Completed Test Sequence**

3 of 12 sc

INTOXILYZER 8000  
Location: ARIZONA  
Serial Number : 80-002201  
01/19/2010 14:41:28

154

Standard Lot#: 826001I  
Last Changed By: S. CIONE #50616

QAS: S. CIONE #50616  
TUCSON PD

Operator: S. CIONE #50616  
TUCSON PD

Subject: SQAP, SQAP  
DOB: 11/11/1111  
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	14:42:38
Air Blank	IPS**	14:42:47

\*RFI Detect  
\*\*Improper Sample

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: ARIZONA  
Serial Number : 80-002201  
01/19/2010

14:42:52

4 of 12 sc

AS4

Standard Lot#: 826001I  
Last Changed By: S. CIONE #50616

QAS: S. CIONE #50616  
TUCSON PD

Operator: S. CIONE #50616  
TUCSON PD

Subject: SQAP, SQAP  
DOB: 11/11/1111  
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	14:43:37
Diagnostic Test	Pass	14:44:04
Air Blank	0.000	14:44:32
0.100 Cal Check	0.100	14:44:51
Air Blank	0.000	14:45:20
Subject Test	0.290*	14:48:28
Air Blank	0.000	14:49:01
Air Blank	0.000	14:49:28
Subject Test	XXX**	14:49:47
Air Blank	0.000	14:50:20

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: ARIZONA  
Serial Number: 80-002201  
01/19/2010 14:51:06

5 of 12 sc

-----  
Standard Lot#: 826001I  
Last Changed By: S. CIONE #50616  
-----

154

QAS: S. CIONE #50616  
TUCSON PD  
-----

Test	g/210L	Time
Air Blank	0.000	14:51:33
Cal Check	0.098	14:51:52
Air Blank	0.000	14:52:21

*Sarah Cione*

Operator's Signature

EXHIBIT G-2

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY  
A.A.C. R13-10-104(A)

6 of 12  
sc

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

TSY  
1-22-2010

QA Specialist Sarah Cione AGENCY Tucson Police Department  
DATE 1/19/10 TIME 1435  
Intoxilyzer Serial # 80-002201 Location Crime Lab

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0. 100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 100 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Sarah Cione

DPS Form Exh G-2 (Rev. 05-01)

INTOXILYZER 8000  
Location: ARIZONA  
Serial Number: 80-002201  
01/19/2010 14:35:46

7 of 12 sc

154

-----  
Standard Lot#: 826001I  
Last Changed By: S. CIONE #50616  
-----

QAS: S. CIONE #50616  
TUCSON PD  
-----

Test	g/210L	Time
Air Blank	0.000	14:36:14
Cal Check	0.100	14:36:33
Air Blank	0.000	14:37:02

Sarah Cione  
Operator's Signature



TUCSON POLICE DEPARTMENT  
CRIME LABORATORY

8 of 12<sub>sc</sub>

INTOXILYZER SERVICE RECORD

By  
1-22-2010

QAS: Sarah Gione

Intoxilyzer Serial Number: 80-002201

Date Returned to Service: 1/19/2010



TUCSON POLICE DEPARTMENT  
CRIME LABORATORY

9 of 12 SC

INTOXILYZER SERVICE RECORD  
FOR  
CALIBRATION GAS TANK CHANGE

Date: 1/19/2010  
 QAS: Sarah Cione  
 Intoxilyzer Serial Number: 80-002201  
 New Lot Number: 826001 I  
 New Expiration Date: 9/24/2010

TSC  
1-22-2010

- (✓) 1. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 2. Remove the old gas tank and install the new one.
- (✓) 3. Hit the "ESC" key twice and swipe the QAS card.
- (✓) 4. Go to the "M-Maintenance" menu and select "S-Standard Change." Select "D" for dry gas and enter the following information:  
 Enter Target Value: 0.100  
 Enter Control Lot #: From the tank label  
 Enter Exp Date: From the tank label (mm/dd/yyyy)  
 Edit info Y/N? Y or N  
 Update Std Info? Y ENTER  
 "Please wait saving..."
- (✓) 5. Hit the "ESC" key twice.
- (✓) 6. Go to the "M-Maintenance" menu and select "Control Testing." Select "D" for dry control test. After selection is made press ENTER.
- (✓) 7. Air blank completed.
- (✓) 8. Calibration check completed. Test results 0.100 AC.
- (✓) 9. Air blank completed.
- (✓) 10. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Sarah Cione

10 of 12 sc

INTOXILYZER 8000  
Location: ARIZONA  
Serial Number: 80-002201  
01/19/2010 14:33:58

-----  
Standard Lot#: 826001I  
Last Changed By: S. CIONE #50616  
-----

TS4

-----  
QAS: S. CIONE #50616  
TUCSON PD  
-----

Test	g/210L	Time
Air Blank	0.000	14:34:26
Cal Check	0.100	14:34:44
Air Blank	0.000	14:35:14

*Sarah Cione*

Operator's Signature



12 of 12 sc  
154

## Certificate of Calibration

### **Intoxilyzer**® 8000

*This is to certify that the calibration of **Intoxilyzer**® serial number 80-002201, manufactured by CMI Inc., a subsidiary of MPD, Inc. of Owensboro, Kentucky, was tested and found to be in compliance with the National Highway Traffic Safety Administration Standard for Devices to Measure Breath Alcohol (F.R., No. 179 48705-48710 Sept. 17, 1993) and with International Recommendation OIML R126 Edition 1998 (E). Calibration solutions are traceable to NIST standard material 1828.*

Date 1/5/10 Signed Kristi Robinson  
Technician



316 East 9<sup>th</sup> Street  
Owensboro, KY 42303  
USA

Part No. 650519 CMI 7/9/09

Pursuant to Section 28.1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.