

1076
40

EXHIBIT Q

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 5000*

AGENCY Tucson Police QA SPECIALIST Irene Doran (46331)
INTOXILYZER SERIAL #04-001226 (#8) LOCATION DUI HQ
(Print Name)
DATE February 4 192004 TIME 1632

✓ Display reads "READY TO START" or "PUSH BUTTON TO START TEST."

DIAGNOSTIC TESTS

- ✓ DVM Test check. Setting should be between .200 and .600. Mode selection switch S2 on, S1 and S3 off or keyboard menu selection "H". Reading is .543.
- ✓ Display Test check. Mode selection switch S1 on and S2 and S3 off or keyboard menu selection "V".
- ✓ Printer Test check. Mode selection switch S1, S2, S3 off or keyboard menu selection "P".
- ✓ Clock time check. Mode selection switch S10 on or keyboard menu selection "E".
- ✓ Date check. Mode selection switch S11 on or keyboard menu selection "E".

OPERATIONAL TESTS

- ✓ Alcohol free subject Test result 0.000 AC
- ✓ Error recognition logic system functioning
✓ Invalid test printed
- ✓ Proper sample recognition system
✓ Invalid sample printed
✓ Deficient sample printed
- N/A Completeness of sample purge with collection tube**
- ✓ Calibration standard 0.080 AC
- Results: 0.083 AC + .0861C ✓
- Instrument operating properly and accurately. YES ✓ NO _

COMMENTS:

SIGNATURE

Irene Doran

0. N/A alcohol concentration standard collected for subsequent analysis.**

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD

**THIS STEP IS ONLY REQUIRED IF THIS DEVICE IS BEING USED FOR SAMPLE CAPTURE.

SN 64-061236

647.02

62/04/04

16:42

#8

EEEEE

EEEEEE

EEEEE

FEDERAL SIGNAL CORP. / CMI INC
INTOXILYZER - ALCOHOL ANALYZER
HZ MODEL 5000 SN 64-001226
02/04/84

#8

TEST	BAC	VALUE	TIME
AIR BLANK	.000	16:32	
SUBJECT TEST	.000	16:32	
AIR BLANK	.000	16:32	

DW = .543

TEST	BAC	VALUE	TIME
AIR BLANK	.000	16:32	
SUBJECT TEST	.000	16:32	
AIR BLANK	.000	16:32	

TEST	BAC	VALUE	TIME
DUT HQ	.000	16:32	
SUBJECT NAME	Death	46331	
INSTRUMENT LOCATION	DUT HQ		
OPERATOR			

TEST	BAC	VALUE	TIME
DUT HQ	.000	16:32	
SUBJECT NAME	Death	46331	
INSTRUMENT LOCATION	DUT HQ		
OPERATOR			

TEST	BAC	VALUE	TIME
DUT HQ	.000	16:32	
SUBJECT NAME	Death	46331	
INSTRUMENT LOCATION	DUT HQ		
OPERATOR			

INTOXILYZER INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

296

296

© 1986 by CMI INC.

INTOXILYZER INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

296

INTOXILYZER INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

296

Tucson Police Department Crime Laboratory

386
84

FEDERAL SIGNAL CORP. / CMI INC
INTOXILYZER - ALCOHOL ANALYZER #8
AZ MODEL 5000 SN 64-001226
02/04/84

TEST	BAC	VALUE	TIME
AIR BLANK	.000	16:34	
CAL. CHECK	.083	16:34	
AIR, BLANK	.000	16:35	

TEST	BAC	VALUE	TIME
AIR BLANK			AIR BLANK
CAL. CHECK			CAL. CHECK
AIR, BLANK			AIR, BLANK

SUBJECT'S NAME	DUT HQ
INSTRUMENT LOCATION	INSTRUMENT LOCATION
TIME FIRST OBSERVED	TIME FIRST OBSERVED
OPERATOR	OPERATOR
ADDITIONAL INFORMATION AND/OR REMARKS	

SUBJECT'S NAME	DUT HQ
INSTRUMENT LOCATION	INSTRUMENT LOCATION
TIME FIRST OBSERVED	TIME FIRST OBSERVED
OPERATOR	OPERATOR
ADDITIONAL INFORMATION AND/OR REMARKS	

CMI Inc.

INTOXILYZER INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

INTOXILYZER INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

105

40g 6
JW

SL 64-031326
247.02
INVALID TEST

02/04/84
16:45

#8

SUBJECT'S NAME

DUI HQ

TIME FIRST OBSERVED

INSTRUMENT LOCATION

OPERATOR

DORAN (40331)

ADDITIONAL INFORMATION AND/OR REMARKS

pushed button during
acc blank

INTOXILYZER[®] INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

CMI[®]
INC.

JW

EXHIBIT P

586
95

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 5000*

STANDARD CALIBRATION CHECK PROCEDURE

Agency Tucson Police Department Date 2/4/04 Time 1632
Intoxilyzer Serial# 64-00126(8) Location DUI HQ
QA Specialist Irene Doran (46331)
(Print Name)

- (✓) 1. Pour a standard alcohol solution of known value into a clean dry simulator jar and assemble the simulator. Insure that a tight seal has been made.
Standard value: 0.080 AC
- (✓) 2. Turn on the simulator and allow the temperature to reach $34^{\circ}\text{C} \pm .2^{\circ}\text{C}$.
- (✓) 3. Set Intoxilyzer mode selection in the ACA mode by switching mode selection switch #9 on or selecting "C" on keyboard menu.
- (✓) 4. Attach simulator to the simulator entrance port on the Intoxilyzer.
- (✓) 5. Intoxilyzer 5000 display reads "READY TO START" or "PUSH BUTTON".
- (✓) 6. Push Start Test button or press enter on keyboard.
- (✓) 7. Insert card in response to display.
- (✓) 8. Air blank completed.
- (✓) 9. Calibration check completed.
Test results 0.083 AC + 0.086AC
- (✓) 10. Air blank completed.
- (✓) 11. When display reads Test Complete remove evidence card. Attach the card to the completed checklist.
- (✓) 12. Return mode selection switch #9 to off position after all calibration checks are complete or type Q and enter on keyboard.

SIGNATURE

Irene Doran

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD

DHS/DSLS/Form C116 (Rev. 7-93)

NOTE: Exhibit P last performed on
12/26/03. Due to scheduling conflicts,
the Exhibit P was run 9 days past its
due date of within 31 days.

FLORIDA SIGNAL CORP., / CMI INC
INTOXILYZER - ALCOHOL ANALYZER
A2 MODEL 5000 SN 64-001226
02/04/84

Feb 04, 1984 SIGNAL CORP., / CMI INC
INTOXILYZER - ALCOHOL ANALYZER #8
A2 MODEL 5000 SN 64-001226
02/04/84

TEST	BAC	VALUE	TIME
AIR BLANK	.000	16:44	
INHALID SAMPLE	.XXX	16:44	
AIR BLANK	.000	16:45	

* DEFICIENT SAMPLE - VALUE
PRINTED WAS HIGHEST OBTAINED.

TEST	BAC	VALUE	TIME
AIR BLANK	.000	16:44	
*SUBJECT TEST			16:48
AIR BLANK	.000	16:49	

* DEFICIENT SAMPLE - VALUE
PRINTED WAS HIGHEST OBTAINED.

DoeAN

INSTRUMENT LOCATION

TIME FIRST OBSERVED

SUBJECT'S NAME

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

INTOXILYZER INSTRUMENT PRINTER CARD

DoeAN (46331)

INSTRUMENT LOCATION

TIME FIRST OBSERVED

SUBJECT'S NAME

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

INTOXILYZER INSTRUMENT PRINTER CARD

696
#8

CMI

INTOXILYZER INSTRUMENT PRINTER CARD

© 1986 by CMI INC

© 1986 by CMI INC

© 1986 by CMI INC