

EXHIBIT Q

1 of 6
JLC

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES *fh*
INTOXILYZER MODEL 5000*

AGENCY Tucson Police QA SPECIALIST Jennifer Coombs 46160
 (Print Name)
 INTOXILYZER SERIAL # 64-001226 (#8) LOCATION LAB
 DATE July 22 19 2004 TIME 1240

Display reads "READY TO START" or "PUSH BUTTON TO START TEST."

DIAGNOSTIC TESTS

- DVM Test check. Setting should be between .200 and .600. Mode selection switch S2 on, S1 and S3 off or keyboard menu selection "H". Reading is .406.
- Display Test check. Mode selection switch S1 on and S2 and S3 off or keyboard menu selection "V".
- Printer Test check. Mode selection switch S1, S2, S3 off or keyboard menu selection "P".
- Clock time check. Mode selection switch S10 on or keyboard menu selection "E".
- Date check. Mode selection switch S11 on or keyboard menu selection "E".

OPERATIONAL TESTS

- Alcohol free subject Test result 0.000 AC
- Error recognition logic system functioning
 - Invalid test printed
- Proper sample recognition system
 - Invalid sample printed
 - Deficient sample printed
- N/A Completeness of sample purge with collection tube**
 - Calibration standard 0.080 AC
 - Results: 0.075 AC & 0.073 AC
- Instrument operating properly and accurately. YES NO

COMMENTS: _____

SIGNATURE

Jennifer L. Coombs

Q. N/A alcohol concentration standard collected for subsequent analysis.**

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD

**THIS STEP IS ONLY REQUIRED IF THIS DEVICE IS BEING USED FOR SAMPLE CAPTURE.

206
file

#8

SH 64-001226	07/22/04
E47, Q2	12:47
EEEEEEEEE	EEEEEEEEE
HHHHHHHHH	HHHHHHHHH
Q 307	-----
Q 306	-----
Q 305	-----
Q 302	-----
Q 303	-----
Q 304	-----
Q 301	-----

FEDERAL SIGNAL CDRP. / CMI INC
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 64-001226
07/22/04

TEST

	BAC VALUE	TIME
AIR BLANK	.000	12:50
SUBJECT TEST	.000	12:51
AIR BLANK	.000	12:51

$$\text{D/H} = .406$$

SUBJECT'S NAME _____

LAB

Hello

OPERATOR

INSTRUMENT LOCATION

TIME FIRST OBSERVED

ADDITIONAL INFORMATION AND/OR REMARKS

INTOXILYZER™ INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

CMI

CMI

4 of 6
file

SN 64-001226
E47.02
INVALID TEST

07/22/04
12:52

#8

SUBJECT'S NAME

LAB

TIME FIRST OBSERVED

INSTRUMENT LOCATION

46160

OPERATOR

ADDITIONAL INFORMATION AND / OR REMARKS

pushed start + test button at
wrong time

INTOXILYZER[®] INSTRUMENT PRINTER CARD

CMI

© 1986 by CMI INC.

5 of
6
file

EXHIBIT P

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 5000*

STANDARD CALIBRATION CHECK PROCEDURE

Agency Tucson Police Department Date 7/22/04 Time 1240
Intoxilyzer Serial# 64-001226 (#8) Location LAB
QA Specialist Jennifer Combs 460160
(Print Name)

- (✓) 1. Pour a standard alcohol solution of known value into a clean dry simulator jar and assemble the simulator.
Insure that a tight seal has been made.
Standard value: 0.080 AC
- (✓) 2. Turn on the simulator and allow the temperature to reach 34°C±.2°C.
- (✓) 3. Set Intoxilyzer mode selection in the ACA mode by switching mode selection switch #9 on or selecting "C" on keyboard menu.
- (✓) 4. Attach simulator to the simulator entrance port on the Intoxilyzer.
- (✓) 5. Intoxilyzer 5000 display reads "READY TO START" or "PUSH BUTTON".
- (✓) 6. Push Start Test button or press enter on keyboard.
- (✓) 7. Insert card in response to display.
- (✓) 8. Air blank completed.
- (✓) 9. Calibration check completed.
Test results 0.075AC & 0.073 AC
- (✓) 10. Air blank completed.
- (✓) 11. When display reads Test Complete remove evidence card.
Attach the card to the completed checklist.
- (✓) 12. Return mode selection switch #9 to off position after all calibration checks are complete or type Q and enter on keyboard.

SIGNATURE

Jennifer L. Combs

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD

6 of
6

FEDERAL SIGNAL CORP. / CMI INC #8
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 64-001226
07/22/84

TEST	BAC	VALUE	TIME
AIR BLANK	.000	12:53	
INVALID SAMPLE	.XXX	12:54	
AIR BLANK	.000	12:54	

* DEFICIENT SAMPLE - VALUE
PRINTED WAS HIGHEST OBTAINED.

FEDERAL SIGNAL CORP. / CMI INC #8
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 64-001226
07/22/84

TEST	BAC	VALUE	TIME
AIR BLANK	.000	12:53	
INVALID SAMPLE	.XXX	12:54	
AIR BLANK	.000	12:54	

* DEFICIENT SAMPLE - VALUE
PRINTED WAS HIGHEST OBTAINED.

SUBJECT'S NAME _____
TIME FIRST OBSERVED _____
INSTRUMENT LOCATION _____
OPERATOR _____
ADDITIONAL INFORMATION AND/OR REMARKS
46160
Sweet Breath

SUBJECT'S NAME _____
TIME FIRST OBSERVED _____
INSTRUMENT LOCATION _____
OPERATOR _____
ADDITIONAL INFORMATION AND/OR REMARKS
46160
Sweet Breath

INTOXILYZER™ INSTRUMENT PRINTER CARD
© 1986 by CMI INC.

CMI Inc

CMI Inc

INTOXILYZER™ INSTRUMENT PRINTER CARD
© 1986 by CMI INC.