

EXHIBIT Q

1 of 6
file

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. B9-14-494 (A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

**STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 5000***

AGENCY Tucson Police QA SPECIALIST Jennifer Coombs 766160
INTOXILYZER SERIAL # 64-001226 (#8) LOCATION Dil Hq
(Print Name)
DATE December 29 TIME 0844

✓ Display reads "READY TO START" or "PUSH BUTTON TO START TEST."

DIAGNOSTIC TESTS

- DVM Test check. Setting should be between .200 and .600. Mode selection switch S2 on, S1 and S3 off or keyboard menu selection "H". Reading is .411.

Display Test check. Mode selection switch S1 on and S2 and S3 off or keyboard menu selection "V".

Printer Test check. Mode selection switch S1, S2, S3 off or keyboard menu selection "P".

Clock time check. Mode selection switch S10 on or keyboard menu selection "E".

Date check. Mode selection switch S11 on or keyboard menu selection "E".

OPERATIONAL TESTS

- ✓ Alcohol free subject Test result 0.000 AC
✓ Error recognition logic system functioning
✓ Invalid test printed
✓ Proper sample recognition system
✓ Invalid sample printed✓
✓ Deficient sample printed✓
N/A Completeness of sample purge with collection tube**
✓ Calibration standard 0.080 AC
Results: 0.084AC & 0.084AC
Instrument operating properly and accurately. YES✓ NO__

COMMENTS:

SIGNATURE John F. Coddington 46160

O. N/A alcohol concentration standard collected for subsequent analysis.**

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD
**THIS STEP IS ONLY REQUIRED IF THIS DEVICE IS BEING USED FOR
SAMPLE CAPTURE.

SH 64-661226

卷之三

12/29/64 #8

FED
INTC
42 12/21

INTOXICATOR STAINLESS STEEL CO., LTD. / CHI LIN
INTOXICATOR - ALCOHOL ANALYZER
HQ2 MODEL 5000 SN 64-061226

FEDERAL SIGNAL CORP. / CHI INC
INTOXILYZER - ALCOHOL ANALYZER
#2 MODEL 5000 SN 64-061226
12/29/04

TEST	BAC	VALUE	TIME
AIR BLANK	.000		08:50
SUBJECT TEST	.000		08:50
AIR BLANK	.000		08:51

$$\Delta V/H = 4/11$$

DW H0
INSTRUMENT LOCATION

SUBJECT'S NAME	4611
TIME FIRST OBSERVED	OPERATOR
ADDITIONAL INFORMATION AND COMMENTS	

SUBJECT'S NAME	Dell HQ	INSTRUMENT LOCATION
TIME FIRST OBSERVED	46160	OPERATOR
ADDITIONAL INFORMATION AND/OR REMARKS		

INTOXILYZER® INSTRUMENT PRINTER CARD

© 1996 by **Cengage**

105

INTOXILYZER® INSTRUMENT PRINTER CARD

206
file

卷之三

3 of 6
file

Tucson Police Department Crime Laboratory

#8

FEDERAL SIGNAL CORP. / CMI INC
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 64-001226
12/29/84

TEST	BAC	VALUE	TIME
AIR BLANK	.000	09:05	
CAL. CHECK	.084	09:05	
AIR BLANK	.000	09:06	

TEST	BAC	VALUE	TIME
AIR BLANK	.000	09:05	
CAL. CHECK	.084	09:05	
AIR BLANK	.000	09:06	

TEST	BAC	VALUE	TIME
AIR BLANK	.000	09:07	
CAL. CHECK	.084	09:07	
AIR BLANK	.000	09:08	

DUT HQ

SUBJECT'S NAME	INSTRUMENT LOCATION

TIME FIRST OBSERVED	OPERATOR	ADDITIONAL INFORMATION AND/OR REMARKS

DUT HQ

SUBJECT'S NAME	INSTRUMENT LOCATION

TIME FIRST OBSERVED	OPERATOR	ADDITIONAL INFORMATION AND/OR REMARKS

DUT HQ

SUBJECT'S NAME	INSTRUMENT LOCATION

TIME FIRST OBSERVED	OPERATOR	ADDITIONAL INFORMATION AND/OR REMARKS

DUT HQ

SUBJECT'S NAME	INSTRUMENT LOCATION

TIME FIRST OBSERVED	OPERATOR	ADDITIONAL INFORMATION AND/OR REMARKS

CMI

INTOXILYZER INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

CMI

INTOXILYZER INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

4616
file

SN 64-001226
E+7.02
INVALID TEST
SAMPLE INTRODUCED
AT IMPROPER TIME.

12/29/84 #8
08:52

SUBJECT'S NAME

DUI HQ

TIME FIRST OBSERVED

INSTRUMENT LOCATION

46160

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

INTOXILYZER® INSTRUMENT PRINTER CARD

CMI

© 1986 by CMI INC.

EXHIBIT P

5 of 6
file

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 5000*STANDARD CALIBRATION CHECK PROCEDURE *AFM*

Agency Tucson Police Department Date 12/29/04 Time 0844
Intoxilyzer Serial# 164-001226(#8) Location DUI HQ
QA Specialist Jennifer Combs 460160

(Print Name)

- (✓) 1. Pour a standard alcohol solution of known value into a clean dry simulator jar and assemble the simulator. Insure that a tight seal has been made.
Standard value: 0.080 AC
- (✓) 2. Turn on the simulator and allow the temperature to reach $34^{\circ}\text{C} \pm .2^{\circ}\text{C}$.
- (✓) 3. Set Intoxilyzer mode selection in the ACA mode by switching mode selection switch #9 on or selecting "C" on keyboard menu.
- (✓) 4. Attach simulator to the simulator entrance port on the Intoxilyzer.
- ✓(✓) 5. Intoxilyzer 5000 display reads "READY TO START" or "PUSH BUTTON".
- ✓(✓) 6. Push Start Test button or press enter on keyboard.
- ✓(✓) 7. Insert card in response to display.
- ✓(✓) 8. Air blank completed.
- ✓(✓) 9. Calibration check completed.
Test results 0.084 AC & 0.084 AC
- ✓(✓) 10. Air blank completed.
- ✓(✓) 11. When display reads Test Complete remove evidence card. Attach the card to the completed checklist.
- (✓) 12. Return mode selection switch #9 to off position after all calibration checks are complete or type Q and enter on keyboard.

SIGNATURE

Jennifer L. Combs

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD

FEDERAL SIGNAL CORP., CMII INC
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 64-001226
12/29/84

FEDERAL SIGNAL CORP., CMII INC
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 64-001226
12/29/84

TEST	BAC	VALUE	TIME
AIR BLANK	.000	08:57	
INVALID SAMPLE	XXX	08:58	
AIR BLANK	.000	08:58	

* DEFICIENT SAMPLE - VALUE
PRINTED WAS HIGHEST OBTAINED.

TEST

TEST	BAC	VALUE	TIME
AIR BLANK	.000	08:59	
SUBJECT TEST	.084	09:01	
AIR BLANK	.000	09:01	

* DEFICIENT SAMPLE - VALUE
PRINTED WAS HIGHEST OBTAINED.

TEST	BAC	VALUE	TIME
AIR BLANK	.000	08:59	
SUBJECT TEST	.084	09:01	
AIR BLANK	.000	09:01	

DUI HQ

SUBJECT'S NAME _____
INSTRUMENT LOCATION _____

TIME FIRST OBSERVED _____

OPERATOR _____

ADDITIONAL INFORMATION AND/OR REMARKS _____

DUI HQ

SUBJECT'S NAME _____
INSTRUMENT LOCATION _____

TIME FIRST OBSERVED _____

OPERATOR _____

ADDITIONAL INFORMATION AND/OR REMARKS _____

Hollo

Sweet Breath

6 of 6 file

INTOXILYZER INSTRUMENT PRINTER CARD
© 1986 by CMII INC

INTOXILYZER INSTRUMENT PRINTER CARD
© 1986 by CMII INC

CMII INC

CMII INC