

EXHIBIT Q

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES *SM*
INTOXILYZER MODEL 5000*

AGENCY Tucson Police QA SPECIALIST Jennifer Coombs 46160
(Print Name)
INTOXILYZER SERIAL #64-002505 (#11) LOCATION DUI HQ
DATE October 7 19 2004 TIME 1233
jlc

☒ Display reads "READY TO START" or "PUSH BUTTON TO START TEST."

DIAGNOSTIC TESTS

☐ DVM Test check. Setting should be between .200 and .600. Mode selection switch S2 on, S1 and S3 off or keyboard menu selection "H". Reading is _____.

☐ Display Test check. Mode selection switch S1 on and S2 and S3 off or keyboard menu selection "V".

☒ Printer Test check. Mode selection switch S1, S2, S3 off or keyboard menu selection "P".

☒ Clock time check. Mode selection switch S10 on or keyboard menu selection "E".

☒ Date check. Mode selection switch S11 on or keyboard menu selection "E".

OPERATIONAL TESTS

☒ Alcohol free subject Test result 0. AC

☐ Error recognition logic system functioning
Invalid test printed

☐ Proper sample recognition system
Invalid sample printed

☐ Deficient sample printed

☐ N/A Completeness of sample purge with collection tube**

☒ Calibration standard 0. AC

Results: 0.

Instrument operating properly and accurately. YES ☒ NO

COMMENTS: *Printer pulling cards into instrument.
** Instrument reading Ambient failed during all first air blanks.
Pigled for repair

SIGNATURE Jennifer L. Coombs

Q. N/A alcohol concentration standard collected for subsequent analysis.**

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD
**THIS STEP IS ONLY REQUIRED IF THIS DEVICE IS BEING USED FOR
SAMPLE CAPTURE.

1 of 2
jlc

EXHIBIT P

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 5000*

STANDARD CALIBRATION CHECK PROCEDURE ¹³⁴

Agency Tucson Police Department Date 10/7/04 Time 1233
Intoxilyzer Serial# 04-002505 (#11) Location DUI HQ
QA Specialist Jennifer Combs

(Print Name)

- (✓) 1. Pour a standard alcohol solution of known value into a clean dry simulator jar and assemble the simulator. Insure that a tight seal has been made. Standard value: 0.080 AC
- (✓) 2. Turn on the simulator and allow the temperature to reach $34^{\circ}\text{C} \pm .2^{\circ}\text{C}$.
- (✓) 3. Set Intoxilyzer mode selection in the ACA mode by switching mode selection switch #9 on or selecting "C" on keyboard menu.
- (✓) 4. Attach simulator to the simulator entrance port on the Intoxilyzer.
- (✓) 5. Intoxilyzer 5000 display reads "READY TO START" or "PUSH BUTTON".
- (✓) 6. Push Start Test button or press enter on keyboard.
- (N/A) 7. Insert card in response to display.
- (*) 8. Air blank completed.
- () 9. Calibration check completed. Test results 0. AC
- () 10. Air blank completed.
- () 11. When display reads Test Complete remove evidence card. Attach the card to the completed checklist.
- () 12. Return mode selection switch #9 to off position after all calibration checks are complete or type Q and enter on keyboard.

SIGNATURE

Jennifer L. Combs

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD

DHS/DSLS/Form C116 (Rev.7-93)

* Ambient Failed

1 of 1
JHC



TUCSON POLICE DEPARTMENT
CRIME LABORATORY

INTOXILYZER SERVICE RECORD

Date: 10/7/04
Intoxilyzer Serial Number: 64-002505 (#11)
QAS: Jennifer Coombs 46160

B4

Problem Description:

Printer pulling cards all the way into instrument. Also, during first air blank of any test, reading "Ambient Failed"

If the instrument is taken out of service, complete the following:

Date out of service 10/7/04

Calibration Check completed prior to removal from service:

☐ Yes Date: _____ ☒ No

SQAP completed prior to removal from service:

☐ Yes Date: _____ ☒ No

Date returned to service: 11/10/04
25 #46331

Corrective Action Taken or Repairs: Sent to Applied Electronics.



MEMORANDUM

DATE: October 19, 2004

TO: Whom it may concern

FROM: Jennifer Coombs
Intoxilyzer 5000 QAS
Tucson Police Department
Crime Lab
270 S. Stone Avenue
Tucson, AZ 85701-1917
(520) 791-4494

SUBJECT: Intoxilyzer Out of Service

Intoxilyzer **64-002505 (#11)** was taken out of service for repair on **October 7, 2004**.

I cannot testify that it was working properly and accurately between the dates of **September 9, 2004** and **October 7, 2004**.

This instrument was located DUI HQ.