

1086
ls

EXHIBIT Q

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES *AS*
INTOXILYZER MODEL 5000*

AGENCY Tucson Police QA SPECIALIST Irene Doran 40331
(Print Name)

INTOXILYZER SERIAL # 64-002516 (H12) LOCATION ODS
DATE December 20 *ls* 2004 TIME 0957

ls Display reads "READY TO START" or "PUSH BUTTON TO START TEST."

DIAGNOSTIC TESTS

- ls* DVM Test check. Setting should be between .200 and .600. Mode selection switch S2 on, S1 and S3 off or keyboard menu selection "H". Reading is .386.
- ls* Display Test check. Mode selection switch S1 on and S2 and S3 off or keyboard menu selection "V".
- ls* Printer Test check. Mode selection switch S1, S2, S3 off or keyboard menu selection "P".
- ls* Clock time check. Mode selection switch S10 on or keyboard menu selection "E".
- ls* Date check. Mode selection switch S11 on or keyboard menu selection "E".

OPERATIONAL TESTS

- ls* Alcohol free subject Test result 0.000 AC
 - ls* Error recognition logic system functioning
Invalid test printed *ls*
 - ls* Proper sample recognition system
Invalid sample printed *ls*
Deficient sample printed *ls*
 - N/A Completeness of sample purge with collection tube**
 - ls* Calibration standard 0.080 AC
Results: 0.0754C + 0.0754C
- Instrument operating properly and accurately. YES NO

COMMENTS: 0954-0958

SIGNATURE Irene Doran

0. N/A alcohol concentration standard collected for subsequent analysis.**

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD
**THIS STEP IS ONLY REQUIRED IF THIS DEVICE IS BEING USED FOR SAMPLE CAPTURE.

386
25

Tucson Police Department Crime Laboratory

24

FEDERAL SIGNAL CORP. / CMI INC
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 64-002516
12/20/04

#12

TEST	BAC VALUE	TIME
AIR BLANK	.000	10:14
CAL. CHECK	.075	10:14
AIR BLANK	.000	10:15

SUBJECT'S NAME

ODS

INSTRUMENT LOCATION

TIME FIRST OBSERVED

40331

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS



INTOXILYZER® INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

FEDERAL SIGNAL CORP. / CMI INC
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 64-002516
12/20/04

#12

TEST	BAC VALUE	TIME
AIR BLANK	.000	10:13
CAL. CHECK	.075	10:13
AIR BLANK	.000	10:13

SUBJECT'S NAME

ODS

INSTRUMENT LOCATION

TIME FIRST OBSERVED

40331

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS



INTOXILYZER® INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

496
lv

SN 64-002516
E47.02
INVALID TEST
SAMPLE INTRODUCED
AT IMPROPER TIME.

12/20/04
10:03

#12

SUBJECT'S NAME	
TIME FIRST OBSERVED	INSTRUMENT LOCATION
OPERATOR	
ADDITIONAL INFORMATION AND / OR REMARKS	

Handwritten entries in the form:
SUBJECT'S NAME: ODS
OPERATOR: 46331

INTOXILYZER[®] INSTRUMENT PRINTER CARD

© 1986 by CMI INC.



PA

EXHIBIT P

5 2 6
ls

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 5000*

STANDARD CALIBRATION CHECK PROCEDURE

104

Agency Tucson Police Department Date 12/20/04 Time 0957
Intoxilyzer Serial# 64-002516 (#12) Location ODS
QA Specialist Gene Dora 46331
(Print Name)

- (A) 1. Pour a standard alcohol solution of known value into a clean dry simulator jar and assemble the simulator. Insure that a tight seal has been made. Standard value: 0.080 AC
- (A) 2. Turn on the simulator and allow the temperature to reach $34^{\circ}\text{C} \pm .2^{\circ}\text{C}$.
- (A) 3. Set Intoxilyzer mode selection in the ACA mode by switching mode selection switch #9 on or selecting "C" on keyboard menu.
- (A) 4. Attach simulator to the simulator entrance port on the Intoxilyzer.
- (A) 5. Intoxilyzer 5000 display reads "READY TO START" or "PUSH BUTTON".
- (A) 6. Push Start Test button or press enter on keyboard.
- (A) 7. Insert card in response to display.
- (A) 8. Air blank completed.
- (A) 9. Calibration check completed. Test results 0.075 + 0.075 AC
- (A) 10. Air blank completed.
- (A) 11. When display reads Test Complete remove evidence card. Attach the card to the completed checklist.
- (A) 12. Return mode selection switch #9 to off position after all calibration checks are complete or type Q and enter on keyboard.

SIGNATURE Gene Dora

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD

FEDERAL SIGNAL CORP. / CMI INC
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 64-002516
12/20/04

#12

TEST	BAC VALUE	TIME
AIR BLANK	.000	10:03
INVALID. SAMPLE	.XXX	10:04
AIR BLANK	.000	10:04

SUBJECT'S NAME

ODS

INSTRUMENT LOCATION

TIME FIRST OBSERVED

46331

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

sweet breath

CMI INC

INTOXILYZER INSTRUMENT PRINTER CARD

© 1986 by CMI INC

PS

FEDERAL SIGNAL CORP. / CMI INC
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 64-002516
12/20/04

#12

TEST	BAC VALUE	TIME
AIR BLANK	.000	10:05
*SUBJECT TEST	.117	10:07
AIR BLANK	.000	10:07

* DEFICIENT SAMPLE - VALUE
PRINTED WAS HIGHEST OBTAINED.

SUBJECT'S NAME

ODS

INSTRUMENT LOCATION

TIME FIRST OBSERVED

46331

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

CMI INC

INTOXILYZER INSTRUMENT PRINTER CARD

© 1986 by CMI INC

PS

686
LS