

EXHIBIT Q

10/6
JLc

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
 A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES *fm*
 INTOXILYZER MODEL 5000*

AGENCY Tucson Police QA SPECIALIST Jennifer Combs 46660
 (Print Name)
 INTOXILYZER SERIAL # 66-001759 LOCATION Dui HQ
 DATE October 7 16 2004 TIME 1150

Display reads "READY TO START" or "PUSH BUTTON TO START TEST."

DIAGNOSTIC TESTS

- DVM Test check. Setting should be between .200 and .600. Mode selection switch S2 on, S1 and S3 off or keyboard menu selection "H". Reading is .357.
- Display Test check. Mode selection switch S1 on and S2 and S3 off or keyboard menu selection "V".
- Printer Test check. Mode selection switch S1, S2, S3 off or keyboard menu selection "P".
- ~~*~~ Clock time check. Mode selection switch S10 on or keyboard menu selection "E".
- Date check. Mode selection switch S11 on or keyboard menu selection "E".

OPERATIONAL TESTS

Alcohol free subject Test result 0.000 AC

Error recognition logic system functioning

Invalid test printed

Proper sample recognition system

Invalid sample printed

Deficient sample printed

N/A Completeness of sample purge with collection tube**

Calibration standard 0.080 AC

Results: 0.081 AC & 0.079 AC

Instrument operating properly and accurately. YES NO

COMMENTS: * moved ahead 1147 → 1150

SIGNATURE

Jennifer L. Combs

Q. N/A alcohol concentration standard collected for subsequent analysis.**

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD

**THIS STEP IS ONLY REQUIRED IF THIS DEVICE IS BEING USED FOR SAMPLE CAPTURE.

3 of 6
file

Tucson Police Department Crime Laboratory

TUCSON POLICE DEPARTMENT
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-001753
10/07/84 11:58

TEST	BRAC VALUE	TIME
AIR BLANK	:000	11:57
CAL. CHECK	:081	11:57
AIR BLANK	:000	11:58

TUCSON POLICE DEPARTMENT
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-001759
10/07/84 11:58

TEST	BRAC VALUE	TIME
AIR BLANK	:000	11:57
CAL. CHECK	:079	11:57
AIR BLANK	:000	11:58

SUBJECT'S NAME	Du Ha
INSTRUMENT LOCATION	INSTRUMENT LOCATION
TIME FIRST OBSERVED	TIME FIRST OBSERVED
OPERATOR	OPERATOR
ADDITIONAL INFORMATION AND/OR REMARKS	

SUBJECT'S NAME	Du Ha
INSTRUMENT LOCATION	INSTRUMENT LOCATION
TIME FIRST OBSERVED	TIME FIRST OBSERVED
OPERATOR	OPERATOR
ADDITIONAL INFORMATION AND/OR REMARKS	

406
file

SN 66-001759
E767.14
INVALID TEST

19/07/84
12:04

SAMPLE INTRODUCED
AT IMPROPER TIME.

SUBJECT'S NAME

Dui HQ

TIME FIRST OBSERVED

INSTRUMENT LOCATION

46160

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

INTOXILYZER INSTRUMENT PRINTER CARD

CMI INC

© 1986 by CMI INC

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 5000*STANDARD CALIBRATION CHECK PROCEDURE *BS*

Agency Tucson Police Department Date 10/7/04 Time 1150
 Intoxilyzer Serial# 16-001759 Location DUI HQ
 QA Specialist Jennifer Coombs 4601160
 (Print Name)

- (✓) 1. Pour a standard alcohol solution of known value into a clean dry simulator jar and assemble the simulator. Insure that a tight seal has been made.
Standard value: 0.080 AC
- (✓) 2. Turn on the simulator and allow the temperature to reach $34^{\circ}\text{C} \pm .2^{\circ}\text{C}$.
- (✓) 3. Set Intoxilyzer mode selection in the ACA mode by switching mode selection switch #9 on or selecting "C" on keyboard menu.
- (✓) 4. Attach simulator to the simulator entrance port on the Intoxilyzer.
- (✓) 5. Intoxilyzer 5000 display reads "READY TO START" or "PUSH BUTTON".
- (✓) 6. Push Start Test button or press enter on keyboard.
- (✓) 7. Insert card in response to display.
- (✓) 8. Air blank completed.
- (✓) 9. Calibration check completed.
Test results 0.081 AC & 0.079 AC
- (✓) 10. Air blank completed.
- (✓) 11. When display reads Test Complete remove evidence card. Attach the card to the completed checklist.
- (✓) 12. Return mode selection switch #9 to off position after all calibration checks are complete or type Q and enter on keyboard.

SIGNATURE

Jennifer L. Coombs

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD

TUCSON POLICE DEPARTMENT
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-001759
12:05
10/07/84

TEST BRAC VALUE TIME
AIR BLANK .000 12:05
INVALID SAMPLE .000 12:05
AIR BLANK .000 12:05

BEGIN NEW DEPRIVATION PERIOD

TUCSON POLICE DEPARTMENT
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-001759
12:09
10/07/84

TEST BRAC VALUE TIME
AIR BLANK .000 12:05
*SUBJECT TEST .065 12:09
AIR BLANK .000 12:09

* DEFICIENT SAMPLE - VALUE
PRINTED WAS HIGHEST OBTAINED.

SUBJECT'S NAME	Dui HQ
INSTRUMENT LOCATION	40000
TIME FIRST OBSERVED	
OPERATOR	
ADDITIONAL INFORMATION AND/OR REMARKS	
Sweet Breath	

INTOXILYZER[®] INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

CMI^{inc}

INTOXILYZER[®] INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

CMI^{inc}

6 of 6
file