

186
D5

EXHIBIT Q

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES *154*
INTOXILYZER MODEL 5000*

AGENCY Tucson Police QA SPECIALIST Irene Doran (46331)
(Print Name)
INTOXILYZER SERIAL #66-00404 (#17) LOCATION DUI HQ
DATE April 22 192007 TIME 1348

JK Display reads "READY TO START" or "PUSH BUTTON TO START TEST."

DIAGNOSTIC TESTS

- JK* DVM Test check. Setting should be between .200 and .600. Mode selection switch S2 on, S1 and S3 off or keyboard menu selection "H". Reading is .294.
- JK* Display Test check. Mode selection switch S1 on and S2 and S3 off or keyboard menu selection "V".
- JK* Printer Test check. Mode selection switch S1, S2, S3 off or keyboard menu selection "P".
- JK* Clock time check. Mode selection switch S10 on or keyboard menu selection "E".
- JK* Date check. Mode selection switch S11 on or keyboard menu selection "E".

OPERATIONAL TESTS

- JK* Alcohol free subject Test result 0.000 AC
JK Error recognition logic system functioning
 Invalid test printed
JK Proper sample recognition system
 Invalid sample printed *JK*
 Deficient sample printed *JK*
N/A Completeness of sample purge with collection tube**
JK Calibration standard 0.080 AC
Results: 0.077 AC + 0.077 AC
Instrument operating properly and accurately. YES NO

COMMENTS:

SIGNATURE

Irene Doran

ON/A alcohol concentration standard collected for subsequent analysis.**

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD

**THIS STEP IS ONLY REQUIRED IF THIS DEVICE IS BEING USED FOR
SAMPLE CAPTURE.

CHI INC
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-004614
004/22/94

TEST	BAC	VALUE	TIME
HAIR BLANK	.000		13:56
SUBJECT TEST	.000		13:56
HAIR BLANK	.000		13:57

DvM = 296

SUBJECT'S NAME	DUST HOG
TIME FIRST OBSERVED	INSTRUMENTATION

11/11/10 INSTRUMENT LOCATION

(46-331)

OPERATOR _____

ADDITIONAL INFORMATION AND / OR REMARKS

INTOXILYZER® INSTRUMENT PRINTER CARD
© 1986 by CMI INC.

INTOXILYZER® INSTRUMENT PRINTER CARD
© 1986 by CMI, INC.

24

86
so

Tucson Police Department Crime Laboratory *fm*

386
S

CHI INC
#17
94/22/04
SH 5-004014
#2 MODEL 5000 - ALCOHOL ANALYZER
#1 TOXILIZER - ALCOHOL ANALYZER
SH 66-004014

TEST	BAC	VALUE	TIME	TEST	BAC	VALUE	TIME
AIR BLANK	.000		13:49	AIR BLANK	.000		13:50
CAL. CHECK	.077		13:49	CAL. CHECK	.977		13:51
AIR BLANK	.000		13:50	AIR BLANK	.000		13:51

SUBJECT'S NAME	DUE HQ	INSTRUMENT LOCATION
TIME FIRST OBSERVED	46331	OPERATOR
ADDITIONAL INFORMATION AND/OR REMARKS		
SUBJECT'S NAME	DUE HQ	INSTRUMENT LOCATION
TIME FIRST OBSERVED	46331	OPERATOR
ADDITIONAL INFORMATION AND/OR REMARKS		

© 1986 by CMN.

© 1986 by CMC Inc.

© 1986 by CMC.

486
PS

SN 66-004014 #17 04/22/84
E647.85 14:04
INVALID TEST
SAMPLE INTRODUCED
AT IMPROPER TIME.

SUBJECT'S NAME _____

TIME FIRST OBSERVED _____ INSTRUMENT LOCATION _____
DUT HQ

OPERATOR _____

ADDITIONAL INFORMATION AND/OR REMARKS

INTOXILYZER® INSTRUMENT PRINTER CARD

CMI
INC

© 1986 by CMI INC

PS

586
80

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 5000*

184

STANDARD CALIBRATION CHECK PROCEDURE

Agency Tucson Police Department Date 4/22/04 Time 1348
Intoxilyzer Serial # 66-004014 (#17) Location DUI HQ
QA Specialist Irene Doran (44331)
(Print Name)

- (1. Pour a standard alcohol solution of known value into a clean dry simulator jar and assemble the simulator. Insure that a tight seal has been made.
Standard value: 0.080 AC
- (2. Turn on the simulator and allow the temperature to reach $34^{\circ}\text{C} \pm .2^{\circ}\text{C}$.
- (3. Set Intoxilyzer mode selection in the ACA mode by switching mode selection switch #9 on or selecting "C" on keyboard menu.
- (4. Attach simulator to the simulator entrance port on the Intoxilyzer.
- (5. Intoxilyzer 5000 display reads "READY TO START" or "PUSH BUTTON".
- (6. Push Start Test button or press enter on keyboard.
- (7. Insert card in response to display.
- (8. Air blank completed.
- (9. Calibration check completed.
Test results 0.077 & 0.077 AC
- (10. Air blank completed.
- (11. When display reads Test Complete remove evidence card. Attach the card to the completed checklist.
- (12. Return mode selection switch #9 to off position after all calibration checks are complete or type Q and enter on keyboard.

SIGNATURE

Irene Doran

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD

CBL 100
 INTOXILYZER - ALCOHOL ANALYZER
 AZ MODEL 5000 SN 66-0001014
 04/22/04

TEST	BAC	TIME
AIR BLANK	.000	14:01
INVALID SAMPLE	.XXX	14:02
AIR BLANK	.000	14:02

#17

CHILLING
INTOXILIZER - ALCOHOL ANALYZER
EZ MODEL 5000 SH 6-004014
814/22764

TEST BAC VALUE TIME

AIR BLANK	000	14:01	AIR BLANK	000	13:58
INVALID SAMPLE	XXX	14:02	*SUBJECT TEST	068	14:00
AIR BLANK	000	14:02	AIR BLANK	000	14:01

* DEFICIENT SAMPLE - VALUE PRINTED WAS HIGHEST OBTAINED.

40.331

OPERATOR _____

~ ADDITIONAL INFORMATION AND / OR REMARKS
Bianca

(44-331)
OPERATOR
ADDITIONAL INFORMATION AND/OR RE

(44331) **OPERATOR**
ADDITIONAL INFORMATION AND OR REMARKS

40.331

OPERATOR _____

~ ADDITIONAL INFORMATION AND / OR REMARKS
Bianca

(44331) **OPERATOR**
ADDITIONAL INFORMATION AND OR REMARKS

40.331

OPERATOR _____

~ ADDITIONAL INFORMATION AND / OR REMARKS
Bianca

INTOXLYZER[®] INSTRUMENT PRINTER CARD
© 1996 LSI CLAIMS

110

INTOXILYZER® INSTRUMENT PRINTER CARD

三

© 1986 by CM INC

© 1986 by CM INC