

EXHIBIT Q

P. 1 of 6
JL

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
 A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES
 INTOXILYZER MODEL 5000*

AGENCY Tucson Police QA SPECIALIST Jennifer Coombs 46160
 (Print Name)
 INTOXILYZER SERIAL #106-00404 (#17) LOCATION DUI HQ
 DATE July 14 19 2004 TIME 1057

Display reads "READY TO START" or "PUSH BUTTON TO START TEST."

DIAGNOSTIC TESTS

- DVM Test check. Setting should be between .200 and .600. Mode selection switch S2 on, S1 and S3 off or keyboard menu selection "H". Reading is .342.
- Display Test check. Mode selection switch S1 on and S2 and S3 off or keyboard menu selection "V".
- Printer Test check. Mode selection switch S1, S2, S3 off or keyboard menu selection "P".
- Clock time check. Mode selection switch S10 on or keyboard menu selection "E".
- Date check. Mode selection switch S11 on or keyboard menu selection "E".

OPERATIONAL TESTS

- Alcohol free subject Test result 0.000 AC
- Error recognition logic system functioning
- Invalid test printed
- Proper sample recognition system
- Invalid sample printed✓
- Deficient sample printed✓
- N/A Completeness of sample purge with collection tube**
- Calibration standard 0.080 AC
- Results: 0.016AC & 0.018AC

Instrument operating properly and accurately. YES NO

COMMENTS:

SIGNATURE

Jennifer L. Coombs

O. N/A alcohol concentration standard collected for subsequent analysis.**

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD

**THIS STEP IS ONLY REQUIRED IF THIS DEVICE IS BEING USED FOR
 SAMPLE CAPTURE.

#17

SN 66-0014014

07/14/04

E647.05 1:01

EEEEEEEEEEEEEEEEEEEEEEEE

EEEEEEEEEEEEEEEEEEEEEEEE

EEEEEEEEEEEEEEEEEEEEEEEE

EEEEEEEEEEEEEEEEEEEEEEEE

EEEEEEEEEEEEEEEEEEEEEEEE

EEEEEEEEEEEEEEEEEEEEEEEE

HHHHHHHHHHHHHHHHHHHHHH

CMI INC
 INTOXILYZER - ALCOHOL ANALYZER
 AZ MODEL 5000 SN 66-004014
 07/14/04

TEST	BAC	VALUE	TIME
AIR BLANK	.000		11:08
SUBJECT TEST	.000		11:08
AIR BLANK	.000		11:09

DVM = .342

DUI HQ

INSTRUMENT LOCATION

DUI HQ

INSTRUMENT LOCATION

419100

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

419100

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

CMI

INTOXILYZER INSTRUMENT PRINTER CARD

© 1986 by CMI INC

INTOXILYZER INSTRUMENT PRINTER CARD

© 1986 by CMI INC

CMI

INTOXILYZER INSTRUMENT PRINTER CARD

© 1986 by CMI INC

P. 3 of 6
JLC

Tucson Police Department Crime Laboratory

#17

#17

CMI INC
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-004014
07/14/84

CMI INC
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-004014
07/14/84

TEST BAC VALUE TIME

AIR BLANK	.000	11:03
CAL. JCHECK	.076	11:03
AIR BLANK	.000	11:03

TEST BAC VALUE TIME

AIR BLANK	.000	11:04
CAL. CHECK	.078	11:04
AIR BLANK	.000	11:05

SUBJECT'S NAME	DUY HU
TIME FIRST OBSERVED	
INSTRUMENT LOCATION	

SUBJECT'S NAME	DUY HU
TIME FIRST OBSERVED	
INSTRUMENT LOCATION	

SUBJECT'S NAME	DUY HU
TIME FIRST OBSERVED	
INSTRUMENT LOCATION	

SUBJECT'S NAME	DUY HU
TIME FIRST OBSERVED	
INSTRUMENT LOCATION	

SUBJECT'S NAME	
TIME FIRST OBSERVED	
INSTRUMENT LOCATION	

SUBJECT'S NAME	
TIME FIRST OBSERVED	
INSTRUMENT LOCATION	

INTOXILYZER® INSTRUMENT PRINTER CARD
© 1986 by CMI INC.

cmi

INTOXILYZER® INSTRUMENT PRINTER CARD
© 1986 by CMI INC.

© 1986 by CMI INC.

cmi

P. 4 of 6
JLC

SN 66-004014
E647.05
INVALID TEST
SAMPLE INTRODUCED
AT IMPROPER TIME.

07/14/84
11:10

#17

SUBJECT'S NAME

DUI HQ

TIME FIRST OBSERVED

INSTRUMENT LOCATION

461160

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

INTOXILYZER® INSTRUMENT PRINTER CARD

CMI

© 1986 by CMI INC.

EXHIBIT P

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 5000*

STANDARD CALIBRATION CHECK PROCEDURE *AP*

Agency Tucson Police Department Date 7/14/04 Time 1057
Intoxilyzer Serial# 16-00404 (#1) Location DUI HQ
QA Specialist Jennifer Coombs 461160
(Print Name)

- (✓) 1. Pour a standard alcohol solution of known value into a clean dry simulator jar and assemble the simulator. Insure that a tight seal has been made.
Standard value: 0.080 AC
- (✓) 2. Turn on the simulator and allow the temperature to reach 34°C ± .2°C.
- (✓) 3. Set Intoxilyzer mode selection in the ACA mode by switching mode selection switch #9 on or selecting "C" on keyboard menu.
- (✓) 4. Attach simulator to the simulator entrance port on the Intoxilyzer.
- (✓) 5. Intoxilyzer 5000 display reads "READY TO START" or "PUSH BUTTON".
- (✓) 6. Push Start Test button or press enter on keyboard.
- (✓) 7. Insert card in response to display.
- (✓) 8. Air blank completed.
- (✓) 9. Calibration check completed.
Test results 0.0764 ± 0.078 AC
- (✓) 10. Air blank completed.
- (✓) 11. When display reads Test Complete remove evidence card. Attach the card to the completed checklist.
- (✓) 12. Return mode selection switch #9 to off position after all calibration checks are complete or type Q and enter on keyboard.

SIGNATURE

Jennifer L. Coombs

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD

#17
CMI INC.
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-004014
07/14/84

TEST	BAC	VALUE	TIME
AIR BLANK	.000	11:11	
INVALID SAMPLE	.XXX	11:11	
AIR BLANK	.000	11:11	

* DEFICIENT SAMPLE - VALUE
PRINTED WAS HIGHEST OBTAINED.

CMI INC.
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-004014
07/14/84

TEST	BAC	VALUE	TIME
AIR BLANK	.000	11:12	
*SUBJECT TEST	.066	11:14	
AIR BLANK	.000	11:15	

* DEFICIENT SAMPLE - VALUE
PRINTED WAS HIGHEST OBTAINED.

SUBJECT'S NAME _____
TIME FIRST OBSERVED _____
INSTRUMENT LOCATION _____
OPERATOR _____
ADDITIONAL INFORMATION AND/OR REMARKS
Du Ha HQ
46160
Sweet Breath

SUBJECT'S NAME _____
TIME FIRST OBSERVED _____
INSTRUMENT LOCATION _____
OPERATOR _____
ADDITIONAL INFORMATION AND/OR REMARKS
Du Ha HQ
46160

INTOXILYZER® INSTRUMENT PRINTER CARD
© 1986 by CMI INC.

CMI
INC.

INTOXILYZER® INSTRUMENT PRINTER CARD

CMI
INC.

© 1986 by CMI INC.

© 1986 by CMI INC.