

EXHIBIT Q

1 of 6
jlc

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 5000*

SM

AGENCY Tucson Police QA SPECIALIST Jennifer Coombs 46160
(Print Name)

INTOXILYZER SERIAL # 66-004814 LOCATION 000
DATE September 1 15 2004 TIME 1543

Display reads "READY TO START" or "PUSH BUTTON TO START TEST."

DIAGNOSTIC TESTS

- DVM Test check. Setting should be between .200 and .600. Mode selection switch S2 on, S1 and S3 off or keyboard menu selection "H". Reading is .358.
- Display Test check. Mode selection switch S1 on and S2 and S3 off or keyboard menu selection "V".
- Printer Test check. Mode selection switch S1, S2, S3 off or keyboard menu selection "P".
- Clock time check. Mode selection switch S10 on or keyboard menu selection "E".
- Date check. Mode selection switch S11 on or keyboard menu selection "E".

OPERATIONAL TESTS

- Alcohol free subject Test result 0.000 AC
 - Error recognition logic system functioning
Invalid test printed
 - Proper sample recognition system
Invalid sample printed ✓
Deficient sample printed ✓
 - N/A Completeness of sample purge with collection tube**
 - Calibration standard 0.050 AC
Results: 0.080AC & 0.080AC
- Instrument operating properly and accurately. YES NO

COMMENTS: _____

SIGNATURE Jennifer L. Coombs

0. N/A alcohol concentration standard collected for subsequent analysis.**

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD
**THIS STEP IS ONLY REQUIRED IF THIS DEVICE IS BEING USED FOR SAMPLE CAPTURE.

3 of 6
jlc

Tucson Police Department Crime Laboratory

TUCSON POLICE DEPARTMENT
 INTOXILYZER - ALCOHOL ANALYZER
 AZ MODEL 5000 SN 66-004814
 09/01/04 16:00

TEST BRAC VALUE TIME
 AIR BLANK .000 15:59
 CAL. CHECK .080 16:00
 AIR BLANK .000 16:00

SUBJECT'S NAME
ODD

INSTRUMENT LOCATION

TIME FIRST OBSERVED

OPERATOR
46160

ADDITIONAL INFORMATION AND/OR REMARKS



INTOXILYZER INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

TUCSON POLICE DEPARTMENT
 INTOXILYZER - ALCOHOL ANALYZER
 AZ MODEL 5000 SN 66-004814
 09/01/04 15:58

TEST BRAC VALUE TIME
 AIR BLANK .000 15:58
 CAL. CHECK .080 15:58
 AIR BLANK .000 15:58

SUBJECT'S NAME

ODD

INSTRUMENT LOCATION

TIME FIRST OBSERVED

OPERATOR
46160

ADDITIONAL INFORMATION AND/OR REMARKS



INTOXILYZER INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

4 of 6
file

SN 66-004814
E767.14
INVALID TEST

09/01/04
16:04

SAMPLE INTRODUCED
AT IMPROPER TIME.

SUBJECT'S NAME

0DD

TIME FIRST OBSERVED

INSTRUMENT LOCATION

46160

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

INTOXILYZER[®] INSTRUMENT PRINTER CARD



© 1986 by CMI INC.

EXHIBIT P

5 of 6
fle

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 5000*

STANDARD CALIBRATION CHECK PROCEDURE *BM*

Agency Tucson Police Department Date 9/1/04 Time 1543
Intoxilyzer Serial# 166-004814 Location ODD
QA Specialist Jennifer Coombs 46160
(Print Name)

- (✓) 1. Pour a standard alcohol solution of known value into a clean dry simulator jar and assemble the simulator. Insure that a tight seal has been made. Standard value: 0.080 AC
- (✓) 2. Turn on the simulator and allow the temperature to reach 34°C±.2°C.
- (✓) 3. Set Intoxilyzer mode selection in the ACA mode by switching mode selection switch #9 on or selecting "C" on keyboard menu.
- (✓) 4. Attach simulator to the simulator entrance port on the Intoxilyzer.
- ✓ ✓ 5. Intoxilyzer 5000 display reads "READY TO START" or "PUSH BUTTON".
- ✓ ✓ 6. Push Start Test button or press enter on keyboard.
- ✓ ✓ 7. Insert card in response to display.
- ✓ ✓ 8. Air blank completed.
- ✓ ✓ 9. Calibration check completed. Test results 0.080AC & 0.080 AC
- ✓ ✓ 10. Air blank completed.
- ✓ ✓ 11. When display reads Test Complete remove evidence card. Attach the card to the completed checklist.
- (✓) 12. Return mode selection switch #9 to off position after all calibration checks are complete or type Q and enter on keyboard.

SIGNATURE Jennifer L. Coombs

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD

TUCSON POLICE DEPARTMENT
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-004814
09/01/04 16:05

TEST BRAC VALUE TIME
AIR BLANK .000 16:04
INVALID SAMPLE .MOA 16:05
AIR BLANK .000 16:05

BEGIN NEW DEPRIVATION PERIOD

SUBJECT'S NAME

TIME FIRST OBSERVED

INSTRUMENT LOCATION

0DD

46160

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

Sweet Breath

INTOXILYZER[®] INSTRUMENT PRINTER CARD

© 1986 by CMI INC.



TUCSON POLICE DEPARTMENT
INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-004814
09/01/04 16:10

TEST BRAC VALUE TIME
AIR BLANK .000 16:05
*SUBJECT TEST .079 16:09
AIR BLANK .000 16:09

* DEFICIENT SAMPLE - VALUE
PRINTED WAS HIGHEST OBTAINED.

SUBJECT'S NAME

TIME FIRST OBSERVED

INSTRUMENT LOCATION

0DD

46160

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

INTOXILYZER[®] INSTRUMENT PRINTER CARD

© 1986 by CMI INC.



606
fle