

EXHIBIT Q

1 of 6
file

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
 A.A.C. R9-14-404 (A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES
 INTOXILYZER MODEL 5000*

AGENCY Tucson Police QA SPECIALIST Jennifer Coombs 46160
 (Print Name)
 INTOXILYZER SERIAL # 160-005004 LOCATION LTB
 DATE October 10 19 2004 TIME 1440

- Display reads "READY TO START" or "PUSH BUTTON TO START TEST."

DIAGNOSTIC TESTS

- DVM Test check. Setting should be between .200 and .600. Mode selection switch S2 on, S1 and S3 off or keyboard menu selection "H". Reading is .413.
- Display Test check. Mode selection switch S1 on and S2 and S3 off or keyboard menu selection "V".
- Printer Test check. Mode selection switch S1, S2, S3 off or keyboard menu selection "P".
- Clock time check. Mode selection switch S10 on or keyboard menu selection "E".
- Date check. Mode selection switch S11 on or keyboard menu selection "E".

OPERATIONAL TESTS

- Alcohol free subject Test result 0.000 AC
- Error recognition logic system functioning
 Invalid test printed
- Proper sample recognition system
 Invalid sample printed
 Deficient sample printed
- N/A Completeness of sample purge with collection tube**
 Calibration standard 0.080 AC
 Results: 0.071AC & 0.074AC
- Instrument operating properly and accurately. YES NO

COMMENTS: _____

SIGNATURE

Jennifer L. Coombs

0. N/A alcohol concentration standard collected for subsequent analysis.**

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD

**THIS STEP IS ONLY REQUIRED IF THIS DEVICE IS BEING USED FOR
 SAMPLE CAPTURE.

TEST	BRAC VALUE	TIME
AIR BLANK	.000	14:49
SUBJECT TEST	.000	14:49
AIR BLANK	.000	14:50

$$DV/H = .4/3$$

SUBJECT'S NAME	LAb	INSTRUMENT LOCATION
TIME FIRST OBSERVED		OPERATOR
ADDITIONAL INFORMATION AND/OR REMARKS Yellow		

CMI INC.
INTOXILYZER[®] INSTRUMENT PRINTER CARD
© 1986 by CMI INC.

INTOXILYZER® INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

三

Zoof
Jec

Tucson Police Department Crime Laboratory

3 of 6
JLC

INTOXILYZER - ALCOHOL ANALYZER	
AZ MODEL 5000	SH 66-005004
10/10/84	14:46
TEST	BRAC VALUE TIME
AIR BLANK	:000 14:45
CAL. CHECK	:077 14:46
AIR BLANK	:000 14:46
TIME FIRST OBSERVED	INSTRUMENT LOCATION
SUBJECT'S NAME	
OPERATOR	ADDITIONAL INFORMATION AND/OR REMARKS
4Hello 1 AB	

4 of 6
JLc

SN 66-005004 E767.14 INVALID TEST	10/10/84 14:51
SAMPLE INTRODUCED AT IMPROPER TIME.	
<hr/> <hr/>	
SUBJECT'S NAME	
LAB	
TIME FIRST OBSERVED	INSTRUMENT LOCATION
<hr/>	
4/11/80	
OPERATOR	
ADDITIONAL INFORMATION AND/OR REMARKS	
<hr/> <hr/> <hr/>	
INTOXILYZER [®] INSTRUMENT PRINTER CARD	
© 1986 by CMI INC.	

CMI
INC

5 of 6
file

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
 A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES
 INTOXILYZER MODEL 5000*STANDARD CALIBRATION CHECK PROCEDURE 104

Agency Tucson Police Department Date 10/10/04 Time 1440
 Intoxilyzer Serial# 1dp-005004 Location LAB
 QA Specialist Jennifer Combs 46660
 (Print Name)

- (✓) 1. Pour a standard alcohol solution of known value into a clean dry simulator jar and assemble the simulator. Insure that a tight seal has been made.
 Standard value: 0.080 AC
- (✓) 2. Turn on the simulator and allow the temperature to reach $34^{\circ}\text{C} \pm .2^{\circ}\text{C}$.
- (✓) 3. Set Intoxilyzer mode selection in the ACA mode by switching mode selection switch #9 on or selecting "C" on keyboard menu.
- (✓) 4. Attach simulator to the simulator entrance port on the Intoxilyzer.
- ✓(✓) 5. Intoxilyzer 5000 display reads "READY TO START" or "PUSH BUTTON".
- ✓(✓) 6. Push Start Test button or press enter on keyboard.
- ✓(✓) 7. Insert card in response to display.
- ✓(✓) 8. Air blank completed.
- ✓(✓) 9. Calibration check completed.
 Test results 0.071 AC & 0.076 AC
- ✓(✓) 10. Air blank completed.
- ✓(✓) 11. When display reads Test Complete remove evidence card. Attach the card to the completed checklist.
- (✓) 12. Return mode selection switch #9 to off position after all calibration checks are complete or type Q and enter on keyboard.

SIGNATURE

Jennifer L. Combs

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD

INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SH 66-065604
10/10/84 14:52

TEST BRAC VALUE TIME
AIR BLANK : 000 14:51
INVALID SAMPLE : .004 14:52
AIR BLANK : .000 14:52

BEGIN NEW DEPRIVATION PERIOD

INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SH 66-065604
10/10/84 14:56

TEST BRAC VALUE TIME
AIR BLANK : 000 14:53
*SUBJECT TEST : .073 14:56
AIR BLANK : .000 14:56

*DEFICIENT SAMPLE - VALUE
PRINTED WAS HIGHEST OBTAINED.

SUBJECT'S NAME _____

Lab

INSTRUMENT LOCATION

INSTRUMENT LOCATION

TIME FIRST OBSERVED

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

410160

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

Sweet Breath

INTOXILYZER[®] INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

CMI[®] Inc

INTOXILYZER[®] INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

CMI[®] Inc

6 of 6
JLc