



CERTIFICATION by CUSTODIAN of RECORDS  
Tucson Police Department Crime Laboratory

154

I hereby certify that these documents are a true and correct copy of the periodic maintenance\* for the Intoxilyzer 5000 Serial Number 166-005004 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.

A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R9-14-404.A.3, for a 0.080AC solution are 0.070AC to 0.090AC.

Date: 12/29/05 Cal 1 0.079 AC Cal 2 0.078 AC

Date: 1/26/06 Cal 1 0.079 AC Cal 2 0.080 AC

Date: N/A Cal 1                      AC Cal 2                      AC

Between the dates of 12/29/05 and 1/26/06 the instrument was operating properly and accurately.

The instrument was taken out of service for repair on N/A.  
The instrument was returned to service on N/A.

Jennifer L. Coombs  
Jennifer L. Coombs

1/27/06

ADPS QAS Permit # QAS 29886

SUBSCRIBED AND SWORN to before me this 27, day of

January, 2006.



OFFICIAL SEAL  
KATHY A. KINNEY  
NOTARY PUBLIC-ARIZONA  
PIMA COUNTY

My Comm. Exp. Sept. 4, 2008

Kathy A. Kinney  
Notary Public

\*Includes documents specified in A.A.C. R9-14-404.A.6

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

## EXHIBIT Q

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY  
A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 5000\*AGENCY Tucson Police QA SPECIALIST Jennifer Coombs 4/6/60  
(Print Name)INTOXILYZER SERIAL # 66-005004 LOCATION DL11 HQ  
DATE January 26 19 2006 TIME 1150

Display reads "READY TO START" or "PUSH BUTTON TO START  
TEST."

## DIAGNOSTIC TESTS

- ☒ DVM Test check. Setting should be between .200 and .600. Mode selection switch S2 on, S1 and S3 off or keyboard menu selection "H". Reading is .366.
- ☒ Display Test check. Mode selection switch S1 on and S2 and S3 off or keyboard menu selection "V".
- ☒ Printer Test check. Mode selection switch S1, S2, S3 off or keyboard menu selection "P".
- ☒ Clock time check. Mode selection switch S10 on or keyboard menu selection "E".
- ☒ Date check. Mode selection switch S11 on or keyboard menu selection "E".

## OPERATIONAL TESTS

- ☒ Alcohol free subject Test result 0.000 AC
- ☒ Error recognition logic system functioning  
Invalid test printed
- ☒ Proper sample recognition system  
Invalid sample printed ☒  
Deficient sample printed ☒
- N/A Completeness of sample purge with collection tube\*\*
- ☒ Calibration standard 0.080 AC  
Results: 0.079 AC & 0.080 AC
- Instrument operating properly and accurately. YES ☒ NO ☐

COMMENTS:

SIGNATURE

Jennifer L. Coombs

0. N/A alcohol concentration standard collected for subsequent analysis.\*\*

\*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD  
\*\*THIS STEP IS ONLY REQUIRED IF THIS DEVICE IS BEING USED FOR  
SAMPLE CAPTURE.

DHS/DSLS/Form C117 (Rev. 7-93)

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

TEST	BRAC VALUE	TIME
HAIR	0.0000	12:00
SUB	0.0000	12:00
HAIR	0.0000	12:00

$$DVM = .3666$$

SUBJECT'S NAME

TIME FIRST OBSERVED

INSTRUMENT LOCATION

OPERATOR

 ADDITIONAL INFORMATION AND/OR REMARKS |

ANTOXILYZER<sup>®</sup> INSTRUMENT PRINTER CARD

© 1986 by CMI INC.



2 of 6  
He

[illegible]

SUBJECT'S NAME:

TIME FIRST OBSERVED

INSTRUMENT LOCATION

## OPERATOR

DATE	TIME	LOCATION	WIND DIRECTION	WIND SPEED	SEA STATE	WEATHER	REMARKS
							ADDITIONAL INFORMATION AND/OR REMARKS

**INTOXILYZER® INSTRUMENT PRINTER CARD**

© 1986 by CMI INC.

ANTOXILYZER<sup>®</sup> INSTRUMENT PRINTER CARD

© 1986 by CMI INC.



3 of 6  
file

# Tucson Police Department Crime Laboratory

TUCSON POLICE DEPARTMENT  
INTOXILYZER - ALCOHOL ANALYZER  
AZ MODEL 5000 SN 66-005304  
01/26/05 12:54

TEST	BRAC VALUE	TIME
AIR BLANK	.000	12:53
CAL. CHECK	.000	12:53
AIR BLANK	.000	12:54

SUBJECT'S NAME

TIME FIRST OBSERVED

INSTRUMENT LOCATION

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

INTOXILYZER® INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

CMI<sup>INC</sup>

TUCSON POLICE DEPARTMENT  
INTOXILYZER - ALCOHOL ANALYZER  
AZ MODEL 5000 SN 66-005304  
01/26/05 12:52

TEST	BRAC VALUE	TIME
AIR BLANK	.000	12:51
CAL. CHECK	.079	12:52
AIR BLANK	.000	12:52

SUBJECT'S NAME

TIME FIRST OBSERVED

INSTRUMENT LOCATION

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

INTOXILYZER® INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

CMI<sup>INC</sup>

4 of 6  
file

SN 66-000004  
E767.14  
INVALID TEST

01/29/86  
12:06

PRESSURE START TEST AT  
IMPROPER TIME

SUBJECT'S NAME

INSTRUMENT LOCATION

TIME FIRST OBSERVED

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

INTOXILYZER® INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

**CMI** INC.

EXHIBIT P

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY  
A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 5000\*

STANDARD CALIBRATION CHECK PROCEDURE

Agency Tucson Police Department Date 1/26/06 Time 1150  
Intoxilyzer Serial # 66-005004 Location DUI HQ  
QA Specialist Jennifer Coombs 46160  
(Print Name)

- (✓) 1. Pour a standard alcohol solution of known value into a clean dry simulator jar and assemble the simulator. Insure that a tight seal has been made. Standard value: 0.080 AC
- (✓) 2. Turn on the simulator and allow the temperature to reach 34°C±.2°C.
- (✓) 3. Set Intoxilyzer mode selection in the ACA mode by switching mode selection switch #9 on or selecting "C" on keyboard menu.
- (✓) 4. Attach simulator to the simulator entrance port on the Intoxilyzer.
- ✓(✓) 5. Intoxilyzer 5000 display reads "READY TO START" or "PUSH BUTTON".
- ✓(✓) 6. Push Start Test button or press enter on keyboard.
- ✓(✓) 7. Insert card in response to display.
- ✓(✓) 8. Air blank completed.
- ✓(✓) 9. Calibration check completed. Test results 0.079 AC & 0.080 AC
- ✓(✓) 10. Air blank completed.
- ✓(✓) 11. When display reads Test Complete remove evidence card. Attach the card to the completed checklist.
- (✓) 12. Return mode selection switch #9 to off position after all calibration checks are complete or type Q and enter on keyboard.

SIGNATURE

Jennifer L. Coombs

\*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD

DHS/DSLS/Form C116 (Rev. 7-93)

TUCSON POLICE DEPARTMENT  
INTOXILYZER - ALCOHOL ANALYZER  
AZ MODEL 5002 SH 66-095004  
01/26/08 12:08

TEST BRAC VALUE TIME  
AIR BLANK .000 12:07  
INVALID SAMPLE .000 12:07  
AIR BLANK .000 12:08

BEGIN NEW DEPRIVATION PERIOD

SUBJECT'S NAME

DUI HQ  
INSTRUMENT LOCATION

TIME FIRST OBSERVED

46160  
OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

Sweet Breath

INTOXILYZER® INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

CMI

TUCSON POLICE DEPARTMENT  
INTOXILYZER - ALCOHOL ANALYZER  
AZ MODEL 5002 SH 66-095004  
01/26/08 12:15

TEST BRAC VALUE TIME  
AIR BLANK .000 12:13  
\*SUSPECT TEST .131 12:15  
AIR BLANK .000 12:15

\* DEFICIENT SAMPLE - VALUE  
PRINTED HAS NOT BEEN OBTAINED.

SUBJECT'S NAME

DUI HQ  
INSTRUMENT LOCATION

TIME FIRST OBSERVED

46160  
OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

INTOXILYZER® INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

CMI

6 of 6  
file