

186
ST

EXHIBIT Q

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 5000*

BM

AGENCY Tucson Police QA SPECIALIST Irene Doran 46331
(Print Name)

INTOXILYZER SERIAL # 66-005005 LOCATION DUI HQ
DATE December 2 18 2004 TIME 0935

ST Display reads "READY TO START" or "PUSH BUTTON TO START TEST."

DIAGNOSTIC TESTS

- ST DVM Test check. Setting should be between .200 and .600. Mode selection switch S2 on, S1 and S3 off or keyboard menu selection "H". Reading is .466.
- ST Display Test check. Mode selection switch S1 on and S2 and S3 off or keyboard menu selection "V".
- ST Printer Test check. Mode selection switch S1, S2, S3 off or keyboard menu selection "P".
- ST Clock time check. Mode selection switch S10 on or keyboard menu selection "E".
- ST Date check. Mode selection switch S11 on or keyboard menu selection "E".

OPERATIONAL TESTS

- ST Alcohol free subject Test result 0.000 AC
 - ST Error recognition logic system functioning
Invalid test printed
 - ST Proper sample recognition system
Invalid sample printed ST
Deficient sample printed ST
 - N/A Completeness of sample purge with collection tube**
 - ST Calibration standard 0.080 AC
Results: 0.079AC + 0.079AC
- Instrument operating properly and accurately. YES NO

COMMENTS: _____

SIGNATURE Irene Doran

0. N/A alcohol concentration standard collected for subsequent analysis.**

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD
**THIS STEP IS ONLY REQUIRED IF THIS DEVICE IS BEING USED FOR
SAMPLE CAPTURE.

Tucson Police Department Crime Laboratory

396
ps

21

INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-005005
12/02/04 09:57

TEST BRAC VALUE TIME
AIR BLANK .000 09:57
CAL. CHECK .079 09:57
AIR BLANK .000 09:57

SUBJECT'S NAME

DUI #10

INSTRUMENT LOCATION

TIME FIRST OBSERVED

46331

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

INTOXILYZER[®] INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

ps



INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-005005
12/02/04 09:56

TEST BRAC VALUE TIME
AIR BLANK .000 09:56
CAL. CHECK .079 09:56
AIR BLANK .000 09:56

SUBJECT'S NAME

DUI #10

INSTRUMENT LOCATION

TIME FIRST OBSERVED

46331

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

INTOXILYZER[®] INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

ps



4 of 6
PS

SN 66-005005
E767.14
INVALID TEST

12/02/04
09:54

SAMPLE INTRODUCED
AT IMPROPER TIME.

SUBJECT'S NAME

TIME FIRST OBSERVED

INSTRUMENT LOCATION

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

DUT HQ

46331

INTOXILYZER[®] INSTRUMENT PRINTER CARD

© 1986 by CMI INC.



PS

EXHIBIT P

586
05

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 5000*

STANDARD CALIBRATION CHECK PROCEDURE *pm*

Agency Tucson Police Department Date 12/2/04 Time 0935
Intoxilyzer Serial# 06-005005 Location DUI HQ
QA Specialist IRVINE DORAN 46331
(Print Name)

- (OK)* 1. Pour a standard alcohol solution of known value into a clean dry simulator jar and assemble the simulator. Insure that a tight seal has been made. Standard value: 0.080 AC
- (OK)* 2. Turn on the simulator and allow the temperature to reach 34°C ± .2°C.
- (OK)* 3. Set Intoxilyzer mode selection in the ACA mode by switching mode selection switch #9 on or selecting "C" on keyboard menu.
- (OK)* 4. Attach simulator to the simulator entrance port on the Intoxilyzer.
- (OK)* 5. Intoxilyzer 5000 display reads "READY TO START" or "PUSH BUTTON".
- (OK)* 6. Push Start Test button or press enter on keyboard.
- (OK)* 7. Insert card in response to display.
- (OK)* 8. Air blank completed.
- (OK)* 9. Calibration check completed. Test results 0.079 + 0.079 AC
- (OK)* 10. Air blank completed.
- (OK)* 11. When display reads Test Complete remove evidence card. Attach the card to the completed checklist.
- (OK)* 12. Return mode selection switch #9 to off position after all calibration checks are complete or type Q and enter on keyboard.

SIGNATURE Irvine Doran

*WITH OR WITHOUT VAPOR RECIRCULATION AND WITH OR WITHOUT KEYBOARD

INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-005005
12/02/04 09:53

TEST BRAC VALUE TIME
AIR BLANK .000 09:52
INVALID SAMPLE .MDA 09:53
AIR BLANK .000 09:53

BEGIN NEW DEPRIVATION PERIOD

SUBJECT'S NAME

DUT HQ

INSTRUMENT LOCATION

TIME FIRST OBSERVED

40331

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

Sweet breath

INTOXILYZER® INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

CMI INC.

PA

INTOXILYZER - ALCOHOL ANALYZER
AZ MODEL 5000 SN 66-005005
12/02/04 09:44

TEST BRAC VALUE TIME
AIR BLANK .000 09:40
*SUBJECT TEST .096 09:43
AIR BLANK .000 09:44

* DEFICIENT SAMPLE - VALUE
PRINTED WAS HIGHEST OBTAINED.

SUBJECT'S NAME

DUT HQ

INSTRUMENT LOCATION

TIME FIRST OBSERVED

40331

OPERATOR

ADDITIONAL INFORMATION AND/OR REMARKS

INTOXILYZER® INSTRUMENT PRINTER CARD

© 1986 by CMI INC.

CMI INC.

PA

626
PA