

CERTIFICATION by CUSTODIAN of RECORDS

Tucson Police Department Crime Laboratory

270 S. Stone Avenue Tucson, AZ 85701-1917

I hereby certify that these documents are a true and correct copy of the periodic maintenance* for the Intoxilyzer 8000 Serial Number 80-555323 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.
A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.
Date: 8.27-89 Cal 1 0.096 AC Cal 2 0.094 AC
Between the dates of 7.27-09 and 8-27-09 the instrument was operating properly and accurately.
The instrument was taken out of service for repair on N/A. The instrument was returned to service on N/A.
Clave Conly
Clare L. Conley
9-1-09
ADPS QAS Permit # 33314
SUBSCRIBED AND SWORN to before me this /2+, day of, 20/19.
AMHUNANDIE OFFICIAL SEAL
Notary Public A.M. HERNÁNDEZ NOTARY PUBLIC-ARIZONA
*Includes documents specified in A.A.C. R13-10-104.A.5

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

EXHIBIT G-4

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

1.68 Chy.09

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Ware Conley AGENCY Tucson Police Department		
DATE 8-27-09 TIME 1238		
INTOXILYZER SERIAL # 80-000323 LOCATION Crime Lab		
(/) 1. Display reads "PUSH BUTTON TO START"		
DIAGNOSTIC TESTS		
<pre>(/) 1. Clock time check. (/) 2. Date check.</pre>		
OPERATIONAL TESTS		
 (V) 1. Alcohol-free subject test result 0. 600 AC. (V) 2. Error recognition logic system functioning. Not a Successfully Completed Test Sequence printed (V) 3. Proper sample recognition system. Not a Successfully Completed Test Sequence printed Deficient sample printed. (V) 4. Standard Calibration Check standard 0. 600 AC Result: 0. 694 AC 		
Instrument is operating properly and accurately. YES V NO		
COMMENTS None.		

DPS Form Exh G-4 (Rev. 05-01)

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number: 80-000323

08/27/2009

12:39:48

Standard Lot#: 811201I
Last Changed By: C. CONLEY #34517

QAS: C. CONLEY #34517

TUCSON PD

Operator: C. CONLEY #34517

TUCSON PD

Subject: TEST, SQAP

DOB: 11/11/1111

Sex: M

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	12:40:42
Diagnostic Test	Pass	12:41:09
Air Blank	0.000	12:41:36
0.100 Cal Check	0.096	12:41:55
Air Blank	0.000	12:42:24
Subject Test	0.000	12:42:42
Air Blank	ABT*	12:43:06
Air Blank	0.000	12:43:34

^{*}Sequence Aborted

388

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number: 80-000323

08/27/2009 12:43:38

Standard Lot#: 811201I
Last Changed By: C. CONLEY #34517

QAS: C. CONLEY #34517

TUCSON PD

Operator: C. CONLEY #34517
TUCSON PD

Subject: TEST, SQAP DOB: 11/11/1111

Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	12:44:11
Air Blank	0.000	12:44:40

*RFI Detect

488

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number: 80-000323

08/27/2009 12:44:4

12:44:48

Standard Lot#: 811201I

Last Changed By: C. CONLEY #34517

QAS: C. CONLEY #34517 TUCSON PD

Operator: C. CONLEY #34517

TUCSON PD

Subject: TEST, SQAP DOB: 11/11/1111

Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	12:45:20
Air Blank	0.000	12:45:49

^{*}Improper Sample

56/8 Un

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number: 80-000323

08/27/2009

12:45:53

Standard Lot#: 811201I Last Changed By: C. CONLEY	#34517
QAS: C. CONLEY TUCSON PD	#34517
Operator: C. CONLEY TUCSON PD	#34517

Subject: TEST, SQAP

DOB: 11/11/1111

Sex: M

Weight: 111

15 Minute Deprivation Period? Yes

	/240=	·
Test	g/210L	Time
Air Blank	0.000	12:46:44
Diagnostic Test	Pass	12:47:11
Air Blank	0.000	12:47:38
0.100 Cal Check	0.096	12:47:57
Air Blank	0.000	12:48:26
Subject Test	0.124*	12:51:35
Air Blank	0.000	12:52:08
Air Blank	0.000	12:52:35
Subject Test	XXX**	12:52:56
Air Blank	0.000	12:53:30

^{*}Deficient Sample

Begin new deprivation period

^{**}Invalid Sample

INTOXILYZER 8000

Location: TUCSON POLICE DEPT Serial Number: 80-000323

08/27/2009

12:53:51

Standard Lot#: 811201I

Last Changed By: C. CONLEY #34517

QAS: C. CONLEY #34517

TUCSON PD

g/210L Time Test 12:54:19 0.000 Air Blank 0.094 Cal Check 12:54:37 Air Blank 0.000 12:55:07

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA Special:	ist Mare Conley	AGENCY Tucson Police Department
DATE 8.3	•	TIME 1238
Intoxilyze	r Serial # 80-000323	Location Crine Lab
(✓) 1.	Ensure that gas tank is at standard alcohol concentrations	tached to instrument and contains a tion solution 0. AC.
	Pour a standard alcohol into a clean dry simulator	and assemble the simulator. Ensure en made. Turn on the simulator and
(/) 2.	Intoxilyzer 8000 displa	ay reads "PUSH BUTTON TO START"
(_V) 3.		ting Menu". Select "D" for dry for wet control test. After ENTER.
(🗸) 4.	Air blank completed.	
(🗸) 5.	Calibration check compl	leted. Test results 0. 09 AC.
(v) 6.	Air blank completed.	
(🗸) 7.	Remove printed record. checklist.	Attach the record to the completed
SIGNATURE	Clare Conley	

DPS Form Exh G-2 (Rev. 05-01)

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number: 80-000323

08/27/2009

12:38:28

Standard Lot#: 811201I

Last Changed By: C. CONLEY #34517

QAS: C. CONLEY #34517

TUCSON PD

Test	g/210L	Time
Air Blank	0.000	12:38:55
Cal Check	0.096	12:39:14
Air Blank	0 000	12-39-43