



CERTIFICATION by CUSTODIAN of RECORDS
 Tucson Police Department Crime Laboratory
 270 S. Stone Avenue Tucson, AZ 85701-1917

*TSy
12-17-09*

I hereby certify that these documents are a true and correct copy of the periodic maintenance* for the Intoxilyzer 8000 Serial Number 80-000323 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.

A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.

Date: 12/17/09 Cal 1 0.089 AC Cal 2 0.098 AC

Between the dates of 4/19/09 and 12/17/09 the instrument was operating properly and accurately.

The instrument was taken out of service for repair on N/A.
 The instrument was returned to service on N/A.

Cheryl K. Straub

Cheryl K. Straub

12/17/09

ADPS QAS Permit # 33408

SUBSCRIBED AND SWORN to before me this 17th, day of December, 2009.

A.M. Hernandez

Notary Public



*Includes documents specified in A.A.C. R13-10-104.A.5

187
CS

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD QUALITY ASSURANCE PROCEDURE

TSY
12-17-09

QA SPECIALIST Cheryl STRAUB AGENCY Tucson Police Department
DATE 12/17/09 TIME 10:28
INTOXILYZER SERIAL # 80-000323 LOCATION CRIME Lab

() 1. Display reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

() 1. Clock time check.
() 2. Date check.

OPERATIONAL TESTS

() 1. Alcohol-free subject test result 0. 000 AC.
() 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
() 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
() 4. Standard Calibration Check standard 0. 100 AC
Result: 0. 098 AC

Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE Cheryl Straub

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

217
CS

TS4

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number: 80-000323
12/17/2005 10:38:18

Standard Lot#: 8112011
Last Changed By: C. CONLEY #34517

OAG: C. STRAUB #50772
TUCSON PD

Operator: C. STRAUB #50772
TUCSON PD

Subject: SOAP, SOAP
DOB: 11/11/1111
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	10:39:09
Diagnostic Test	Pass	10:39:38
Air Blank	0.000	10:40:05
0.100 Cal Check	0.099	10:40:24
Air Blank	0.000	10:40:53
Subject Test	0.000	10:41:18
Air Blank	0.000	10:41:45
Five Minute Wait	ABF*	
Air Blank	0.000	10:42:20

*Sequence Aborted

Not a Successfully
Completed
Test Sequence

377
CS
154

INTOXICOR 8000
Location: TUCSON POLICE DEPT
Serial Number: 80-010323
12/17/2009 10:42:35

Standard Lot#: 8112011
Last Changed By: C. CONLEY #34517

Off: C. STRUB #50772
TUCSON PD

Operator: C. STRUB #50772
TUCSON PD

Subject: SOAP SOAP
DOB: 11/11/1111
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RF1*	10:43:19
Air Blank	IP5**	10:43:29

*RF1 Detect
**Improper Sample

Not a Successfully
Completed
Test Sequence

497
CS

TS4

INSTRUMENT #000
Location: TUCSON POLICE DEPT
Serial Number: 00-100323
12/17/2009 10:44:25

Standard Lots: 811201
Last Changed By: C. COMLEY #34517

Off: C. STRAUB #50772
TUCSON PD

Operator: C. STRAUB #50772
TUCSON PD

Subject: SOAP_SOW
DOB: 11/11/1111
Sex: F Weight: 111

15 Minute Deprivation Period? yes

Test	g/210	Time
Air Blank	0.000	10:45:15
Diagnostic Test	Pass	10:45:45
Air Blank	0.000	10:46:12
0.100 Cal Check	1.099	10:46:31
Air Blank	0.000	10:47:00
Subject Test	0.344*	10:50:09
Air Blank	0.000	10:50:41
Air Blank	0.000	10:51:09
Subject Test	0.00**	10:51:38
Air Blank	0.000	10:52:11

*Deficient Sample
**Invalid Sample
Begin new deprivation period

Not a Successfully
Completed
Test Sequence

517
CS

10A

INTOXLYZER 8100
Location: TUCSON POLICE DEPT
Serial Number: 80-803323
12/17/2009 10:50:05

Standard Load: 8112811
Last Changed By: C. COMBY #34517

OAS: C. STRAUB #50772
TUCSON PD

Test	g/210L	Time
Air Blank	0.000	10:50:02
Cal Check	0.000	10:50:01
Air Blank	0.000	10:50:00

Chap. Straub
Operator's Signature
|

EXHIBIT G-2

677
CS

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

BY
12-17-09

QA Specialist Cheryl STRAUB AGENCY Tucson Police Department
DATE 12/17/09 TIME 10:28
Intoxilyzer Serial # 80-000323 Location CRIME LAB

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0. 100 AC.
OR
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 099 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Cheryl Straub

DPS Form Exh G-2 (Rev. 05-01)

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

797
CS

154

INTOXILYZER 8011
Location: TUCSON POLICE DEPT
Serial Number: 90-100323
12/17/2009 10:29:46

Standard Lock: 6112811
Last Changed By: C. CONLEY #54517

Off: C. STRAUB #50772
TUCSON PD

Test	g/210L	Time
Air Blank	0.000	10:29:13
Cal Check	0.099	10:29:32
Air Blank	0.000	10:30:31


Operator's Signature