

## CERTIFICATION by CUSTODIAN of RECORDS

Tucson Police Department Crime Laboratory 270 S. Stone Avenue Tucson, AZ 85701-1917

| I hereby certify that these documents are a true and correct copy of the periodic maintenance* for the Intoxilyzer 8000 Serial Number  80-000323 maintained by the Tucson Police Department  Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program. |
|--|
| A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.   |
| Date: 05/06/10 Cal 1 0.098 AC Cal 2 0.098 AC   |
| Between the dates of 04/08/10 and 05/06/10 the instrument was operating properly and accurately.   |
| The instrument was taken out of service for repair on $N/A$ .  The instrument was returned to service on $N/A$ .   |
| Chyld. Stan  |
| Cheryl K. Straub   |
| 05/06/10   |
| ADPS QAS Permit # 33408  |
| SUBSCRIBED AND SWORN to before me this 645, day of   |
| - May ,20/0.  ANH ornanded   |
| Notary Public  Notary Public  Notary Public  Notary Public  Notary Public-ARIZON  PIMA COUNTY  |

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the

information contained in computer storage devices of the Tucson Police Department.

My Comm. Exp. Mar. 17, 2012

\*Includes documents specified in A.A.C. R13-10-104.A.5

## THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

754 5.6.2010

STANDARD QUALITY ASSURANCE PROCEDURE

| 1 / 2   | <u> </u> |
|---|----------|
| QA SPECIALIST Charyl Straub AGENCY Tucson Police Department DATE 05/06/10 TIME 10:29  |          |
| INTOXILYZER SERIAL # 80-000323 LOCATION CRIME LAb   |          |
| 1. Display reads "PUSH BUTTON TO START"   |          |
| DIAGNOSTIC TESTS  |          |
| ( ) 1. Clock time check. ( ) 2. Date check.   |          |
| OPERATIONAL TESTS   |          |
| 1. Alcohol-free subject test result 0. 000 AC.  2. Error recognition logic system functioning.  Not a Successfully Completed Test Sequence printed  3. Proper sample recognition system.  Not a Successfully Completed Test Sequence printed  Deficient sample printed.  4. Standard Calibration Check standard 0. 100 AC  Result: 0. 098 AC  Instrument is operating properly and accurately. YES NO |          |
| COMMENTS  |          |
|   |          |
|   |          |
|   |          |
|   |          |

DPS Form Exh G-4 (Rev. 05-01)

Pursuant to section 28-1327. Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

INTOXILYZER 8000 Location: TUCSON POLICE DEPT

Serial Number: 80-000323 05/06/2010

10:41:03

Standard Lot#: 642845

Last Changed By: C. STRAUB #50772

QAS: C. STRAUB #50772

TUCSON PD

Operator: C. STRAUB #50772

TUCSON PD

Subject: SOAP, SOAP DOB: 11/11/1111

Sex: F Weight: 111

15 Minute Deprivation Period? Yes

| Test   | g/210L   | Time   |
|--|--|--|
| Air Blank Diagnostic Test Air Blank 0.100 Cal Check Air Blank Subject Test Air Blank Air Blank | 0.000<br>Pass<br>0.000<br>0.098<br>0.000<br>0.000<br>ABT*<br>0.000 | 10:41:54<br>10:42:23<br>10:42:50<br>10:43:09<br>10:43:38<br>10:44:11<br>10:44:18 |

\*Sequence Aborted

Not a Successfully Completed Test Sequence 184

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number: 86-888323

Standard Lot#: 642845 `
Last Changed By: C. STRAUB #50772

QAS: C. STRAUB #50772

TUCSON PD

Operator: C. STRAUB #50772

TUCSON PD

Subject: SQAP,SQAP DOB: 11/11/1111

Weight: 111

15 Minute Deprivation Period? Yes

| Test                   | g/210L        | Time                 |
|------------------------|---------------|----------------------|
| Air Blank<br>Air Blank | RF[*<br>[PS** | 10:46:14<br>10:46:22 |
| *RFI Detect            | nie           |                      |

Not a Successfully Completed Test Sequence

INTOXILYZER 8000 Location: TUCSON POLICE DEPT

Serial Number: 80-000323

10:47:09

Standard Lot#: 642845

Last Changed By: C. STRAUB #50772

QAS: C. STRAUB #50772

TUCSON PD

Operator: C. STRAUB #50772

TUCSON PD

Subject: SQAP,SQAP DOB: 11/11/1111 Weight: 111

Sex: F

15 Minute Deprivation Period? Yes

| Test  | g/210L   | Time  |
|---|--|---|
| Air Blank Diagnostic Test Air Blank 0.100 Cal Check Air Blank Subject Test Air Blank Air Blank Subject Test Air Blank | 0.000<br>Pass<br>0.000<br>0.098<br>0.000<br>0.340*<br>0.000<br>0.000<br>XXX**<br>0.000 | 10: 47: 59 10: 48: 28 10: 48: 55 10: 49: 14 10: 49: 43 10: 52: 52 10: 53: 25 10: 53: 52 10: 54: 66 10: 54: 40 |
|   |  |   |

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: TUCSON POLICE DEPT Serial Number: 80-003323

05/06/2010

10:55:10

Standard Lot#: 642845

Last Changed By: C. STRAUB #50772

QAS: C. STRAUB #50772

TUCSON PD

 Test
 g/210L
 Time

 Air Blank
 0.000
 10:55:38

 Cal Check
 0.098
 10:55:56

 Air Blank
 0.000
 10:56:26

Oly Ita
Operator's stgrature

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

B4 2010

STANDARD CALIBRATION CHECK PROCEDURE

| QA         | Spe        | ciali | ist ChERUL STRAUB   | _ AGENCY              | Tucson Poli          | ce Department    |
|------------|------------|-------|---|-----------------------|----------------------|------------------|
| DA'        | TE _       | 05    | 5/06/10   | TIME                  | 10:29                |                  |
| In         | toxi       | lyze: | ist <u>Cheryl StrauB</u><br>5/06/10<br>r Serial # 80-000323   | _ Locatio             | n CRIME              | LAb              |
|            |            |       |   |                       |                      |                  |
| ( <b>v</b> | )          | 1.    | Ensure that gas tank is a standard alcohol concentr   | 301011 3010           | 101011 0. <u>/00</u> |                  |
|            |            |       | Pour a standard alcoholinto a clean dry simulato that a tight seal has be allow temperature to reac | een made.<br>h 34°C + | Turn on the 0.2°C    | simulator and    |
| ( 1        | <b>/</b> ) | 2.    | Intoxilyzer 8000 displ  | ay reads              | "PUSH BUTT           | ON TO START"     |
|            |            |       | Intoxilyzer 8000 displ<br>Go to the "Control Te<br>control test or "W"<br>selection is made pres    | S ENTER.              | CONCEON              | Jest, Arter      |
| ( •        | )          | 4.    | Air blank completed.  |                       |                      |                  |
| ( i        | <b>/</b> } | 5.    | Air blank completed.  Calibration check comp  | oleted. T             | est results          | 0. <u>098</u> AC |
| ( L        | <b>/</b> ) | 6.    | Air blank completed.  |                       |                      |                  |
| ( <b>v</b> | <u></u>    | 7.    | Remove printed record checklist.  | . Attach              | the record to        | the completed    |
|            |            |       | Clay Star   |                       |                      |                  |
| SI         | GNAI       | 'URE  | My Jan  |                       |                      |                  |

DPS Form Exh G-2 (Rev. 05-01)

INTOXILYZER 8000 Location: TUCSON POLICE DEPT Serial Number: 80-000323 05/06/2010

10:29:56

Standard Lot#: 642845

Last Changed By: C. STRAUB #50772

OAS: C. STRAUB #50772

TUCSON PD

| Test                                | g/210L                  | Time                             |  |
|-------------------------------------|-------------------------|----------------------------------|--|
| Air Blank<br>Cal Check<br>Air Blank | 0.000<br>0.098<br>0.000 | 10:30:24<br>10:30:42<br>10:31:12 |  |