

#### CERTIFICATION by CUSTODIAN of RECORDS

Tucson Police Department Crime Laboratory 270 S. Stone Avenue Tucson, AZ 85701-1917

6/7/10 I hereby certify that these documents are a true and correct copy of the periodic maintenance\* for the Intoxilyzer 8000 Serial Number MO-000323 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program. A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC. Date: 6.3.2010 Cal 1 0.098 AC Cal 2 0.096 AC Between the dates of 5.6.2010 and 6.3.2010the instrument was operating properly and accurately. The instrument was taken out of service for repair on The instrument was returned to service on Seen (allego Terry S. Gallegos 6.7-2010 ADPS QAS Permit # 33169 SUBSCRIBED AND SWORN to before me this **7** , day of m//wnanded Notary Public PIMA COUNTY \*Includes documents specified in A.A.C. R13-10-104.A.5 My Comm. Exp. Mar. 17, 2012

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

## THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Terry Gallegos 28847	AGENCY To	icson Police D	epartment
DATE 6.3.2010	TIME	1531	
DATE 6.3.2010 INTOXILYZER SERIAL # 80-000323	LOCATION	CrimeLab	
1. Display reads "PUSH BUTTON TO		en e	
DIAGNOSTIC TESTS			
( ) 1. Clock time check. 2. Date check.			e de la companya de l
OPERATIONAL TESTS			
1. Alcohol-free subject test rest 2. Error recognition logic system Not a Successfully Completed 5 3. Proper sample recognition system Not a Successfully Completed 5 Deficient sample printed. 4. Standard Calibration Check standard Calibration Check standard control of the complete of t	Pest Sequenc tem. Test Seguenc	ce printed	
Instrument is operating properly and acc	curately. YE	ES NO	
COMMENTS			
			·
SIGNATURE Sunfallego			* ·

DPS Form Exh G-4 (Rev. 05-01)

Location: TUCSON POLICE DEPT

Serial Number: 80-000323

06/03/2010

15:32:31

Standard Lot#: 642845

Last Changed By: C. STRAUB #50772

QAS: T. GALLEGOS #28847

TUCSON PD CRIME LAB

Operator: T. GALLEGOS #28847

TUCSON PD CRIME LAB

Subject: TEST, SQAP DOB: 11/11/1111

Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L		Time	
Air Blank	0.000		15:33:21	•
Diagnostic Test	Pas <b>s</b>		15:33:48	
Air Blank	0.000		15:34:15	
0.100 Cal Check	0.097		15:34:34	
Air Blank	0.000		15:35:03	
Subject Test	0.000	*	15:35:18	
Air Blank	ABT*		15:35:26	
Air Blank	0.000		15:35:55	

<sup>\*</sup>Sequence Aborted

### Not a Successfully Completed Test Sequence

Pursuant to section 28-1327. Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

2057

Location: TUCSON POLICE DEPT

Serial Number: 80-000323

06/03/2010

15:35:59

Standard Lot#: 642845

Last Changed By: C. STRAUB

#50772

QAS: T. GALLEGOS

TUCSON PD CRIME LAB

#28847

Operator: T. GALLEGOS

#28847

TUCSON PD CRIME LAB

Subject: TEST, SQAP

DOB: 11/11/1111

Sex: F

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	15:36:42
Air Blank	IPS**	15:36:51

<sup>\*</sup>RFI Detect

### Not a Successfully Completed Test Sequence

<sup>\*\*</sup>Improper Sample

Location: TUCSON POLICE DEPT

Serial Number: 80-000323

06/03/2010

15:37:29

Standard Lot#: 642845

Last Changed By: C. STRAUB

#50772

QAS: T. GALLEGOS

#28847

TUCSON PD CRIME LAB

Operator: T. GALLEGOS

#28847

TUCSON PD CRIME LAB

Subject: TEST, SQAP DOB: 11/11/1111

Sex: F

Weight: 11

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	15:38:19
Diagnosti <b>c</b> Te <b>st</b>	Pas <b>s</b>	15:38:46
Air Blank	0.000	15:39:14
0.100 Cal Check	0.097	15:39:32
Air Blank	0.000	15:40:02
Subject Test	0.178*	15:43:10
Air Blank	0.000	15:43:42
Air Blank	0.000	15:44:10
Subject Test	XXX**	15:44:27
Air Blank	0.000	15:44:59

<sup>\*</sup>Deficient Sample

Begin new deprivation period

#### Not a Successfully Completed Test Sequence

<sup>\*\*</sup>Invalid Sample

Location: TUCSON POLICE DEPT

Serial Number: 80-000323

06/03/2010 15:45:17

Standard Lot#: 642845

Last Changed By: C. STRAUB #50772

·----

QAS: T. GALLEGOS #28847

TUCSON PD CRIME LAB

Test	ì	g/210L	Time
Air Blank		0.000 0.096	15:45:44
Air Blank		0.000	15:46:32

Sun Callege Operators Signature

# THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.O. R13-10-104(A,

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

6057

5C

STANDARD CALIBRATION CHECK PROCEDURE

DATE .	6	ist <u>Terry Gallegos 2884.</u> · 3·2010	TIME	1531	
Intoxi	llyze	-3.2010 r Serial # 80-000323	Location	Crime Lab	
	÷.	Ensure that gas tank is a standard alcohol concentra OR	tached to tion solut:	instrument and coion 0.100 AC.	ontains a
		Pour a standard alcohol into a clean dry simulator that a tight seal has be allow temperature to reach	r and assem en made. T	able the simulator Turn on the simul	AC, . Ensure ator and
	2.	Intoxilyzer 8000 displa	ay reads '	PUSH BUTTON TO	START"
	3.	Go to the "Control Test control test or "W" selection is made press	ting Men	u". Select "D"	for drv
	4.	Air blank completed.			
<b>O</b> )	5.	Calibration check compl	eted. Tes	st results 0. <u>0</u> 9	8_AC.
<b>/</b> )	6.	Air blank completed.			
<b>(</b> )	7.	Remove printed record. checklist.	Attach th	e record to the o	ompleted:
SIGNAT	TURE	Sempeallegn			

DPS Form Exh G-2 (Rev. 05-01)

Location: TUCSON POLICE DEPT

Serial Number: 80-000323

06/03/2010

15:31:00

Standard Lot#: 642845

Last Changed By: C. STRAUB #50772

QAS: T. GALLEGOS #28847

TUCSON PD CRIME LAB

Test	g/210L	Time
Air Blank	0.000	15:31:28
Cal Check	<b>0.</b> 098	15:31:46
Air Blank	0.000	15:32:15

Operator's Signature