



CERTIFICATION by CUSTODIAN of RECORDS

Tucson Police Department Crime Laboratory
270 S. Stone Avenue Tucson, AZ 85701-1917

I hereby certify that these documents are a true and correct copy of the periodic maintenance* for the Intoxilyzer 8000 Serial Number 80-006331 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.

154
10-2-09

A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.

Date: 09/17/09 Cal 1 0.098 AC Cal 2 0.096 AC

Between the dates of 08/20/09 and 09/17/09 the instrument was operating properly and accurately.

The instrument was taken out of service for repair on N/A.
The instrument was returned to service on N/A.

Cheryl K. Straub

Cheryl K. Straub

09/29/09

ADPS QAS Permit # 33408

SUBSCRIBED AND SWORN to before me this 29th, day of September, 20 09.

A.M. Hernández

Notary Public



*Includes documents specified in A.A.C. R13-10-104.A.5

EXHIBIT G-4

197
cs

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD QUALITY ASSURANCE PROCEDURE

154
10.2.09

QA SPECIALIST Cheryl Straub AGENCY Tucson Police Department
DATE 09/17/09 TIME 10:51
INTOXILYZER SERIAL # 80-000331 LOCATION DUI/HQ

() 1. Display reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- () 1. Clock time check.
- () 2. Date check.

OPERATIONAL TESTS

- () 1. Alcohol-free subject test result 0. 000 AC.
- () 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
- () 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
- () 4. Standard Calibration Check standard 0. 100 AC
Result: 0. 094 AC

Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE

Cheryl Straub

287
CS

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number : 90-000331
09/17/2009

10:54:45

SM

Standard Lot#: 811201I
Last Changed By: C. CONLEY #34517

QAS: C. STRAUB #50772
TUCSON PD

Operator: C. STRAUB #50772
TUCSON PD

Subject: SQAP, SQAP
DOB: 11/11/1111
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	10:55:43
Diagnostic Test	Pass	10:56:10
Air Blank	0.000	10:56:38
0.100 Cal Check	0.098	10:56:56
Air Blank	0.000	10:57:26
Subject Test	0.000	10:57:44
Air Blank	ABT*	10:57:51
Air Blank	0.000	10:58:19

*Sequence Aborted

Not a Successfully Completed Test Sequence

387
CS

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number : 80-000331
09/17/2009 10:59:21

134

Standard Lot#: 811201I
Last Changed By: C. CONLEY #34517

QAS: C. STRAUB #50772
TUCSON PD

Operator: C. STRAUB #50772
TUCSON PD

Subject: SQAP, SQAP
DOB: 11/11/1111
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	11:00:53
Air Blank	IPS**	11:01:07

*RFI Detect
**Improper Sample

Not a Successfully Completed Test Sequence

477
CS

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number : 80-000331
09/17/2009 11:02:22

Standard Lot#: 811201I
Last Changed By: C. CONLEY #34517

B4

QAS: C. STRAUB #50772
TUCSON PD

Operator: C. STRAUB #50772
TUCSON PD

Subject: SQAP, SQAP
DOB: 11/11/1111
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	11:03:20
Diagnostic Test	Pass	11:03:47
Air Blank	0.000	11:04:14
0.100 Cal Check	0.097	11:04:33
Air Blank	0.000	11:05:03
Subject Test	0.260*	11:08:11
Air Blank	0.000	11:08:44
Air Blank	0.000	11:09:11
Subject Test	XXX**	11:09:32
Air Blank	0.000	11:10:06

*Deficient Sample
**Invalid Sample
Begin new deprivation period

Not a Successfully Completed Test Sequence

Pursuant to section 26-1527, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

587
CS

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number: 80-000331
09/17/2009 11:10:44

Standard Lot#: 8112011
Last Changed By: C. CONLEY #34517

TS4

QAS: C. STRAUB #50772
TUCSON PD

Test	g/210L	Time
Air Blank	0.000	11:11:11
Cal Check	0.096	11:11:30
Air Blank	0.000	11:12:00



Operator's Signature

EXHIBIT G-2

677
cg

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

TS4
10-2-09

QA Specialist Cheryl Straub AGENCY Tucson Police Department
DATE 09/17/09 TIME 10:51
Intoxilyzer Serial # 80-000331 Location DUI/HQ

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0. 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 098 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Cheryl Straub

DPS Form Exh G-2 (Rev. 05-01)

787
es

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number: 80-000331
09/17/2009 10:51:53

154

Standard Lot#: 811201I
Last Changed By: C. CONLEY #34517

QAS: C. STRAUB #50772
TUCSON PD

Test	g/210L	Time
Air Blank	0.000	10:52:20
Cal Check	0.098	10:52:39
Air Blank	0.000	10:53:09



Operator's Signature