

CERTIFICATION by CUSTODIAN of RECORDS

Tucson Police Department Crime Laboratory
270 S. Stone Avenue Tucson, AZ 85701-1917

I hereby certify that these documents are a true and correct copy of the periodic maintenance* for the Intoxilyzer 8000 Serial Number  80 - 000 331  maintained by the Tucson Police Department
Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.
A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.
Date: 5-26-10 Cal 1 0 100 AC Cal 2 0.099 AC
Between the dates of 4.29-10 and 5-26-10 the instrument was operating properly and accurately.
The instrument was taken out of service for repair on N/A.  The instrument was returned to service on N/A.
<u>Clare Cooling</u>
Clare L. Conley
5.26-10
ADPS QAS Permit # 33314
SUBSCRIBED AND SWORN to before me this 26th, day of
$\frac{1/1/Ay}{2010}$
Notary Public OFFICIAL SEAL
*Includes documents specified in A.A.C. R13-10-104.A.5  *Includes documents specified in A.A.C. R13-10-104.A.5  *Includes documents specified in A.A.C. R13-10-104.A.5  *Includes documents specified in A.A.C. R13-10-104.A.5

# THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

54.2010 6.1.2010

QA SPECIALIST <u>Clare Con ley</u>	AGENCY Tucson Police Department
DATE 5.26-10	TIME 0826
INTOXILYZER SERIAL # 80-000331	LOCATION ODW DW HQ
( ✔ ) 1. Display reads "PUSH BUTTON TO	START"
DIAGNOSTIC TESTS	
( ) 1. Clock time check. ( ) 2. Date check.  OPERATIONAL TESTS	
<ul> <li>( ) 1. Alcohol-free subject test res</li> <li>( ) 2. Error recognition logic systemot a Successfully Completed</li> <li>( ) 3. Proper sample recognition sysmot a Successfully Completed Deficient sample printed.</li> <li>( ) 4. Standard Calibration Check standard: 0.099 AC</li> </ul>	em functioning. Test Sequence printed stem. Test Sequence printed
Instrument is operating properly and ac	curately. YES / NO
COMMENTS Nme.	
SIGNATURE Une Coly	

DPS Form Exh G-4 (Rev. 05-01)

Location: TUCSON POLICE DEPT

Serial Number: 80-000331

05/26/2010

08:28:07

Standard Lot#: 654682

Last Changed By: C. CONLEY #34517

QAS: C. CONLEY

#34517

TUCSON PD

Operator: C. CONLEY

#34517

TUCSON PD

Subject: TEST, SOAP DOB: 11/11/1111

Sex: M

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:29:01
Diagnosti <b>c Tes</b> t	Pass	08:29:28
Air Blank	0.000	08:29:56
0.100 Cal Check	0.099	08:30:14
Air Blank	0.000	08:30:44
Subject Test	0.000	08:30:59
Air Blank	ABT*	08:31:06
Air Blank	0.000	08:31:34

<sup>\*</sup>Sequence Aborted

## Not a Successfully Completed Test Sequence

Location: TUCSON POLICE DEPT

Serial Number: 80-000331

05/26/2010 Standard Lot#: 654682

Last Changed By: C. CONLEY #34517

08:31:38

QAS: C. CONLEY #34517

TUCSON PD

Operator: C. CONLEY #34517

TUCSON PD

Subject: TEST, SQAP

DOB: 11/11/1111
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

		•
Test	g/210L	Time
Air Blank	0.000	08:32:26
Diagnostic <b>Tes</b> t	Pass	08:32:54
Air Blank	0.000	08:33:21
0.100 Cal Check	0.099	08:33:40
Air Blank	RFI*	08:33:49
Air Blank	IPS**	08:34:00

<sup>\*</sup>RFI Detect

## Not a Successfully Completed Test Sequence

347

<sup>\*\*</sup>Improper Sample

Location: TUCSON POLICE DEPT

Serial Number: 80-000331

05/26/2010

08:34:27

Standard Lot#: 654682

Last Changed By: C. CONLEY #34517

QAS: C. CONLEY

#34517

TUCSON PD

Operator: C. CONLEY

#3451**7** 

TUCSON PD

Subject: TEST, SQAP

DOB: 11/11/1111

Sex: M

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:35:12
Diagnostic Test	Pass	08:35:40
Air Blank	0.000	08:36:07
0.100 Cal Check	0.099	08:36:26
Air Blank	0.000	08:36:55
Subject Test	0.198*	08:40:03
Air Blank	0.000	08:40:36
Air Blank	0.000	08:41:03
Subject Test	XXX**	08:41:49
Air Blank	0.000	08:42:24

<sup>\*</sup>Deficient Sample

Begin new deprivation period

### Not a Successfully Completed Test Sequence

lu Un

34

<sup>\*\*</sup>Invalid Sample

Location: TUCSON POLICE DEPT

Serial Number: 80-000331

05/26/2010

08:42:39

Standard Lot#: 654682

Last Changed By: C. CONLEY #34517

QAS: C. CONLEY #34517

TUCSON PD

\_\_\_\_\_

Test	•	g/210L	Time
Air	Blank	0.000	08:43:06
Cal	Check	0.099	08:43:25
Air	Blank	0.000	08:43:55

Operator's Signature

# THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

667

ARIZONA DEPARTMENT OF PUBLIC SAFETY

al C

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

(10,7010

STANDARD CALIBRATION CHECK PROCEDURE

QA Special	ist Clare Con ley AGENCY Tucson Police Department
DATE 5	
Intoxilyze	r Serial # 80 000331 Location ODW DWHQ
	en e
( ) 1.	Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0. 100 AC.
	Pour a standard alcohol concentration solution AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
( ~ ) 2.	Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
( <sub>V</sub> ) 3.	Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
( / ) 4.	Air blank completed.
( <b>/</b> ) . 5.	Calibration check completed. Test results 0. 100 AC.
(V) 6.	Air blank completed.
( V ) 7.	Remove printed record. Attach the record to the completed checklist.
SIGNATURE	Clare Conley
DDG ===================================	

Location: TUCSON POLICE DEPT

Serial Number: 80-000331

05/26/2010

08:26:02

Standard Lot#: 654682

Last Changed By: C. CONLEY #34517

QAS: C. CONLEY #34517

grib. c. conder "3

TUCSON PD

\_\_\_\_\_

Test	·	g/210L	Time
Air	Blank	0.000	08:26:29
Cal	Check	0.100	08:26:48
Air	Blank	0.000	08:27:17

Operator's Signature

756 lu

pol