



CERTIFICATION by CUSTODIAN of RECORDS  
Tucson Police Department Crime Laboratory

TS4  
8-8-07

I hereby certify that these documents are a true and correct copy of the periodic maintenance\* for the Intoxilyzer 8000 Serial Number 80-000332 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.

A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.

Date: 8-2-07 Cal 1 0.098 AC Cal 2 0.098 AC

Between the dates of 7-6-07 and 8-2-07 the instrument was operating properly and accurately.

The instrument was taken out of service for repair on N/A.  
The instrument was returned to service on N/A.

Clare Conley


Clare L. Conley

8-3-07

ADPS QAS Permit # 33314

SUBSCRIBED AND SWORN to before me this 3, day of

August, 2007.



Kathy A Kinney

Notary Public

\*Includes documents specified in A.A.C. R13-10-104.A.5

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY  
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000B4  
8-8-07

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Clare Conley AGENCY Tucson Police DepartmentDATE 8-2-07 TIME 08:21INTOXILYZER SERIAL # 80-000332 LOCATION ODE

(✓) 1. Display reads "PUSH BUTTON TO START"

## DIAGNOSTIC TESTS

- (✓) 1. Clock time check.
- 
- (✓) 2. Date check.

## OPERATIONAL TESTS

- (✓) 1. Alcohol-free subject test result 0.
- 000
- AC.
- 
- (✓) 2. Error recognition logic system functioning.
- 
- Not a Successfully Completed Test Sequence printed
- 
- (✓) 3. Proper sample recognition system.
- 
- Not a Successfully Completed Test Sequence printed
- 
- Deficient sample printed.
- 
- (✓) 4. Standard Calibration Check standard 0.
- 100
- AC
- 
- Result: 0.
- 098
- AC

Instrument is operating properly and accurately. YES ☒ NO ☐COMMENTS None.  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_SIGNATURE Clare Conley

268 cu

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number : 80-000332

08/02/2007

08:29:53

B4

Standard Lot#: 630601I

Last Changed By: J. COOMBS

#46160

QAS: C. CONLEY

#34517

TUCSON PD

Operator: C. CONLEY

#34517

TUCSON PD

Subject: TEST, SQAP

DOB: 11/11/1111

Sex: M

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:30:43
Diagnostic Test	Pass	08:31:11
Air Blank	0.000	08:31:38
0.100 Cal Check	0.099	08:31:57
Air Blank	0.000	08:32:26
Subject Test	0.103*	08:35:35
Air Blank	0.000	08:36:08
Air Blank	0.000	08:36:35
Subject Test	XXX**	08:37:00
Air Blank	0.000	08:37:34

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

### Not a Successfully Completed Test Sequence

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

3 of 8 cu

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number : 80-000332

08/02/2007

08:28:47

TSY

Standard Lot#: 630601I

Last Changed By: J. COOMBS

#46160

QAS: C. CONLEY

#34517

TUCSON PD

Operator: C. CONLEY

#34517

TUCSON PD

Subject: TEST, SQAP

DOB: 11/11/1111

Sex: M

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	08:29:20
Air Blank	0.000	08:29:49

\*Improper Sample

### Not a Successfully Completed Test Sequence

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

4 of 8 cu

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number : 80-000332

08/02/2007

08:27:26

by

Standard Lot#: 630601I

Last Changed By: J. COOMBS

#46160

QAS: C. CONLEY

#34517

TUCSON PD

Operator: C. CONLEY

#34517

TUCSON PD

Subject: TEST, SQAP

DOB: 11/11/1111

Sex: M

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	08:28:15
Air Blank	0.000	08:28:43

\*RFI Detect

### Not a Successfully Completed Test Sequence

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

5 of 8 cu

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number : 80-000332

08/02/2007

08:23:53

134

Standard Lot#: 630601I

Last Changed By: J. COOMBS

#46160

QAS: C. CONLEY

#34517

TUCSON PD

Operator: C. CONLEY

#34517

TUCSON PD

Subject: TEST, SQAP

DOB: 11/11/1111

Sex: M

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:24:46
Diagnostic Test	Pass	08:25:13
Air Blank	0.000	08:25:41
0.100 Cal Check	0.098	08:25:59
Air Blank	0.000	08:26:29
Subject Test	0.000	08:26:45
Air Blank	ABT*	08:26:54
Air Blank	0.000	08:27:22

\*Sequence Aborted

### Not a Successfully Completed Test Sequence

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

6 of 8 cu

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number: 80-000332

08/02/2007

08:38:01

PH

-----  
Standard Lot#: 6306011

Last Changed By: J. COOMBS #46160  
-----

QAS: C. CONLEY #34517

TUCSON PD  
-----

Test	g/210L	Time
Air Blank	0.000	08:38:28
Cal Check	0.098	08:38:47
Air Blank	0.000	08:39:17

Clare Conley

Operator's Signature

7 of 8u

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY  
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

PH

## STANDARD CALIBRATION CHECK PROCEDURE

QA Specialist Clare Conley AGENCY Tucson Police Department  
DATE 8-2-07 TIME 08:21  
Intoxilyzer Serial # 80-000332 Location ODE

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0. 100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 098 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Clare Conley

DPS Form Exh G-2 (Rev. 05-01)



8 of 8 cu

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number: 80-000332

08/02/2007

08:21:40

SH

-----  
Standard Lot#: 630601I

Last Changed By: J. COOMBS #46160  
-----

QAS: C. CONLEY #34517

TUCSON PD  
-----

Test	g/210L	Time
Air Blank	0.000	08:22:08
Cal Check	0.098	08:22:26
Air Blank	0.000	08:22:56

Clare Conley

Operator's Signature