

CERTIFICATION by CUSTODIAN of RECORDS

Tucson Police Department Crime Laboratory 1306 W. Miracle Mile, Tucson, AZ 85705

I hereby certify that these documents are a true and correct copy of the periodic maintenance* for the Intoxilyzer 8000 Serial Number 1.12.101 20-000332		
of Public Safety for a quality assurance program. A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.		
Date: 12-8-11 Cal 1 0.096 AC Cal 2 0.094 AC		
Between the dates of and and and the instrument was operating properly and accurately.		
The instrument was taken out of service for repair on N/A The instrument was returned to service on N/A		
Clau Comby		
Clare L. Conley		
12-8-11		
ADPS QAS Permit # 33314		
SUBSCRIBED AND SWORN to before me this 8, day of December, 2011. And L. Brecoff		
Notary Public OFFICIAL SEAL ANN L. BEECROFT NOTARY PUBLIC-ARIZONA		
Includes documents specified in A.A.C. R13-10-104.A.5 PIMA COUNTY My Comm. Exp. Jan. 4, 2013		

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

EXHIBIT G-4

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Clare Conley	AGENCY Tucson Police Department
1 -	TIME 6704
INTOXILYZER SERIAL # 80-600332	LOCATION ODW DU HO
(/) 1. Display reads "PUSH BUTTON TO	START"
DIAGNOSTIC TESTS	
(/) 1. Clock time check. (/) 2. Date check.	
OPERATIONAL TESTS	
 1. Alcohol-free subject test rest 2. Error recognition logic system Not a Successfully Completed 3. Proper sample recognition system Not a Successfully Completed Deficient sample printed. 4. Standard Calibration Check standard: 0. 094 AC 	em functioning. Test Sequence printed etem. Test Sequence printed
Instrument is operating properly and ac	curately. YESNO
COMMENTS Nme,	
SIGNATURE Care Conly	

DPS Form Exh G-4 (Rev. 05-01)

Location: TUCSON POLICE DEPT

Serial Number: 80-000332

12/08/2011

07:05:47

Standard Lot#: 04911100A1 Last Changed By: S. CIONE

#50616

QAS: C. CONLEY

#34517

TUCSON PD

Operator: C. CONLEY

#34517

TUCSON PD

Subject: TEST, SQAP

DOB: 11/11/1111

Sex: M

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	07:06:32
Diagnostic Test	Pass	07:06:59
Air Blank	0.000	07:07:27
0.100 Cal Check	0.096	07:07:46
Air Blank	0.000	07:08:15
Subject Test	0.000	07:08:44
Air Blank	ABT*	07:08:50
Air Blank	RFI**	07:08:58

^{*}Sequence Aborted

Not a Successfully Completed Test Sequence

^{**}RFI Detect

Location: TUCSON POLICE DEPT

Serial Number: 80-000332

12/08/2011 07:09:28

Standard Lot#: 04911100A1 Last Changed By: S. CIONE #50616

QAS: C. CONLEY #34517

TUCSON PD

Operator: C. CONLEY #34517
TUCSON PD

Subject: TEST, SQAP

DOB: 11/11/1111

Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	07:09:53
Air Blank	0.000	07:10:22

^{*}Improper Sample

Not a Successfully Completed Test Sequence

3567 Ur

134

Location: TUCSON POLICE DEPT

Serial Number: 80-000332

12/08/2011

07:10:25

Standard Lot#: 04911100A1

Last Changed By: S. CIONE

#50616

QAS: C. CONLEY

#34517

TUCSON PD

Operator: C. CONLEY

#34517

TPD

Subject: TEST, SQAP

DOB: 11/11/1111

Weight: 111

Sex: M

Weight. III

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	07:11:31
Diagnostic Test	Pass	07:11:59
Air Blank	0.000	07:12:26
0.100 Cal Check	0.096	07:12:45
Air Blank	0.000	07:13:15
Subject Test	0.259*	07:16:23
Air Blank	0.000	07:16:55
Air Blank	0.000	07:17:23
Subject Test	XXX**	07:17:39
Air Blank	0.000	07:18:13

^{*}Deficient Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

On 188

184

^{**}Invalid Sample

Location: TUCSON POLICE DEPT

Serial Number: 80-000332

12/08/2011

07:18:37

Standard Lot#: 04911100A1

Last Changed By: S. CIONE #50616

hase changed by. S. Clone #30010

QAS: C. CONLEY #34517

TUCSON PD

Test	g/210L	Time
Air Blank	0.000	07:19:04
Cal Check	0.094	07:19:23
Air Blank	0.000	07:19:53

134

Operator's Signature

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

154 12.12-2011

STANDARD CALIBRATION CHECK PROCEDURE

QA Special	ist Mare Conley AGENCY Tucson Police Department 8-11 TIME 0704
DATE 12	-8-11 TIME 0704
Intoxilyze	r Serial # 80-000332 Location ODW DUM HR
(🗸) 1.	Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
	Pour a standard alcohol concentration solution AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} + 0.2^{\circ}\text{C}$
	Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
(>) 3.	Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
(🗸) 4.	Air blank completed.
	Calibration check completed. Test results 0.096 AC.
	Air blank completed.
(/) 7.	Remove printed record. Attach the record to the completed checklist.
SIGNATURE	Clare Conly
DPS Form Ex	th G-2 (Rev. 05-01)

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

Location: TUCSON POLICE DEPT

Serial Number: 80-000332

12/08/2011

07:04:20

Standard Lot#: 04911100A1

Last Changed By: S. CIONE #50616

QAS: C. CONLEY #34517

TUCSON PD

Test	g/210L	Time
Air Blank	0.000	07:04:47
Cal Check	0.096	07:05:06
Air Blank	0.000	07:05:36