



CERTIFICATION by CUSTODIAN of RECORDS

Tucson Police Department Crime Laboratory
1306 W. Miracle Mile, Tucson, AZ 85705

Tsh
2.11.2016

I hereby certify that these documents are a true and correct copy of the periodic maintenance* for the Intoxilyzer 8000 Serial Number 90-000332 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.

A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.

Date: 2/2/2016 Cal 1 0.101 AC Cal 2 N/A AC

Between the dates of _____ and _____ the instrument was operating properly and accurately.

The instrument was taken out of service for repair on 2/2/2016.
The instrument was returned to service on _____.

William L Loomis

William L. Loomis

2/10/2016

ADPS QAS Permit # 33720

SUBSCRIBED AND SWORN to before me this 10th, day of February, 2016.

Elma Escarcega

Notary Public

*Includes documents specified in A.A.C. R13-10-104.A.5




EXHIBIT G-2

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

1 of 4
WLL
B4 2.11.2016

QA Specialist W. LOOMIS 53545 AGENCY Tucson Police Department
DATE 2/2/2016 TIME 13:45
Intoxilyzer Serial # 80-000332 Location CRIME LAB

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0. 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 101 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE William Loomis

DPS Form Exh G-2 (Rev. 05-01)

INTOXILYZER 8000
Location: TUCSON PD
Serial Number: 80-000332
Core Version: 8105.48
02/02/2016 13:45:41

Standard Lot#: 15615100A1
Last Changed By: W. LOOMIS #53545

QAS: W. LOOMIS #53545
TUCSON PD

Operator: W. LOOMIS #53545
TUCSON PD

Test	g/210L	Time
Air Blank	0.000	13:46:10
0.100 Cal Check	0.101	13:46:31
Air Blank	0.000	13:47:00

Preventative Maintenance
Performed -- Timer Reset

2 of 4
DH



TUCSON POLICE DEPARTMENT
CRIME LABORATORY

INTOXILYZER SERVICE RECORD

3 of 4
Bh 2.11.2016

Date: 2/2/2016

QAS: W. LOOMIS 53545

Intoxilyzer Serial Number: 80-000332

Reason For Service: AFTER SYSTEM START-UP, INSTRUMENT WILL GO TO READY MODE BUT START TEST BUTTON DOES NOT START SEQUENCE. ABLE TO ACCESS OPTIONS MENU, BUT ON EXITING OPTIONS, THE DISPLAY GETS "STUCK" ON LAST MENU OR SCREEN.

Date Out of Service: 2/2/2016

Calibration Check completed prior to removal from service

Date: 2/2/2016

SQAP completed prior to removal from service

Date: UNABLE TO COMPLETE

Corrective Action:

☒ Sent to CMI, Inc.

DOWNLOADED DATA PRIOR TO SENDING TO CMI



Service Evaluation Form

This form MUST be completed and enclosed with instrument to be serviced.

4 of 4 m.c.

Note: Please ship items in their original shipping container.

BH

Contact information: Customer Number 105075 (contact Customer Service)

Name TERRY GALLEGOS Phone: (520) 837-7677

Fax: (520) 791-4642 Email: TERRY.GALLEGOS@TUCSONAZ.GOV

Bill to Address:

FINANCE/ACCOUNTS PAYABLE
PO BOX 27450
TUCSON, AZ 85726-7450

Ship to Address:

TUCSON POLICE DEPARTMENT CRIMELAB
1306 W. MIRAGE MILE
TUCSON, AZ 85705

Instrument Serial Number: 80-000332

Detailed Description of Problem:

AFTER SYSTEM START-UP, INSTRUMENT WILL GO TO READY MODE, BUT START TEST
BUTTON DOES NOT START SEQUENCE. ABLE TO ACCESS OPTIONS MENU, BUT ON EXITING OPTIONS,
THE DISPLAY GETS STUCK ON LAST MENU OR SCREEN.

***Hazardous Material Warning! - DO NOT return gas cylinders with instrument! ***

Note: An evaluation fee of \$39.50 for handheld instruments or \$79.00 for IR evidential instruments will apply to estimates which are not repaired.

☐ I Authorize Repairs Up To: ☐ All ☐ \$250 ☐ \$500 ☐ \$750 ☐ Other \$ _____
Purchase Order Number 37768

Authorized By:

TERRY GALLEGOS TOXICOLOGY COORDINATOR

Name (Please Print)

Title

[Signature]

Signature

2.11.2016

Date

Ship to:

CMI, Inc.
Attn: Service Dept.
316 East Ninth Street
Owensboro, KY 42303

☒ No, please send estimate before repairs are made. An estimate will be faxed before performing any repairs and may cause delays in service.