

CERTIFICATION by CUSTODIAN of RECORDS

Tucson Police Department Crime Laboratory 270 S. Stone Avenue Tucson, AZ 85701-1917

2/0 S. Stone Avenue Tucson, AZ 85/01-191/
I hereby certify that these documents are a true and correct copy of the periodic maintenance* for the Intoxilyzer 8000 Serial Number Storing Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.
A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.
Date: 7.21-11 Cal 1 0.098 AC Cal 2 0.09% AC
Between the dates of
The instrument was taken out of service for repair on NIA. The instrument was returned to service on NIA.
Lare Conley
Clare L. Conley
7-21-11
ADPS QAS Permit # 33314
SUBSCRIBED AND SWORN to before me this 2/11, day of
Notary Public OFFICIAL SEAL SUSAN A. HITCHYE
*Includes documents specified in A.A.C. R13-10-104.A.5

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

My Comm. Exp. April 20, 2015

EXHIBIT G-4

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

1567 Cm By 2011 2.25.2011

QA SPECIALIST Clare Conley	AGENCY Tucson Police Department
DATE 7.21-11	TIME 1229
INTOXILYZER SERIAL # So-000334	LOCATION Crime Lab
(✓) 1. Display reads "PUSH BUTTON TO	
DIAGNOSTIC TESTS	
() 1. Clock time check. () 2. Date check.	
OPERATIONAL TESTS	
() 1. Alcohol-free subject test re () 2. Error recognition logic syst Not a Successfully Completed () 3. Proper sample recognition sy Not a Successfully Completed Deficient sample printed. () 4. Standard Calibration Check s Result: 0. 098 AC	em functioning. Test Sequence printed stem. Test Sequence printed
Instrument is operating properly and as	ccurately. YES NO
COMMENTS Now.	
SIGNATURE Man Conly	

DPS Form Exh G-4 (Rev. 05-01)

287 Cm

INTOXILYZER 8000

Location: TUCSON POLICE DEPT Serial Number: 80-000334

07/21/2011

12:30:50

Standard Lot#: 20910100A1

Last Changed By: S. CIONE

#50616

QAS: C. CONLEY

#34517

TUCSON PD

Operator: C. CONLEY

#34517

TUCSON PD

Subject: TEST, SQAP

DOB: 11/11/1111

Sex: M

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	12:31:45
Diagnostic Test	Pass	12:32:12
Air Blank	0.000	12:32:40
0.100 Cal Check	0.099	12:32:58
Air Blank	0.000	12:33:28
Subject Test	0.000	12:33:48
Air Blank	ABT*	12:33:54
Air Blank	0.000	12:34:23

^{*}Sequence Aborted

Not a Successfully Completed Test Sequence

Location: TUCSON POLICE DEPT

Serial Number: 80-000334

07/21/2011

12:34:26

Standard Lot#: 20910100A1 Last Changed By: S. CIONE

QAS: C. CONLEY #34517

TUCSON PD

#34517 Operator: C. CONLEY

TUCSON PD

Subject: TEST, SQAP

DOB: 11/11/1111

Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	12:34:52
Air Blank	IPS**	12:35:02

^{*}RFI Detect

Not a Successfully Completed Test Sequence

#50616

^{**}Improper Sample

Location: TUCSON POLICE DEPT

Serial Number: 80-000334 07/21/2011

12:35:50

Standard Lot#: 20910100A1
Last Changed By: S. CIONE #50616

QAS: C. CONLEY #34517
TUCSON PD

Operator: C. CONLEY #34517

Subject: TEST, SQAP DOB: 11/11/1111

TUCSON PD

Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	12:36:39
Diagnostic Test	Pass	12:37:06
Air Blank	0.000	12:37:34
0.100 Cal Check	0.098	12:37:52
Air Blank	0.000	12:38:22
Subject Test	0.154*	12:41:30
Air Blank	0.000	12:42:02
Air Blank	0.000	12:42:30
Subject Test	XXX**	12:43:13
Air Blank	0.000	12:43:46

^{*}Deficient Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

4567 Cn

Bh

^{**}Invalid Sample

Location: TUCSON POLICE DEPT

Serial Number: 80-000334

07/21/2011

12:44:04

Standard Lot#: 20910100A1

Last Changed By: S. CIONE #50616

QAS: C. CONLEY #34517

TUCSON PD

Test	g/210L	Time
Air Blank	0.000	12:44:31
Cal Check	0.098	12:44:50
Air Blank	0 000	12.45.19

Operator's Signature

storage devices of the Tucson Police Department.

EXHIBIT G-2

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

6-67 Ch

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

7.25.2011

STANDARD CALIBRATION CHECK PROCEDURE

QA Speciali	st Clare Con bey	AGENCY	Tucson Police Department
DATE 7.	21-11	TIME	1229
Intoxilyzer	Serial # <u>\$0-600334</u>	Location	rume Lab
(🗸) - 1.	Ensure that gas tank is at standard alcohol concentration OR		tion 0. 100 AC.
	into a clean dry simulator	and asseen made.	ation solution AC, emble the simulator. Ensure Turn on the simulator and 0.2°C
(🗸) 2.	Intoxilyzer 8000 displa	ay reads	"PUSH BUTTON TO START"
(,) 3.		for we	nu". Select "D" for dry t control test. After
(/) 4.	Air blank completed.		
(✓) 5.	Calibration check compl	leted. T	est results 0. <u>098</u> AC.
(/) 6.	Air blank completed.		
	Remove printed record. checklist.	Attach	the record to the completed
SIGNATURE	Clare Conly		· .
DPS Form Ex	kh G-2 (Rev. 05-01)		

Location: TUCSON POLICE DEPT

Serial Number: 80-000334

07/21/2011

12:29:28

7587 Cm 154

Standard Lot#: 20910100A1

Last Changed By: S. CIONE #50616

QAS: C. CONLEY #34517

TUCSON PD

Test	g/210L	Time
Air Blank	0.000	12:29:55
Cal Check	0.098	12:30:14
Air Blank	0.000	12:30:43

Uare linky

Operator's Signature

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.