



CERTIFICATION by CUSTODIAN of RECORDS
Tucson Police Department Crime Laboratory

TS4
8-8-07

I hereby certify that these documents are a true and correct copy of the periodic maintenance* for the Intoxilyzer 8000 Serial Number 80-500712 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.

A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.

Date: 8-2-07 Cal 1 0.098 AC Cal 2 0.097 AC

Between the dates of 7-6-07 and 8-2-07 the instrument was operating properly and accurately.

The instrument was taken out of service for repair on N/A.
The instrument was returned to service on N/A.


Clare Conley

Clare L. Conley

8-3-07

ADPS QAS Permit # 33314

SUBSCRIBED AND SWORN to before me this 3, day of August, 2007.



Kathy A Kinney

Notary Public

*Includes documents specified in A.A.C. R13-10-104.A.5

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

134
8-8-07

QA SPECIALIST Clare Conley AGENCY Tucson Police Department

DATE 8-2-07 TIME 09:05

INTOXILYZER SERIAL # 80-000712 LOCATION Traffic

(✓) 1. Display reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- (✓) 1. Clock time check.
- (✓) 2. Date check.

OPERATIONAL TESTS

- (✓) 1. Alcohol-free subject test result 0. 000 AC.
- (✓) 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
- (✓) 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
- (✓) 4. Standard Calibration Check standard 0. 100 AC
Result: 0. 097 AC

Instrument is operating properly and accurately. YES NO

COMMENTS None.

SIGNATURE Clare Conley

248 cu

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number : 80-000712
08/02/2007

09:06:23 **B4**

Standard Lot#: 635301I
Last Changed By: G. EPPLEY #26227

QAS: C. CONLEY #34517
TUCSON PD

Operator: C. CONLEY #34517
TUCSON PD

Subject: TEST, SQAP
DOB: 11/11/1111
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:07:10
Diagnostic Test	Pass	09:07:37
Air Blank	0.000	09:08:04
0.100 Cal Check	0.098	09:08:23
Air Blank	0.000	09:08:52
Subject Test	0.000	09:09:10
Air Blank	ABT*	09:09:17
Air Blank	0.000	09:09:45

*Sequence Aborted

Not a Successfully Completed Test Sequence

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

368cm

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number : 80-000712
08/02/2007

09:09:49

Bh

Standard Lot#: 635301I
Last Changed By: G. EPPLEY #26227

QAS: C. CONLEY #34517
TUCSON PD

Operator: C. CONLEY #34517
TUCSON PD

Subject: TEST, SQAP
DOB: 11/11/1111
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:10:38
Diagnostic Test	Pass	09:11:06
Air Blank	RFI*	09:11:14
Air Blank	0.000	09:11:42

*RFI Detect

Not a Successfully Completed Test Sequence

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

4 of 8 cu

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number : 80-000712
08/02/2007

09:11:47

B4

Standard Lot#: 635301I
Last Changed By: G. EPPLEY #26227

QAS: C. CONLEY #34517
TUCSON PD

Operator: C. CONLEY #34517
TUCSON PD

Subject: TEST, SQAP
DOB: 11/11/1111
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	09:12:17
Air Blank	0.000	09:12:45

*Improper Sample

Not a Successfully Completed Test Sequence

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

Sof8cu

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number : 80-000712
08/02/2007

09:12:50

TS4

Standard Lot#: 635301I
Last Changed By: G. EPPLEY #26227

QAS: C. CONLEY #34517
TUCSON PD

Operator: C. CONLEY #34517
TUCSON PD

Subject: TEST, SQAP
DOB: 11/11/1111
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:13:41
Diagnostic Test	Pass	09:14:09
Air Blank	0.000	09:14:36
0.100 Cal Check	0.098	09:14:55
Air Blank	0.000	09:15:24
Subject Test	0.318*	09:18:32
Air Blank	0.000	09:19:05
Air Blank	0.000	09:19:32
Subject Test	XXX**	09:19:50
Air Blank	0.000	09:20:24

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

6.78 cu

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number: 80-000712
08/02/2007 09:20:44

PH

Standard Lot#: 635301I
Last Changed By: G. EPPLEY #26227

QAS: C. CONLEY #34517
TUCSON PD

Test	g/210L	Time
Air Blank	0.000	09:21:12
Cal Check	0.097	09:21:30
Air Blank	0.000	09:22:00

Chase Conley
Operator's Signature

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

154

STANDARD CALIBRATION CHECK PROCEDURE

QA Specialist Clare Conley AGENCY Tucson Police Department
DATE 8-2-07 TIME 09:05
Intoxilyzer Serial # 80-000712 Location Traffic

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0. 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 098 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Clare Conley

DPS Form Exh G-2 (Rev. 05-01)

8 of 8 cu

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number: 80-000712
08/02/2007 09:05:02

TS4

Standard Lot#: 635301I
Last Changed By: G. EPPLEY #26227

QAS: C. CONLEY #34517
TUCSON PD

Test	g/210L	Time
Air Blank	0.000	09:05:30
Cal Check	0.098	09:05:49
Air Blank	0.000	09:06:18

Clare Conley

Operator's Signature