



CERTIFICATION by CUSTODIAN of RECORDS

Tucson Police Department Crime Laboratory  
270 S. Stone Avenue Tucson, AZ 85701-1917

TS4  
4.10.09

I hereby certify that these documents are a true and correct copy of the periodic maintenance\* for the Intoxilyzer 8000 Serial Number 80-000712 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.

A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.

Date: 3/26/09 Cal 1 N/A AC Cal 2 N/A AC

Between the dates of - and - the instrument was operating properly and accurately.

The instrument was taken out of service for repair on 3/20/09.  
The instrument was returned to service on \_\_\_\_\_.

Sarah Cione

Sarah L. Cione


3/31/09

ADPS QAS Permit # 33355

SUBSCRIBED AND SWORN to before me this 31, day of March, 2009.

A.M. Hernandez

Notary Public



\*Includes documents specified in A.A.C. R13-10-104.A.5



TUCSON POLICE DEPARTMENT  
CRIME LABORATORY

INTOXILYZER SERVICE RECORD

TSY  
4.10.09

Date: 3/26/09

QAS: Sarah Cione

Intoxilyzer Serial Number: 80-000712

Reason For Service: Breath hose does not heat up, will not pass diagnostics (temp reg. failure). Hose is bad (tests "OPEN" on breath hose tester), but a good hose does not heat up either. Also, internal printer paper gets jammed inside of door when printing. Works fine if printer door is removed.

Date Out of Service: 3/20/08

Calibration Check completed prior to removal from service

Date: N/A

SQAP completed prior to removal from service

Date: N/A

Corrective Action:  Sent to CMI, Inc.

downloaded 10 subject tests prior to sending to CMI



# Return Material Authorization Form

By 4-10-09

This form **MUST** be completed and enclosed with item(s) being returned.  
Failure to complete and return this form may cause delays in crediting / repairing.

1. Item(s) Returned:  Instrument  Supplies  Other  
(Note: please ship items in their original shipping container(s) or a similar protective box.)

2. Instrument Model: Intox 8000 Serial Number 80-000712

3. Bill To Address: Finance / Accounts Payable Ship To Address: (NO PO BOXES)  
P.O. Box 27450 Tucson Police Department  
TUCSON, AZ 85726-7450 Crime Lab  
270 S. Stone Ave  
TUCSON, AZ 85701-1917

Customer # 105075  
Name Terry Gallegos  
Phone: (520) 837-7677 Fax: (520) 791-4642  
Email: Terry.Gallegos@tucsonaz.gov

4. Reason for Return: Note - If you are sending an item for repair, please give a detailed description of the problem. Please list any special instructions that you may have concerning this return. Note - Returns for credit may be subject to a 15% restocking charge.

Breath hose does not heat up, will not pass diagnostics. Hose is bad (tested "OPEN" on breath hose tester), but a good hose does not heat up on intox either. Also, internal printer gets jammed inside of door when printing. Works fine if printer door is removed.

5. Purchase Order Number (attach a copy of P.O. if applicable) 20669

6. I Authorize All Repairs:  Yes  No

7. I Authorize Repairs Up To: \$                     

8. I Require An Estimate Regardless of Cost  Yes  No

Note: An evaluation fee will apply to estimates that are not repaired.

9. Authorized By: Terry Gallegos

Name (Please Print)

Terry Gallegos  
Signature

3.26.2009  
Date