

CERTIFICATION by CUSTODIAN of RECORDS
Tucson Police Department Crime Laboratory
270 S. Stone Avenue Tucson, AZ 85701-1917

I hereby certify that these documents are a true and correct copy of the
periodic maintenance* for the Intoxilyzer 8000 Serial Number
80-000712 maintained by the Tucson Police Department
Crime Laboratory pursuant to the requirements of the Arizona Department
of Public Safety for a quality assurance program.
A quality assurance specialist (QAS) tested the instrument on the following
dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.3
for a 0.100AC solution are 0.090AC to 0.110AC.
Date: 5/26/11 Cal 1 0.100 AC Cal 2 0.098 A
Date 5/26/11 Cal 1 0.100 AC Cal 2 0.098 A
D
Between the dates of $\frac{4/28/n}{}$ and $\frac{5/26/l}{}$ the instrument was operating properly and accurately.
the instrument was operating properly and accurately.
The instrument was taken out of service for repair on
The instrument was returned to service on -
Earth Ciant
Sarah Cione
Sarah L. Cione
Saran L. Cione
5/30/11
, .
ADPS QAS Permit # 33355
SUBSCRIBED AND SWORN to before me this 30° , day of
OFFICIAL SEAL ALL 20 11
SUE A. SOZA
NOTARY PUBLIC-ARIZONA
PIMA COUNTY Dee 4 00%
My Comm. Exp. Dec. 30, 2012
Notary Public

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

^{*}Includes documents specified in A.A.C. R13-10-104.A.5

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

154 5.31.2011

QA SPECIALIST Sarah Cione	AGENCY Tucson Police Department
DATE 5/26/11	TIME 826
INTOXILYZER SERIAL # 80-000712	LOCATION Team 2-0DW
(\checkmark) 1. Display reads "PUSH BUTTON TO	START"
DIAGNOSTIC TESTS	
(\checkmark) 1. Clock time check. (\checkmark) 2. Date check.	
OPERATIONAL TESTS	
<pre>(/) 1. Alcohol-free subject test res (/) 2. Error recognition logic syste</pre>	tem. Test Sequence printed andard 0. 100 AC
SIGNATURE Sonah Cionl	

DPS Form Exh G-4 (Rev. 05-01)

INTOXILYZER 8000

Location: TUCSON POLICE DEPT Serial Number: 80-000712

05/26/2011

08:27:40

Standard Lot#: 719919 Last Changed By: S. CLEARY #44524 QAS: S. CIONE #50616 TUCSON PD Operator: S. CIONE #50616 TUCSON PD

Subject: SQAP, SQAP DOB: 11/11/1111

Sex: F

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:28:24
Diagnostic Test	Pass	08:28:51
Air Blank	0.000	08:29:18
0.100 Cal Check	0.099	08:29:37
Air Blank	0.000	08:30:07
Subject Test	0.000	08:30:35
Air Blank	0.000	08:31:02
Five Minute Wait	ABT*	
Air Blank	0.000	08:31:36

^{*}Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number: 80-000712

05/26/2011

08:31:44

Standard Lot#: 719919

Last Changed By: S. CLEARY #44524

QAS: S. CIONE #50616

TUCSON PD

Operator: S. CIONE #50616

TUCSON PD

Subject: SQAP, SQAP DOB: 11/11/1111

Sex: F

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:32:27
Diagnostic Test	Pass	08:32:54
Air Blank	RFI*	08:33:04
Air Blank	IPS**	08:33:12

^{*}RFI Detect

Not a Successfully Completed Test Sequence

^{**}Improper Sample

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number: 80-000712

05/26/2011

08:33:19

154

Standard Lot#: 719919	
Last Changed By: S. CLEARY	#44524
QAS: S. CIONE TUCSON PD	#50616
Operator: S. CIONE TUCSON PD	#50616

Subject: SQAP, SQAP DOB: 11/11/1111

Sex: F

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:34:01
Diagnostic Test	Pass	08:34:28
Air Blank	0.000	08:34:55
0.100 Cal Check	0.099	08:35:14
Air Blank	0.000	08:35:44
Subject Test	0.064*	08:38:52
Air Blank	0.000	08:39:24
Air Blank	0.000	08:39:51
Subject Test	XXX * *	08:40:07
Air Blank	0.000	08:40:40

^{*}Deficient Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

^{**}Invalid Sample

5 of 750

INTOXILYZER 8000

Location: TUCSON POLICE DEPT Serial Number: 80-000712

05/26/2011

08:40:57

Standard Lot#: 719919

Last Changed By: S. CLEARY #44524

QAS: S. CIONE #50616

TUCSON PD

 By

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

6 of 7sc 134 5.31.2011

QA Spec	ciali	st Sarah Cionu	AGENCY Tucson Police Department
DATE _	5/26	2/11	TIME _826
	/	Serial # <u>80-000712</u>	
(1.	OR	
/		into a clean dry simulator that a tight seal has been allow temperature to reach	
(🗸)	2.	Intoxilyzer 8000 display	y reads "PUSH BUTTON TO START"
(🗸)	3.	Go to the "Control Test	ting Menu". Select "D" for dry for wet control test. After
(V)		Air blank completed.	
()	5.	Calibration check comple	eted. Test results 0. 100 AC.
	6.	Air blank completed.	
	7.	Remove printed record. checklist.	Attach the record to the completed
SIGNAT	URE /	Sarah Cione	

DPS Form Exh G-2 (Rev. 05-01)

7 of 7 sc

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number: 80-000712

05/26/2011

08:26:17

Standard Lot#: 719919

Last Changed By: S. CLEARY #44524

QAS: S. CIONE #50616

TUCSON PD

Test	i	g/210L	Time
Air	Blank	0.000	08:26:44
Cal	Check	0.100	08:27:03
Air	Blank	0.000	08:27:32