



CERTIFICATION by CUSTODIAN of RECORDS

Tucson Police Department Crime Laboratory
270 S. Stone Avenue Tucson, AZ 85701-1917

By
1-25-2010

I hereby certify that these documents are a true and correct copy of the periodic maintenance* for the Intoxilyzer 8000 Serial Number 80-001163 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.

A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.

Date: 1/22/2010 Cal 1 0.097 AC Cal 2 0.097 AC

Between the dates of 12/11/09 and 1/22/10 the instrument was operating properly and accurately.

The instrument was taken out of service for repair on -.
The instrument was returned to service on -.

Sarah Cione

Sarah L. Cione

1/25/2010

ADPS QAS Permit # 33355

SUBSCRIBED AND SWORN to before me this 25th, day of January, 2010.

A.M. Hernandez

Notary Public



*Includes documents specified in A.A.C. R13-10-104.A.5

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

By
1-25-2010

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Sarah Cione AGENCY Tucson Police Department
DATE 1/22/2010 TIME 1249
INTOXILYZER SERIAL # 80-001163 LOCATION Crime Lab

(✓) 1. Display reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- (✓) 1. Clock time check.
- (✓) 2. Date check.

OPERATIONAL TESTS

- (✓) 1. Alcohol-free subject test result 0. 000 AC.
- (✓) 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
- (✓) 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
- (✓) 4. Standard Calibration Check standard 0. 100 AC
Result: 0. 097 AC

Instrument is operating properly and accurately. YES NO

COMMENTS More than 31 days b/n Calibration Checks, but no Evidentiary Tests performed since 12/11/09

SIGNATURE Sarah Cione

2 of 9 sc

INTOXILYZER 8000
Location: ARIZONA
Serial Number : 80-001163
01/22/2010

12:51:02

BA

Standard Lot#: 808101G
Last Changed By: S. CIONE #50616

QAS: S. CIONE #50616
TUCSON PD

Operator: S. CIONE #50616
TUCSON PD

Subject: SQAP, SQAP
DOB: 11/11/1111
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	12:51:50
Diagnostic Test	Pass	12:52:17
Air Blank	0.000	12:52:45
0.100 Cal Check	0.097	12:53:03
Air Blank	0.000	12:53:33
Subject Test	0.000	12:54:05
Air Blank	0.000	12:54:32
Five Minute Wait	ABT*	
Air Blank	0.000	12:55:09

*Sequence Aborted

Not a Successfully Completed Test Sequence

3 of 9 sc

INTOXILYZER 8000
Location: ARIZONA
Serial Number : 80-001163
01/22/2010 12:56:12

Standard Lot#: 808101G
Last Changed By: S. CIONE #50616

174

QAS: S. CIONE #50616
TUCSON PD

Operator: S. CIONE #50616
TUCSON PD

Subject: SQAP, SQAP
DOB: 11/11/1111
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	12:56:43
Air Blank	IPS**	12:56:52

*RFI Detect
**Improper Sample

Not a Successfully Completed Test Sequence

4 of 9 sc

INTOXILYZER 8000
 Location: ARIZONA
 Serial Number : 80-001163
 01/22/2010 12:57:45

Standard Lot#: 808101G
 Last Changed By: S. CIONE #50616

QAS: S. CIONE #50616
 TUCSON PD

Operator: S. CIONE #50616
 TUCSON PD

Subject: SQAP, SQAP
 DOB: 11/11/1111
 Sex: F Weight: 111

15 Minute Deprivation Period? Yes

134

Test	g/210L	Time
Air Blank	0.000	12:58:30
Diagnostic Test	Pass	12:58:57
Air Blank	0.000	12:59:24
0.100 Cal Check	0.097	12:59:43
Air Blank	0.000	13:00:13
Subject Test	0.043*	13:03:21
Air Blank	0.000	13:03:53
Air Blank	0.000	13:04:20
Subject Test	XXX**	13:04:42
Air Blank	0.000	13:05:15

*Deficient Sample
 **Invalid Sample
 Begin new deprivation period

Not a Successfully Completed Test Sequence

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

5 of 9 sc

INTOXILYZER 8000
Location: ARIZONA
Serial Number: 80-001163
01/22/2010 13:05:28

Standard Lot#: 808101G
Last Changed By: S. CIONE #50616

TS4

QAS: S. CIONE #50616
TUCSON PD

Test	g/210L	Time
Air Blank	0.000	13:05:55
Cal Check	0.097	13:06:14
Air Blank	0.000	13:06:44

Sarah Cione
Operator's Signature

EXHIBIT G-2

6 of 9 sc

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

By
1-25-2010

QA Specialist Sarah Cione AGENCY Tucson Police Department
DATE 1/22/2010 TIME 1249
Intoxilyzer Serial # 80-00163 Location Crime Lab

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0. 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 097 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE

Sarah Cione

DPS Form Exh G-2 (Rev. 05-01)

7 of 9 SC

INTOXILYZER 8000
Location: ARIZONA
Serial Number: 80-001163
01/22/2010 12:49:40

Standard Lot#: 808101G
Last Changed By: S. CIONE #50616



QAS: S. CIONE #50616
TUCSON PD

Test	g/210L	Time
Air Blank	0.000	12:50:07
Cal Check	0.097	12:50:26
Air Blank	0.000	12:50:56

Sarah Cione
Operator's Signature



TUCSON POLICE DEPARTMENT
CRIME LABORATORY

8 of 9 SC

INTOXILYZER SERVICE RECORD
FOR
CALIBRATION GAS TANK CHANGE

By
1-25-2010

Date: 1/22/2010
 QAS: Sarah Cione
 Intoxilyzer Serial Number: 80-001163
 New Lot Number: 808101G
 New Expiration Date: 3/31/2010

- (✓) 1. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 2. Remove the old gas tank and install the new one.
- (✓) 3. Hit the "ESC" key twice and swipe the QAS card.
- (✓) 4. Go to the "M-Maintenance" menu and select "S-Standard Change." Select "D" for dry gas and enter the following information:
 Enter Target Value: 0.100
 Enter Control Lot #: From the tank label
 Enter Exp Date: From the tank label (mm/dd/yyyy)
 Edit info Y/N? Y or N
 Update Std Info? Y ENTER
 "Please wait saving..."
- (✓) 5. Hit the "ESC" key twice.
- (✓) 6. Go to the "M-Maintenance" menu and select "Control Testing." Select "D" for dry control test. After selection is made press ENTER.
- (✓) 7. Air blank completed.
- (✓) 8. Calibration check completed. Test results 0.097 AC.
- (✓) 9. Air blank completed.
- (✓) 10. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Sarah Cione

9 of 9
SC

INTOXILYZER 8000
Location: ARIZONA
Serial Number: 80-001163
01/22/2010 12:48:18

104

Standard Lot#: 808101G
Last Changed By: S. CIONE #50616

QAS: S. CIONE #50616
TUCSON PD

Test	g/210L	Time
Air Blank	0.000	12:48:45
Cal Check	0.097	12:49:04
Air Blank	0.000	12:49:34

Sarah Cione
Operator's Signature