



CERTIFICATION by CUSTODIAN of RECORDS

Tucson Police Department Crime Laboratory
270 S. Stone Avenue Tucson, AZ 85701-1917

TS4
2-10-2009

I hereby certify that these documents are a true and correct copy of the periodic maintenance* for the Intoxilyzer 8000 Serial Number 80-001165 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.

A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.

Date: 2-6-09 Cal 1 0.100 AC Cal 2 0.099 AC

Between the dates of 01-09-09 and 2-6-09 the instrument was operating properly and accurately.

The instrument was taken out of service for repair on N/A.
The instrument was returned to service on N/A.

Clare Conley

Clare L. Conley

2-10-09

ADPS QAS Permit # 33314

SUBSCRIBED AND SWORN to before me this 10th, day of February, 2009.

A.M. Hernández

Notary Public



OFFICIAL SEAL
A.M. HERNÁNDEZ
NOTARY PUBLIC-ARIZONA
PIMA COUNTY
My Comm. Exp. Mar. 17, 2012

*Includes documents specified in A.A.C. R13-10-104.A.5

EXHIBIT G-4

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

1069
cu
By
2-10-2009

QA SPECIALIST Clare Conley AGENCY Tucson Police Department

DATE 2-6-09 TIME 0751

INTOXILYZER SERIAL # S0-001165 LOCATION ODS

(✓) 1. Display reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- (✓) 1. Clock time check.
- (✓) 2. Date check.

OPERATIONAL TESTS

- (✓) 1. Alcohol-free subject test result 0. 000 AC.
- (✓) 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
- (✓) 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
- (✓) 4. Standard Calibration Check standard 0. 100 AC
Result: 0. 099 AC

Instrument is operating properly and accurately. YES NO

COMMENTS None.

SIGNATURE Clare Conley

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number : 80-001165
02/06/2009 07:53:58

2669
CW
TSY

Standard Lot#: 804302I
Last Changed By: S. CIONE #50616

QAS: C. CONLEY #34517
TUCSON PD

Operator: C. CONLEY #34517
TUCSON PD

Subject: TEST, SQAP
DOB: 11/11/1111
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	07:54:52
Diagnostic Test	Pass	07:55:19
Air Blank	0.000	07:55:46
0.100 Cal Check	0.100	07:56:05
Air Blank	0.000	07:56:36
Subject Test	0.000	07:57:04
Air Blank	ABT*	07:57:11
Air Blank	0.000	07:57:39

*Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number : 80-001165
02/06/2009 07:57:44

369
cu
B4

Standard Lot#: 804302I
Last Changed By: S. CIONE #50616

QAS: C. CONLEY #34517
TUCSON PD

Operator: C. CONLEY #34517
TUCSON PD

Subject: TEST, SQAP
DOB: 11/11/1111
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	07:58:13
Air Blank	IPS**	07:58:22

*RFI Detect
**Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number : 80-001165
02/06/2009

07:58:30

469
CW

Standard Lot#: 804302I
Last Changed By: S. CIONE #50616

QAS: C. CONLEY #34517
TUCSON PD

Operator: C. CONLEY #34517
TUCSON PD

Subject: TEST, SQAP
DOB: 11/11/1111
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	07:59:13
Diagnostic Test	Pass	07:59:41
Air Blank	0.000	08:00:08
0.100 Cal Check	0.101	08:00:27
Air Blank	0.000	08:00:57
Subject Test	0.238*	08:04:06
Air Blank	0.000	08:04:39
Air Blank	0.000	08:05:06
Subject Test	RGE**	08:05:29
Air Blank	0.000	08:06:04

*Deficient Sample

**Range Exceeded

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number : 80-001165
02/06/2009 08:06:08

5-69
104 cu

Standard Lot#: 804302I
Last Changed By: S. CIONE #50616

QAS: C. CONLEY #34517
TUCSON PD

Operator: C. CONLEY #34517
TUCSON PD

Subject: TEST, SQAP
DOB: 11/11/1111
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:06:53
Diagnostic Test	Pass	08:07:20
Air Blank	0.000	08:07:47
0.100 Cal Check	0.098	08:08:06
Air Blank	0.000	08:08:36
Subject Test	XXX*	08:08:59
Air Blank	IPS**	08:09:13

*Invalid Sample
Begin new deprivation period
**Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number: 80-001165
02/06/2009 08:09:32

b-69
Cu
784

Standard Lot#: 804302I
Last Changed By: S. CIONE #50616

QAS: C. CONLEY #34517
TUCSON PD

Test	g/210L	Time
Air Blank	AMB*	08:10:00
Air Blank	0.000	08:10:28

*Ambient Fail

Clare Conley

Operator's Signature

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number: 80-001165

02/06/2009

08:10:33

769
CW

Standard Lot#: 804302I

Last Changed By: S. CIONE #50616

TS4

QAS: C. CONLEY #34517

TUCSON PD

Test g/210L Time

Air Blank 0.000 08:11:00
Cal Check 0.099 08:11:19
Air Blank 0.000 08:11:49

C. Conley

Operator's Signature

EXHIBIT G-2

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

8569
TSY
2-10-09
CU

QA Specialist Clare Conley AGENCY Tucson Police Department
DATE 2-6-09 TIME 0751
Intoxilyzer Serial # 80-001165 Location ODS

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0. 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 100 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Clare Conley

DPS Form Exh G-2 (Rev. 05-01)

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number: 80-001165
02/06/2009 07:51:52

Standard Lot#: 804302I
Last Changed By: S. CIONE #50616

QAS: C. CONLEY #34517
- TUCSON PD

929
cu
TS4

Test	g/210L	Time
Air Blank	0.000	07:52:20
Cal Check	0.100	07:52:38
Air Blank	0.000	07:53:08

Clare Conley

Operator's Signature