



CERTIFICATION by CUSTODIAN of RECORDS

Tucson Police Department Crime Laboratory
270 S. Stone Avenue Tucson, AZ 85701-1917

I hereby certify that these documents are a true and correct copy of the periodic maintenance* for the Intoxilyzer 8000 Serial Number 80-001199 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.

TS4
3-9-2011

A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.

Date: 3/3/11 Cal 1 0.103 AC Cal 2 0.100 AC

Between the dates of 2/3/11 and 3/3/11 the instrument was operating properly and accurately.

The instrument was taken out of service for repair on -
The instrument was returned to service on -

Sarah Cione

Sarah L. Cione


3/9/11

ADPS QAS Permit # 33355

SUBSCRIBED AND SWORN to before me this 9 day of March, 2011.

Am Hernandez

Notary Public



*Includes documents specified in A.A.C. R13-10-104.A.5

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

EXHIBIT G-4

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R13-10-104(A)

1 of 7
sc

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

154
3.9.2011

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Sarah Cione AGENCY Tucson Police Department
DATE 3/3/11 TIME 935
INTOXILYZER SERIAL # 80-001199 LOCATION Vehicle # 942

(✓) 1. Display reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

(✓) 1. Clock time check.
(✓) 2. Date check.

OPERATIONAL TESTS

(✓) 1. Alcohol-free subject test result 0. 000 AC.
(✓) 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
(✓) 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
(✓) 4. Standard Calibration Check standard 0. 100 AC
Result: 0. 100 AC

Instrument is operating properly and accurately. YES ✓ NO

COMMENTS

SIGNATURE Sarah Cione

2 of 7_{SC}

BS4

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number: 80-001199
03/03/2011 09:37:20

Standard Lot#: 719919
Last Changed By: C. MORAND #24219

OAS: S. CIONE #50616
TUCSON PD

Operator: S. CIONE #50616
TUCSON PD

Subject: SQAP, SQAP
DOB: 11/11/1111
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:38:19
Diagnostic Test	Pass	09:38:48
Air Blank	0.000	09:39:15
0.100 Cal Check	0.099	09:39:34
Air Blank	0.000	09:40:03
Subject Test	0.000	09:40:19
Air Blank	0.000	09:40:46
Five Minute Wait	ABT*	
Air Blank	0.000	09:41:21

*Sequence Aborted

Not a Successfully
Completed
Test Sequence

3 of 7sc

TS4

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number: 80-001199
03/03/2011 09:41:30

Standard Lot#: 719919
Last Changed By: C. MORAND #24219

QAS: S. CIONE #50616
TUCSON PD

Operator: S. CIONE #50616
TUCSON PD

Subject: SQAP, SQAP
DOB: 11/11/1111
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	09:42:06
Air Blank	IPS**	09:42:14

*RFI Detect
**Improper Sample

Not a Successfully
Completed
Test Sequence

4 of 7sc

TS4

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number: 80-001199
03/03/2011 09:42:28

Standard Lot#: 719919
Last Changed By: C. MORAND #24219

OAS: S. CIONE #50616
TUCSON PD

Operator: S. CIONE #50616
TUCSON PD

Subject: SQAP, SQAP
DOB: 11/11/1111
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:43:14
Diagnostic Test	Pass	09:43:43
Air Blank	0.000	09:44:10
0.100 Cal Check	0.099	09:44:29
Air Blank	0.000	09:44:58
Subject Test	0.097*	09:48:07
Air Blank	0.000	09:48:38
Air Blank	0.000	09:49:06
Subject Test	XXX**	09:49:23
Air Blank	0.000	09:49:56

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully
Completed
Test Sequence

5 of 7
sc

154

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number: 80-001199
03/03/2011 09:50:14

Standard Lot#: 719919
Last Changed By: C. MORAND #24219

DRS: S. CIONE #50616
TUCSON PD

Test	g/210L	Time
Air Blank	0.000	09:50:41
Cal Check	0.100	09:51:00
Air Blank	0.000	09:51:30

Sarah Cione
Operator's Signature

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

TSY
3.9.2011

QA Specialist Sarah Cione AGENCY Tucson Police Department
DATE 3/3/11 TIME 935
Intoxilyzer Serial # 80-001199 Location Vehicle #942

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0. 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 103 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Sarah Cione

DPS Form Exh G-2 (Rev. 05-01)

7 of 7sc

ph

INTOXILYZER 8000
Location: TUCSON POLICE DEPT
Serial Number: 80-001199
03/03/2011 09:35:57

Standard Lot#: 719919
Last Changed By: C. MORAND #24219

QAS: S. CIONE #50616
TUCSON PD

Test	g/210L	Time
Air Blank	0.000	09:36:24
Cal Check	0.103	09:36:43
Air Blank	0.000	09:37:12

Sarah Cione
Operator's Signature