



CERTIFICATION by CUSTODIAN of RECORDS
 Tucson Police Department Crime Laboratory
 270 S. Stone Avenue Tucson, AZ 85701-1917

I hereby certify that these documents are a true and correct copy of the periodic maintenance* for the Intoxilyzer 8000 Serial Number 80-002201 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.

TSY
10-26-09

A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.

Date: 10/13/09 Cal 1 - AC Cal 2 - AC

Between the dates of - and - the instrument was operating properly and accurately.

The instrument was taken out of service for repair on 10/10/09.
 The instrument was returned to service on _____.

Sarah Cione

Sarah L. Cione


10/14/09

ADPS QAS Permit # 33355

SUBSCRIBED AND SWORN to before me this 14th day of October, 2009.

A.M. Hernandez

Notary Public



*Includes documents specified in A.A.C. R13-10-104.A.5

1 of 2 sc



Return Material Authorization Form

This form MUST be completed and enclosed with item(s) being returned.
Failure to complete and return this form may cause delays in crediting / repairing.

TS4
10-26-09

1. Item(s) Returned: Instrument Supplies Other
(Note: please ship items in their original shipping container(s) or a similar protective box.)

2. Instrument Model: I8000 Serial Number 80-002201

3. Bill To Address: Finance / Accounts Payable P.O. Box 27450 TUCSON, AZ 85726-7450
Ship To Address: (NO PO BOXES) Tucson Police Department Crime Lab 270 S. Stone Ave TUCSON, AZ 85701-1917

Customer # 105075
Name Terry Gallegos
Phone: (520) 837-7677 Fax: (520) 791-4642
Email: Terry.Gallegos@tucsonaz.gov

4. Reason for Return: Note - If you are sending an item for repair, please give a detailed description of the problem. Please list any special instructions that you may have concerning this return. Note - Returns for credit may be subject to a 15% restocking charge.
Breath hose does not heat up. Breath hose is good - checked it on tester.
Fails diagnostics - temperature regulation fails.

5. Purchase Order Number (attach a copy of P.O. if applicable) 20669

6. I Authorize All Repairs: Yes No

7. I Authorize Repairs Up To: \$ _____

8. I Require An Estimate Regardless of Cost Yes No

Note: An evaluation fee will apply to estimates that are not repaired.

9. Authorized By: Terry Gallegos
Name (Please Print)
[Signature] 10-13-09
Signature Date



TUCSON POLICE DEPARTMENT
CRIME LABORATORY

INTOXILYZER SERVICE RECORD

2 of 2 SC

TSY
10-26-09

Date: 10/13/09

QAS: Sarah Cione

Intoxilyzer Serial Number: 80-002201

Reason For Service: Breath hose does not heat up.

Date Out of Service: 10/10/09

Calibration Check completed prior to removal from service

Date: -

SQAP completed prior to removal from service

Date: -

Corrective Action:

Sent to CMI, Inc.

downloaded data prior to sending to CMI

