



CERTIFICATION by CUSTODIAN of RECORDS

Tucson Police Department Crime Laboratory  
270 S. Stone Avenue Tucson, AZ 85701-1917

TSK 5.11.09

I hereby certify that these documents are a true and correct copy of the periodic maintenance\* for the Intoxilyzer 8000 Serial Number 80-002261 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.

A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.

Date: 5/6/09 Cal 1 0.099 AC Cal 2 0.099 AC

Between the dates of N/A and 5/6/09 the instrument was operating properly and accurately.

The instrument was taken out of service for repair on \_\_\_\_\_  
The instrument was returned to service on 5/6/09

Sarah Cione

Sarah L. Cione


5/7/09

ADPS QAS Permit # 33355

SUBSCRIBED AND SWORN to before me this 7<sup>th</sup> day of May, 2009.

A.M. Hernandez

Notary Public



\*Includes documents specified in A.A.C. R13-10-104.A.5

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.



TUCSON POLICE DEPARTMENT  
CRIME LABORATORY

1 of 12 sc

INTOXILYZER SERVICE RECORD

TJH  
5.11.09

QAS: Sarah Cione

Intoxilyzer Serial Number: 80-002261

Date Returned to Service: 5/6/09

2 of 12 sc

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY  
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

TSK  
5.11.09

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Sarah Cione AGENCY Tucson Police Department

DATE 5/6/09 TIME 13:39

INTOXILYZER SERIAL # 80-002261 LOCATION Crime Lab

(✓) 1. Display reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- (✓) 1. Clock time check.
- (✓) 2. Date check.

OPERATIONAL TESTS

- (✓) 1. Alcohol-free subject test result 0. 000 AC.
- (✓) 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed
- (✓) 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.
- (✓) 4. Standard Calibration Check standard 0. 100 AC  
Result: 0. 099 AC

Instrument is operating properly and accurately. YES  NO

COMMENTS  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SIGNATURE Sarah Cione

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

3 of 12 sc

INTOXILYZER 8000  
Location: TUCSON PD  
Serial Number : 80-002261  
05/06/2009 13:41:30

184

Standard Lot#: 808101G  
Last Changed By: S. CIONE #50616

QAS: S. CIONE #50616  
TUCSON PD

Operator: S. CIONE #50616  
TUCSON PD

Subject: SQAP, TEST  
DOB: 11/11/1111  
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	13:42:14
Diagnostic Test	Pass	13:42:41
Air Blank	0.000	13:43:09
0.100 Cal Check	0.099	13:43:28
Air Blank	0.000	13:43:57
Subject Test	0.000	13:44:16
Air Blank	0.000	13:44:43
Five Minute Wait	ABT*	
Air Blank	0.000	13:45:22

\*Sequence Aborted

**Not a Successfully Completed Test Sequence**

Pursuant to Section 28.1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department

4 of 12 sc

INTOXILYZER 8000  
Location: TUCSON PD  
Serial Number : 80-002261  
05/06/2009 13:45:43

BY

Standard Lot#: 808101G  
Last Changed By: S. CIONE #50616

QAS: S. CIONE #50616  
TUCSON PD

Operator: S. CIONE #50616  
TUCSON PD

Subject: SQAP, TEST  
DOB: 11/11/1111  
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	13:46:28
Air Blank	IPS**	13:46:38

\*RFI Detect  
\*\*Improper Sample

**Not a Successfully Completed Test Sequence**

Pursuant to Section 28.1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department

5 of 12 sc

INTOXILYZER 8000  
Location: TUCSON PD  
Serial Number : 80-002261  
05/06/2009 13:46:43

BM

Standard Lot#: 808101G  
Last Changed By: S. CIONE #50616

QAS: S. CIONE #50616  
TUCSON PD

Operator: S. CIONE #50616  
TUCSON PD

Subject: SQAP, TEST  
DOB: 11/11/1111  
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	13:47:33
Diagnostic Test	Pass	13:48:00
Air Blank	0.000	13:48:28
0.100 Cal Check	0.100	13:48:47
Air Blank	0.000	13:49:16
Subject Test	0.069*	13:52:24
Air Blank	0.000	13:52:56
Air Blank	0.000	13:53:23
Subject Test	XXX**	13:53:46
Air Blank	0.000	13:54:18

\*Deficient Sample  
\*\*Invalid Sample  
Begin new deprivation period

**Not a Successfully Completed Test Sequence**

Pursuant to Section 28.1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department

INTOXILYZER 8000  
Location: TUCSON PD  
Serial Number: 80-002261  
05/06/2009

13:54:34

6 of 12<sub>sc</sub>

-----  
Standard Lot#: 808101G  
Last Changed By: S. CIONE #50616  
-----

SM

-----  
QAS: S. CIONE #50616  
TUCSON PD  
-----

Test	g/210L	Time
Air Blank	0.000	13:55:01
Cal Check	0.099	13:55:20
Air Blank	0.000	13:55:49

*Sarah Cione*

Operator's Signature

Pursuant to Section 28.1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department

EXHIBIT G-2

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY  
A.A.C. R13-10-104(A)

7 of 12<sub>sc</sub>

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

DH  
5/11/09

QA Specialist Sarah Cione AGENCY Tucson Police Department  
DATE 5/6/09 TIME 13:39  
Intoxilyzer Serial # 80-002261 Location Crime Lab

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0. 100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 099 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Sarah Cione

DPS Form Exh G-2 (Rev. 05-01)

INTOXILYZER 8000  
Location: TUCSON PD  
Serial Number: 80-002261  
05/06/2009

13:39:40

8 of 12 SC

-----  
Standard Lot#: 808101G  
Last Changed By: S. CIONE #50616  
-----

124

QAS: S. CIONE #50616  
TUCSON PD  
-----

Test	g/210L	Time
Air Blank	0.000	13:40:07
Cal Check	0.099	13:40:26
Air Blank	0.000	13:40:56

*Sarah Cione*

Operator's Signature

Pursuant to Section 28.1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department



TUCSON POLICE DEPARTMENT  
CRIME LABORATORY

9 of 12 SC

INTOXILYZER SERVICE RECORD  
FOR  
CALIBRATION GAS TANK CHANGE

TSY  
5.11.09

Date: 5/4/09  
QAS: Sarah Cione  
Intoxilyzer Serial Number: 80-002261  
New Lot Number: 808101G  
New Expiration Date: 3/31/2010

- (✓) 1. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 2. Remove the old gas tank and install the new one.
- (✓) 3. Hit the "ESC" key twice and swipe the QAS card.
- (✓) 4. Go to the "M-Maintenance" menu and select "S-Standard Change."  
Select "D" for dry gas and enter the following information:  
  - Enter Target Value: 0.100
  - Enter Control Lot #: From the tank label
  - Enter Exp Date: From the tank label (mm/dd/yyyy)
  - Edit info Y/N? Y or N
  - Update Std Info? Y ENTER
  - "Please wait saving..."
- (✓) 5. Hit the "ESC" key twice.
- (✓) 6. Go to the "M-Maintenance" menu and select "T-Control Testing."  
Select "D" for dry control test. After selection is made press ENTER.
- (✓) 7. Air blank completed.
- (✓) 8. Calibration check completed. Test results 0. 100 AC.
- (✓) 9. Air blank completed.
- (✓) 10. Remove and sign printed record. Attach the record to the completed checklist.

SIGNATURE Sarah Cione

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

10 of 12sc

AS4

INTOXILYZER 8000  
Location: TUCSON PD  
Serial Number: 80-002261  
05/04/2009 14:25:08

-----  
Standard Lot#: 808101G  
Last Changed By: S. CIONE #50616

-----  
QAS: S. CIONE #50616  
TUCSON PD  
-----

Test	g/210L	Time
Air Blank	0.000	14:25:35
Cal Check	0.100	14:25:54
Air Blank	0.000	14:26:23

Sarah Cione  
Operator's Signature

Pursuant to Section 28.1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department



12 of 12 sc

Ph

# CERTIFICATE OF CALIBRATION

## INTOXILYZER® 8000

This is to certify that the calibration of INTOXILYZER serial number 80-002261, manufactured by CMI, subsidiary of MPD, Inc. of Owensboro, Kentucky, was tested and found to be in compliance with the National Safety Administration Standard for Devices to Measure Breath Alcohol (F.R., No. 179 48705-48710 Sept. 17, 1993) and with International Recommendation OIML R126 Edition 1998 (E). Calibration solutions are traceable to NIST (NBS) standard material 1828.

Date 4/17/09

Signed Kristi Robinson  
Technician

316 EAST 9TH STREET  
OWENSBORO, KY. 42303

**CMI** INC.

Pursuant to Section 28.1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department