

CERTIFICATION by CUSTODIAN of RECORDS

Tucson Police Department Crime Laboratory 270 S. Stone Avenue Tucson, AZ 85701-1917

154.09

I hereby certify that these documents are a true and correct copy of the periodic maintenance* for the Intoxilyzer 8000 Serial Number	
maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.	
A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2 for a 0.100AC solution are 0.090AC to 0.110AC.	
Date: 6/25/09 Cal 1 0.100 AC Cal 2 0.099 AC	•
Between the dates of and and	
The instrument was taken out of service for repair on The instrument was returned to service on	
Sanah Cione	
Sarah L. Cione	
6/25/09	
ADPS QAS Permit # 33355	
SUBSCRIBED AND SWORN to before me this 25 th , day of 20 09.	
June , 20 09. Anthernancles	• •

Notary Public

*Includes documents specified in A.A.C. R13-10-104.A.5

OFFICIAL SEAL
A.M. HERNÁNDEZ
NOTARY PUBLIC-ARIZONA
PIMA COUNTY
My Comm. Exp. Mar. 17, 2012

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

1047_{5c} 154_{0.09}

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Sarah Cione	AGENCY Tucson Police Department
DATE 6/25/09	
INTOXILYZER SERIAL # 80-002262	LOCATION DUT HQ
() 1. Display reads "PUSH BUTTON TO	START"
DIAGNOSTIC TESTS	
() 1. Clock time check. () 2. Date check.	
OPERATIONAL TESTS	
 () 1. Alcohol-free subject test res () 2. Error recognition logic systemote Not a Successfully Completed () 3. Proper sample recognition sysmote a Successfully Completed Deficient sample printed. () 4. Standard Calibration Check standard: 0.049 AC Instrument is operating properly and accessing the subject test res	m functioning. Test Sequence printed tem. Test Sequence printed andard 0. 100 AC
COMMENTS	
SIGNATURE Sanoh Cional	·

DPS Form Exh G-4 (Rev. 05-01)

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

INTOXILYZER 8000 Location: ARIZONA

Serial Number: 80-002262

06/25/2009

09:23:45

Standard Lot#: 826001I

Last Changed By: C. MORAND #24219

#50616 QAS: S. CIONE

TUCSON PD

Operator: S. CIONE #50616

TUCSON PD

Subject: SQAP, TEST DOB: 11/11/1111

Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:24:36
Diagnostic Test	Pass	09:25:04
Air Blank	0.000	09:25:31
0.100 Cal Check	0.101	09:25:50
Air Blank	0.000	09:26:19
Subject Test	0.000	09:26:52
Air Blank	0.000	09:27:19
Five Minute Wait	ABT*	
Air Blank	0.000	09:27:54

^{*}Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000 Location: ARIZONA

Serial Number: 80-002262

06/25/2009

09:28:55

Standard Lot#: 826001I

Last Changed By: C. MORAND #24219

QAS: S. CIONE #50616

TUCSON PD

#50616 Operator: S. CIONE

TUCSON PD

Subject: SQAP, TEST DOB: 11/11/1111

Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	09:29:46
Air Blank	IPS**	09:29:54

^{*}RFI Detect

Not a Successfully Completed Test Sequence

^{**}Improper Sample

INTOXILYZER 8000 Location: ARIZONA

Serial Number: 80-002262

09:30:39 06/25/2009

Standard Lot#: 826001I #24219 Last Changed By: C. MORAND #50616 OAS: S. CIONE TUCSON PD Operator: S. CIONE #50616 TUCSON PD

Subject: SQAP, TEST DOB: 11/11/1111

Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:31:40
Diagnostic Test	Pass	09:32:08
Air Blank	0.000	09:32:35
0.100 Cal Check	0.100	09:32:54
Air Blank	0.000	09:33:23
Subject Test	0.148*	09:36:32
Air Blank	0.000	09:37:04
Air Blank	0.000	09:37:31
Subject Test	XXX**	09:38:02
Air Blank	0.000	09:38:35

^{*}Deficient Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

^{**}Invalid Sample

5 of 7_{sc}

INTOXILYZER 8000 Location: ARIZONA

Serial Number: 80-002262

06/25/2009

09:38:53

Standard Lot#: 826001I

Last Changed By: C. MORAND #24219

QAS: S. CIONE #50616

TUCSON PD

Test	g/210L	Time
Air Blank	0.000	09:39:20
Cal Check	0.099	09:39:39
Air Blank	0.000	09:40:09



THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

154.09

STANDARD CALIBRATION CHECK PROCEDURE

QA Sp	ecial:	ist Sarah Cione AGENCY Tucson Police Department
DATE	6/2	5/09 TIME 9:21
Intox	ilyze:	r Serial # 80-002262 Location DUT HQ
(1)	1.	Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0. 100 AC.
		Pour a standard alcohol concentration solution AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} + 0.2^{\circ}\text{C}$
()	2.	Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
()	3.	Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
()	4.	Air blank completed.
(🗸)	5.	Calibration check completed. Test results 0. 100 AC.
	6.	Air blank completed.
()	7.	Remove printed record. Attach the record to the completed checklist.
SIGNA	TURE	Sanah Cione
DPS F	orm E:	xh G-2 (Rev. 05-01)

7 of 7sc

INTOXILYZER 8000 Location: ARIZONA

Serial Number: 80-002262

06/25/2009

09:21:48

Standard Lot#: 826001I

Last Changed By: C. MORAND #24219

QAS: S. CIONE #50616

TUCSON PD -

Test	g/210L	Time
Air Blank	0.000	09:22:15
Cal Check	0.100	09:22:34
Air Blank	0.000	09:23:04

