



# CERTIFICATION by CUSTODIAN of RECORDS

*Tucson Police Department Crime Laboratory*

I hereby certify that these documents are a true and correct copy of the periodic maintenance\* for the Intoxilyzer 8000 Serial Number 80-002264 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.

A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.

Date: 2.15.2007 Cal 1 0.101 AC Cal 2 0.101 AC

Between the dates of 1.17.2007 and 2.15.2007 the instrument was operating properly and accurately.


The instrument was taken out of service for repair on —.  
The instrument was returned to service on —.

*Terry S. Gallegos*  
Terry S. Gallegos

2.27.2007

ADPS QAS Permit # 33169

SUBSCRIBED AND SWORN to before me this 27, day of February, 2007.



*Kathy A. Kinney*  
Notary Public

\*Includes documents specified in A.A.C. R13-10-104.A.5

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

1 of 8  
fleTHIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY  
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

BM

QA SPECIALIST Jennifer Coombs 46660 AGENCY Tucson Police DepartmentDATE February 15, 2007 TIME 1228INTOXILYZER SERIAL # 80-002264 LOCATION Traffic

(✓) 1. Display reads "PUSH BUTTON TO START"

## DIAGNOSTIC TESTS

- (✓) 1. Clock time check.  
(✓) 2. Date check.

## OPERATIONAL TESTS

- (✓) 1. Alcohol-free subject test result 0.000 AC.  
(✓) 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
(✓) 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed ✓✓  
Deficient sample printed.  
(✓) 4. Standard Calibration Check standard 0.100 AC  
Result: 0.101 AC

Instrument is operating properly and accurately. YES ☒ NO ☐

COMMENTS + moved clock back 4 minutes 1235 → 1231 after  
improper sample check.

SIGNATURE

Jennifer L. Coombs

DPS Form Exh G-4 (Rev. 05-01)

2 of 8  
file

INTOXILYZER 8000  
Location: TUCSON POLICE DEPT  
Serial Number : 80-002264  
02/15/2007

12:30:00

ph

Standard Lot#: 629202I  
Last Changed By: T. GALLEGOS #28847

QAS: J. COOMBS #46160  
TUCSON PD CRIME LAB

Operator: J. COOMBS #46160  
TUCSON PD CRIME LAB

Subject: SQAP,TEST  
DOB: 11/11/1111  
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	12:30:49
Diagnostic Test	Pass	12:31:16
Air Blank	0.000	12:31:43
0.100 Cal Check	0.101	12:32:02
Air Blank	0.000	12:32:32
Subject Test	0.000	12:33:00
Air Blank	0.000	12:33:27
Five Minute Wait	ABT*	
Air Blank	0.000	12:34:13

\*Sequence Aborted

**Not a Successfully Completed Test Sequence**

3 of 8  
fle

INTOXILYZER 8000  
Location: TUCSON POLICE DEPT  
Serial Number : 80-002264  
02/15/2007

12:34:19

Bh

Standard Lot#: 629202I  
Last Changed By: T. GALLEGOS #28847

QAS: J. COOMBS #46160  
TUCSON PD CRIME LAB

Operator: J. COOMBS #46160  
TUCSON PD CRIME LAB

Subject: SQAP, TEST  
DOB: 11/11/1111  
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	12:35:08
Air Blank	0.000	12:35:36

\*Improper Sample

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number : 80-002264

02/15/2007

12:31:05

408  
jle

154

Standard Lot#: 629202I

Last Changed By: T. GALLEGOS

#28847

QAS: J. COOMBS

#46160

TUCSON PD CRIME LAB

Operator: J. COOMBS

#46160

TUCSON PD CRIME LAB

Subject: SQAP,TEST

DOB: 11/11/1111

Sex: F

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	12:32:02
Diagnostic Test	Pass	12:32:29
Air Blank	0.000	12:32:56
0.100 Cal Check	RFI*	12:33:13
Air Blank	0.000	12:33:43

\*RFI Detect

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: TUCSON POLICE DEPT  
Serial Number : 80-002264  
02/15/2007

12:33:55

BM

5 of 8  
file

Standard Lot#: 629202I  
Last Changed By: T. GALLEGOS #28847

QAS: J. COOMBS #46160  
TUCSON PD CRIME LAB

Operator: J. COOMBS #46160  
TUCSON PD CRIME LAB

Subject: SQAP, TEST  
DOB: 11/11/1111  
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	12:34:50
Diagnostic Test	Pass	12:35:17
Air Blank	0.000	12:35:44
0.100 Cal Check	0.101	12:36:03
Air Blank	0.000	12:36:33
Subject Test	0.075*	12:39:41
Air Blank	0.000	12:40:12
Air Blank	0.000	12:40:40
Subject Test	XXX**	12:41:01
Air Blank	0.000	12:41:34

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number: 80-002264

02/15/2007

12:42:05

-----  
Standard Lot#: 629202I

Last Changed By: T. GALLEGOS #28847  
-----

QAS: J. COOMBS #46160

TUCSON PD CRIME LAB  
-----

Test	g/210L	Time
Air Blank	0.000	12:42:33
Cal Check	0.101	12:42:51
Air Blank	0.000	12:43:21

  
Operator's Signature

EXHIBIT G-2

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY  
A.A.C. R13-10-104(A)

7048  
jlc

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

154

STANDARD CALIBRATION CHECK PROCEDURE

QA Specialist Jennifer Coombs 46160 AGENCY Tucson Police Department  
DATE February 15, 2007 TIME 1228  
Intoxilyzer Serial # 80-002264 Location Traffic

- ( ✓ ) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0. 100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- ( ✓ ) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- ( ✓ ) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- ( ✓ ) 4. Air blank completed.
- ( ✓ ) 5. Calibration check completed. Test results 0. 101 AC.
- ( ✓ ) 6. Air blank completed.
- ( ✓ ) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE

Jennifer L. Coombs

DPS Form Exh G-2 (Rev. 05-01)



INTOXILYZER 8000

Location: TUCSON POLICE DEPT

Serial Number: 80-002264

02/15/2007

12:28:16

*8 of 8*  
*jk*


-----  
Standard Lot#: 629202I

Last Changed By: T. GALLEGOS #28847  
-----

QAS: J. COOMBS #46160

TUCSON PD CRIME LAB  
-----

Test	g/210L	Time
Air Blank	0.000	12:28:43
Cal Check	0.101	12:29:02
Air Blank	0.000	12:29:32

  
\_\_\_\_\_  
Operator's Signature