

CERTIFICATION by CUSTODIAN of RECORDS Tucson Police Department Crime Laboratory

270 S. Stone Avenue Tucson, AZ 85701-1917

I hereby certify that these documents are a true and correct copy of the periodic maintenance* for the Intoxilyzer 8000 Serial Number 80-002345 maintained by the Tucson Police Department
Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.
A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.
Date: 07/01/10 Cal 1 0./01 AC Cal 2 0./00 AC
Between the dates of $06/03/10$ and $07/01/10$ the instrument was operating properly and accurately.
The instrument was taken out of service for repair on \(\nu/A\). The instrument was returned to service on \(\nu/A\).
Chyl K. Stu
Cheryl K. Straub
07/01/10 ADPS QAS Permit # 33408
SUBSCRIBED AND SWORN to before me this, day of, 20 /O.
Sul flernandes
Notary Public Notary Public *Includes documents specified in A A C. R13-10-104 A 5

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

My Comm. Exp. Mar. 17, 2012

*Includes documents specified in A.A.C. R13-10-104.A.5

EXHIBIT G-4

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

2.6.1

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

DATE 07/01/10	AGENCY Tucson Police Department
DATE 07/01/10	TIME 08:00
INTOXILYZER SERIAL # 80-002345	LOCATION TEAM 3
() 1. Display reads "PUSH BUTTON TO	O START"
DIAGNOSTIC TESTS	
() 1. Clock time check. () 2. Date check.	
OPERATIONAL TESTS	
1. Alcohol-free subject test res 2. Error recognition logic system Not a Successfully Completed 3. Proper sample recognition system Not a Successfully Completed Deficient sample printed. 4. Standard Calibration Check standard: Result: 0. /// AC	Test Sequence printed stem. Test Sequence printed
Instrument is operating properly and ac	ccurately. YESNO
COMMENTS	
SIGNATURE Chyl Sin	
•	

DPS Form Exh G-4 (Rev. 05-01)

INTOXILYZER 8000

Location: ARIZONA

Serial Number: 80-002345

07/01/2010

08:01:45

Standard Lot#: 642845

Last Changed By: C. STRAUB #50772

QAS: C. STRAUB #50772

TUCSON PD

Operator: C. STRAUB #50772

TUCSON PD

Subject: SQAP, SQAP

DOB: 11/11/1111

Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:02:38
Diagnostic Test	Pass	08:03:06
Air Blank	0.000	08:03:33
0.100 Cal Check	0.101	08:03:52
Air Blank	0.000	08:04:21
Subject Test	0.000	08:04:41
Air Blank	ABT*	08:04:48
Air Blank	0.000	08:05:16

^{*}Sequence Aborted

Not a Successfully Completed Test Sequence

134

INTOXILYZER 8000 Location: ARIZONA

Serial Number: 80-002345

07/01/2010

08:05:22

Standard Lot#: 642845

Last Changed By: C. STRAUB

#50772

QAS: C. STRAUB

#50772

TUCSON PD

Operator: C. STRAUB

#50772

TUCSON PD

Subject: SQAP, SQAP

DOB: 11/11/1111 Sex: F

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time	1
Air Blank	RFI*	08:05:55	
Air Blank	IPS**	08:06:04	

^{*}RFI Detect

Not a Successfully Completed Test Sequence

^{**}Improper Sample

INTOXILYZER 8000

Location: ARIZONA

Serial Number: 80-002345

07/01/2010

08:06:17

Standard Lot#: 642845

Last Changed By: C. STRAUB

#50772

QAS: C. STRAUB

#50772

TUCSON PD

Operator: C. STRAUB

#50772

TUCSON PD

Subject: SQAP, SQAP

DOB: 11/11/1111

Sex: F

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:07:13
Diagnostic Test	Pass	08:07:41
Air Blank	0.000	08:08:08
0.100 Cal Check	0.101	08:08:27
Air Blank	0.000	08:08:56
Subject Test	0.207*	08:12:05
Air Blank	0.000	08:12:37
Air Blank	0.000	08:13:05
Subject Test	XXX**	08:13:38
Air Blank	0.000	08:14:12

^{*}Deficient Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

^{**}Invalid Sample

INTOXILYZER 8000 Location: ARIZONA

Serial Number: 80-002345

07/01/2010

08:15:18

10

Standard Lot#: 642845

Last Changed By: C. STRAUB #50772

QAS: C. STRAUB #50772

TUCSON PD

	1	1.1
Test	g/210L	Time
Air Blank	0.000	08:15:45
Cal Check	0.100	08:16:04
Air Blank	0.000	08:16:34

Operator's Signature

EXHIBIT G-2

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

		SIMBIND CHIBRATION CHECK PROCEDURE
QA Spe	ciali	Ist Chepyl Strauß agency Tucson Police Department 101/10 TIME 08:00
DATE _	07	01/10 TIME 08:00
Intoxi	lyzer	Serial # <u>80-002345</u> Location <u>TEAM 3</u>
()	1.	Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0. $\cancel{100}$ AC.
	12	Pour a standard alcohol concentration solution NA AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
()	2.	Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
		Intoxilyzer 8000 display reads "PUSH BUTTON TO START" Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
()	4.	Air blank completed.
()	5.	Calibration check completed. Test results 0. /0/ AC.
(V)	6.	Air blank completed.
(V)	7.	Air blank completed. Calibration check completed. Test results 0. /0/ AC. Air blank completed. Remove printed record. Attach the record to the completed checklist.
SIGNAT		Chyl Fin

DPS Form Exh G-2 (Rev. 05-01)

INTOXILYZER 8000 Location: ARIZONA

Serial Number: 80-002345

07/01/2010

08:00:20

Standard Lot#: 642845

Last Changed By: C. STRAUB #50772

QAS: C. STRAUB #50772

TUCSON PD

Test	Ė	g/210L	Time
Air	Blank	0.000	08:00:47
Cal	Check	0.101	08:01:06
Air	Blank	0.000	08:01:35

Operator's Signature