

CERTIFICATION by CUSTODIAN of RECORDS
Tucson Police Department Crime Laboratory
270 S. Stone Avenue Tucson, AZ 85701-1917

I hereby certify that these documents are a true and correct copy of the periodic maintenance* for the Intoxilyzer 8000 Serial Number
periodic maintenance* for the Intoxilyzer 8000 Serial Number 80 - 002345 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department
of Public Safety for a quality assurance program.
A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.
Date: $6/23/11$ Cal 1 0.102 AC Cal 2 0.101 AC
Between the dates of $\frac{5/26/11}{26/11}$ and $\frac{6/23/11}{26/11}$ the instrument was operating properly and accurately.
The instrument was taken out of service for repair on
The instrument was returned to service on
Sarah Cione
Sarah L. Cione
6/27/11
ADPS QAS Permit # 33355
SUBSCRIBED AND SWORN to before me this 27, day of, day of, 20_//
AM Vernande &
OFFICIAL SEAL A M. HERNÁNDEZ
Notary Public Notary Public Notary Public NOTARY PUBLIC-ARIZONA PIMA COUNTY
*Includes documents specified in A.A.C. R13-10-104.A.5 My Comm. Exp. Mar. 17, 2012

Pursuant to section 28-1327. Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

EXHIBIT G-4

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

1077sc By 22.2011

QA SPECIALIST Sarah Cione	AGENCY Tucson Police Department	
DATE 6/23/1	TIME <u>1310</u>	
INTOXILYZER SERIAL # 80-002345		
(\checkmark) 1. Display reads "PUSH BUTTON TO		
DIAGNOSTIC TESTS		
() 1. Clock time check. () 2. Date check.		
OPERATIONAL TESTS		
<pre>(/) 1. Alcohol-free subject test result 0.000AC. (/) 2. Error recognition logic system functioning.</pre>		
COMMENTS		
signature Sarah Con		

DPS Form Exh G-4 (Rev. 05-01)

2 of 7₅₀

INTOXILYZER 8000 Location: ARIZONA

Serial Number: 80-002345

06/23/2011

13:14:37

Standard Lot#: 20910100A1 Last Changed By: C. STRAUB #50772 QAS: S. CIONE #50616 TUCSON PD Operator: S. CIONE #50616 TUCSON PD

Subject: SQAP, SQAP DOB: 11/11/1111

Sex: F

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	13:15:22
Diagnostic Test	Pass	13:15:49
Air Blank	0.000	13:16:16
0.100 Cal Check	0.102	13:16:35
Air Blank	0.000	13:17:05
Subject Test	0.000	13:17:49
Air Blank	0.000	13:18:16
Five Minute Wait	ABT*	
Air Blank	0.000	13:19:07

^{*}Sequence Aborted

Not a Successfully Completed Test Sequence

3 of 7_{sc}

134

INTOXILYZER 8000 Location: ARIZONA

Serial Number: 80-002345

06/23/2011 13:19:22

Standard Lot#: 20910100A1
Last Changed By: C. STRAUB #50772

QAS: S. CIONE #50616
TUCSON PD

Operator: S. CIONE #50616
TUCSON PD

Subject: SQAP, SQAP DOB: 11/11/1111

Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	13:19:58
Air Blank	IPS**	13:20:05

^{*}RFI Detect

Not a Successfully Completed Test Sequence

^{**}Improper Sample

INTOXILYZER 8000

Location: ARIZONA

Serial Number: 80-002345

06/23/2011

13:20:38

Standard Lot#: 20910100A1

Last Changed By: C. STRAUB

#50772

QAS: S. CIONE

#50616

TUCSON PD

Operator: S. CIONE

#50616

TUCSON PD

Subject: SQAP, SQAP

DOB: 11/11/1111

Sex: F

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	13:21:27
Diagnostic Test	Pass	13:21:55
Air Blank	0.000	13:22:22
0.100 Cal Check	0.102	13:22:41
Air Blank	0.000	13:23:10
Subject Test	0.042*	13:26:18
Air Blank	0.000	13:26:50
Air Blank	0.000	13:27:17
Subject Test	XXX**	13:27:51
Air Blank	0.000	13:28:24

^{*}Deficient Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

4 of 750

^{**}Invalid Sample

INTOXILYZER 8000

Location: ARIZONA

Serial Number: 80-002345

06/23/2011

13:29:52

Standard Lot#: 20910100A1

Last Changed By: C. STRAUB #50772

QAS: S. CIONE #50616

TUCSON PD

g/210L Test Time Air Blank 0.000 0.101 Cal Check 13:30:38

Air Blank 0.000 13:31:07

5 of 7_{sc}

EXHIBIT G-2

6 of 75c

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

6.27.2011

QA Speciali	list <u>Darah Lion</u> AGENCY Tucsor	n Police Department
DATE 6/2	23/11 TIME 1310	
	er Serial # <u>80-002345</u> Location <u>Crin</u>	re Lab
*		. *
1.	Ensure that gas tank is attached to instract standard alcohol concentration solution 0.	rument and contains a 100 AC.
	Pour a standard alcohol concentration into a clean dry simulator and assemble that a tight seal has been made. Turn allow temperature to reach 34°C \pm 0.2°C	the simulator. Ensure
	Intoxilyzer 8000 display reads "PUSH	
(pV) 3.	Go to the "Control Testing Menu". So control test or "W" for wet con selection is made press ENTER.	Select "D" for dry trol test. After
(✓ ,) 4.	Air blank completed.	
(🗸) 5.	Calibration check completed. Test re	sults 0. <u>102</u> AC.
(/) 6.	Air blank completed.	
() 7.	Air blank completed. Remove printed record. Attach the record.	ord to the completed
SIGNATURE	SanahCion	

DPS Form Exh G-2 (Rev. 05-01)

INTOXILYZER 8000 Location: ARIZONA

Serial Number: 80-002345

06/23/2011

Standard Lot#: 20910100A1

Last Changed By: C. STRAUB #50772

QAS: S. CIONE #50616

TUCSON PD

Test		g/210L	Time
Air I	Blank	0.000	13:10:45
Cal	Check	0.102	13:11:04
Air I	Blank	0.000	13:11:34

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.