



CERTIFICATION by CUSTODIAN of RECORDS

Tucson Police Department Crime Laboratory  
270 S Stone Avenue Tucson, AZ 85701-1917

I hereby certify that these documents are a true and correct copy of the periodic maintenance\* for the Intoxilyzer 8000 Serial Number 80-002759 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.

TS4  
9.19.2011

A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.

Date: 08/25/11 Cal 1 0.098 AC Cal 2 0.097 AC

Between the dates of 07/28/11 and 08/25/11 the instrument was operating properly and accurately.

The instrument was taken out of service for repair on N/A  
The instrument was returned to service on N/A

Cheryl K. Straub

Cheryl K. Straub


09/08/2011

ADPS QAS Permit # 33408

SUBSCRIBED AND SWORN to before me this 8 day of Sept, 20 11.

Ann L. Bieczcroft

Notary Public



\*Includes documents specified in A.A.C. R13-10-104.A.5

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

EXHIBIT G-4

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY  
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

By  
9.19.2011

QA SPECIALIST Cheryl STRAUB AGENCY Tucson Police Department

DATE 08/25/2011 TIME 11:02

INTOXILYZER SERIAL # 80-002759 LOCATION TEAM 4

(✓) 1. Display reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- (✓) 1. Clock time check.  
(✓) 2. Date check.

OPERATIONAL TESTS

- (✓) 1. Alcohol-free subject test result 0. 000 AC.  
(✓) 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
(✓) 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.  
(✓) 4. Standard Calibration Check standard 0. 100 AC  
Result: 0. 097 AC

Instrument is operating properly and accurately. YES ✓ NO     

COMMENTS

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SIGNATURE

Cheryl Straub

277  
C8

INTOXILYZER 8000

Location:

Serial Number : 80-002759

08/25/2011

11:11:05

Standard Lot#: 04911100A1

Last Changed By: C. STRAUB

#50772

QAS: C. STRAUB

#50772

TUCSON PD

Operator: C. STRAUB

#50772

TUCSON PD

Subject: SQAP, SQAP

DOB: 11/11/1111

Sex: F

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	11:11:56
Diagnostic Test	Pass	11:12:23
Air Blank	0.000	11:12:51
0.100 Cal Check	0.098	11:13:09
Air Blank	0.000	11:13:39
Subject Test	0.000	11:13:59
Air Blank	ABT*	11:14:06
Air Blank	0.000	11:14:34

\*Sequence Aborted

### Not a Successfully Completed Test Sequence

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

377  
'8

INTOXILYZER 8000

Location:

Serial Number : 80-002759

08/25/2011

11:14:45

BY

Standard Lot#: 04911100A1

Last Changed By: C. STRAUB

#50772

QAS: C. STRAUB

#50772

TUCSON PD

Operator: C. STRAUB

#50772

TUCSON PD

Subject: SQAP, SQAP

DOB: 11/11/1111

Sex: F

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	11:15:15
Air Blank	IPS**	11:15:25

\*RFI Detect

\*\*Improper Sample

**Not a Successfully Completed Test Sequence**

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

487  
CS

INTOXILYZER 8000

Location:

Serial Number : 80-002759

08/25/2011

11:15:31

134

Standard Lot#: 04911100A1

Last Changed By: C. STRAUB

#50772

QAS: C. STRAUB

#50772

TUCSON PD

Operator: C. STRAUB

#50772

TUCSON PD

Subject: SQAP, SQAP

DOB: 11/11/1111

Sex: F

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	11:16:19
Diagnostic Test	Pass	11:16:47
Air Blank	0.000	11:17:14
0.100 Cal Check	0.098	11:17:33
Air Blank	0.000	11:18:03
Subject Test	0.315*	11:21:11
Air Blank	0.000	11:21:46
Air Blank	0.000	11:22:13
Subject Test	XXX**	11:22:37
Air Blank	0.000	11:23:11

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

**Not a Successfully Completed Test Sequence**

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

527  
08

INTOXILYZER 8000

Location:

Serial Number: 80-002759

08/25/2011

11:23:46


-----  
Standard Lot#: 04911100A1

Last Changed By: C. STRAUB #50772  
-----

QAS: C. STRAUB #50772

TUCSON PD  
-----

Test	g/210L	Time
Air Blank	0.000	11:24:14
Cal Check	0.097	11:24:32
Air Blank	0.000	11:25:02



Operator's Signature

677  
CS

EXHIBIT G-2

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY  
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

BS  
9.19.2011

QA Specialist Cheryl STRAUB AGENCY Tucson Police Department  
DATE 08/25/2011 TIME 11:02  
Intoxilyzer Serial # 80-002759 Location TEAM 4

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0. 100 AC.  
OR  
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 098 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE

Cheryl Straub

DPS Form Exh G-2 (Rev. 05-01)

787  
08

134

INTOXILYZER 8000  
Location:  
Serial Number: 80-002759  
08/25/2011 11:02:38

Standard Lot#: 04911100A1  
Last Changed By: C. STRAUB #50772

OAS: C. STRAUB #50772  
TUCSON PD

Test	g/210L	Time
Air Blank	0.000	11:03:05
Cal Check	0.098	11:03:24
Air Blank	0.000	11:03:54

  
Operator's Signature