



WLL
2/11/16

CERTIFICATION by CUSTODIAN of RECORDS

Tucson Police Department Crime Laboratory
1306 W. Miracle Mile, Tucson, AZ 85705

I hereby certify that these documents are a true and correct copy of the periodic maintenance* for the Intoxilyzer 8000 Serial Number 80-006054 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.

A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.

Date: 1-28-2016 Cal 1 — AC Cal 2 — AC

Between the dates of — and — the instrument was operating properly and accurately.

The instrument was taken out of service for repair on 1-28-2016.
The instrument was returned to service on —.

Terry S. Gallegos

Terry S. Gallegos

2-9-2016

ADPS QAS Permit # 33169

SUBSCRIBED AND SWORN to before me this 9th, day of February, 2016.

Elma Escarcega

Notary Public

*Includes documents specified in A.A.C. R13-10-104.A.5




EXHIBIT G-2

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA Specialist Terry Callegos 28847 AGENCY Tucson Police Department
DATE 1-28-2016 TIME —
Intoxilyzer Serial # 80-006054 Location DUI-HQ

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0. 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "~~PUSH BUTTON TO START~~" "Fault #018 Contact CMI" TS4
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. _____ AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE

Terry Callegos

DPS Form Exh G-2 (Rev. 05-01)

No calibration check completed.
TS4



TUCSON POLICE DEPARTMENT
CRIME LABORATORY

INTOXILYZER SERVICE RECORD

2 of 3
BX
WLL

Date: 1-28-2016

QAS: Terry Gallegos 28847

Intoxilyzer Serial Number: 80-006054

Reason For Service: "Fault #018, Contact CMI" on display.
and failed these Diagnostics: EEPROM, DSP,
Analytical Test.

Date Out of Service: 1-28-2016

Calibration Check completed prior to removal from service

Date: No

SQAP completed prior to removal from service

Date: No

Corrective Action: ☒ Sent to CMI, Inc.

Records downloaded

30f3
TS4
WLL

EMI INC.
Service Evaluation Form

This form MUST be completed and enclosed with instrument to be serviced.

Note: Please ship items in their original shipping container.

Contact information: Customer Number 105075 (contact Customer Service)
Name TERRY GALLEGOS Phone: (520) 837-7677
Fax: (520) 791-4642 Email: TERRY.GALLEGOS@TUCSONAZ.GOV

Bill to Address:	Ship to Address:
<u>FINANCE/ACCOUNTS PAYABLE</u>	<u>TUCSON POLICE DEPARTMENT CRIMELAB</u>
<u>PO BOX 27450</u>	<u>1306 W. MIRAGE MILE</u>
<u>TUCSON, AZ 85726-7450</u>	<u>TUCSON, AZ 85705</u>

Instrument Serial Number: 80-006054

Detailed Description of Problem:


Display "System Fault #018 Contact CMT" and it failed
Diagnostics: EEPROM, DSP, Analytical Test.

***Hazardous Material Warning! - DO NOT return gas cylinders with instrument! ***

Note: An evaluation fee of \$39.50 for handheld instruments or \$79.00 for IR evidential instruments will apply to estimates which are not repaired.

☐ I Authorize Repairs Up To: ☐ All ☐ \$250 ☐ \$500 ☐ \$750 ☐ Other \$ _____
Purchase Order Number 37768

Authorized By:

<u>TERRY GALLEGOS TOXICOLOGY COORDINATOR</u>	Ship to:
Name (Please Print)	Title
<u></u>	<u>1-28-2016</u>
Signature	Date

CMI, Inc.
Attn: Service Dept.
316 East Ninth Street
Owensboro, KY 42303

☒ No, please send estimate before repairs are made. An estimate will be faxed before performing any repairs and may cause delays in service.

316 EAST NINTH STREET • OWENSBORO, KENTUCKY 42303 • 1-866-835-0690

02/12