



CERTIFICATION by CUSTODIAN of RECORDS  
Tucson Police Department Crime Laboratory

154

I hereby certify that these documents are a true and correct copy of the periodic maintenance\* for the Intoxilyzer 5000 Serial Number 666-004814 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.

A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R9-14-404.A.3, for a 0.080AC solution are 0.070AC to 0.090AC.

Date: 1/20/05 Cal 1 0.075 AC Cal 2 0.075 AC

Date: 2/4/05 Cal 1 N/A AC Cal 2 N/A AC

Date: N/A Cal 1                      AC Cal 2                      AC

Between the dates of N/A and N/A the instrument was working properly and accurately.

The instrument was taken out of service for repair on 2/4/05.


The instrument was returned to service on 3/18/05.

Jennifer L. Coombs  
Jennifer L. Coombs

3/23/05

ADPS QAS Permit # QAS 29886

SUBSCRIBED AND SWORN to before me this 23, day of March, 2005.



Kathy A. Kinney  
Notary Public

\*Includes documents specified in A.A.C. R9-14-404.A.6

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.




# MEMORANDUM

---

**DATE:** February 9, 2005

**TO:** Whom it may concern

**FROM:** Jennifer Coombs  
Intoxilyzer 5000 QAS  
Tucson Police Department  
Crime Lab  
270 S. Stone Avenue  
Tucson, AZ 85701-1917  
(520) 791-4494



**SUBJECT:** Intoxilyzer Out of Service

Intoxilyzer 5000, serial number **66-004814**, was taken out of service for repair on **February 4, 2005**.

I cannot testify that it was working properly and accurately between the dates of **January 20, 2005** and **February 4, 2005**.

This instrument was located at **ODD**.



TUCSON POLICE DEPARTMENT  
CRIME LABORATORY

INTOXILYZER SERVICE RECORD

Date: February 4, 2005

QAS: Jennifer Coombs 46160

Intoxilyzer Serial Number: 666-004814

Reason For Service: On 2/4/05, Ofc. Nelson left a message stating that the instrument displayed "Out of Range".

Date Out of Service: 2/4/05

Calibration Check completed prior to removal from service

Date: No

SQAP completed prior to removal from service

Date: No

Corrective Action:

☒ Sent to Applied Electronics

Instrument brought to lab. During diagnostic test, instrument displayed "Processor Error 5". Plugged pump into fan connector and ran for 1.25 hours. Diagnostic checks still produced "Processor Error 5".