



CERTIFICATION by CUSTODIAN of RECORDS
Tucson Police Department Crime Laboratory

BR

I hereby certify that these documents are a true and correct copy of the periodic maintenance* for the Intoxilyzer 5000 Serial Number 166-004814 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.

A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R9-14-404.A.3, for a 0.080AC solution are 0.070AC to 0.090AC.


Date: 9/1/05 Cal 1 0.077 AC Cal 2 0.077 AC
 Date: 9/23/05 Cal 1 N/A AC Cal 2 N/A AC
 Date: N/A Cal 1 AC Cal 2 AC

Between the dates of N/A and N/A the instrument was operating properly and accurately.

The instrument was taken out of service for repair on 9/23/05.
 The instrument was returned to service on N/A.

Jennifer L. Coombs
 Jennifer L. Coombs
10/4/05

ADPS QAS Permit # QAS 29886
 SUBSCRIBED AND SWORN to before me this 4, day of October, 2005.



Kathy A. Kinney
 Notary Public

*Includes documents specified in A.A.C. R9-14-404.A.6



MEMORANDUM

DATE: September 23, 2005

TO: Whom it may concern

FROM: Jennifer L. Coombs
Intoxilyzer 5000 QAS
Tucson Police Department
Crime Lab
270 S. Stone Avenue
Tucson, AZ 85701-1917
(520) 791-4494

SUBJECT: Intoxilyzer Out of Service

Intoxilyzer 5000, serial number **66-004814**, was taken out of service for repair on **September 23, 2005**.

I cannot testify that it was working properly and accurately between the dates of **September 1, 2005** and **September 23, 2005**.

This instrument was located at **ODM**.



TUCSON POLICE DEPARTMENT
CRIME LABORATORY

INTOXILYZER SERVICE RECORD

DH

Date: September 23, 2005

QAS: Jennifer Coombs

Intoxilyzer Serial Number: 666-004814 @ODM

Reason For Service: Ofc. Bonnie (26131) reported the instrument would not get passed the air blank. Upon testing, the instrument displayed "unstable reference". I then brought it back to the lab where during the diagnostic check a Processor Error 5 occurred.

Date Out of Service: 9/23/05

Calibration Check completed prior to removal from service

Date: N/A

SQAP completed prior to removal from service

Date: N/A

Corrective Action:

Sent to Applied Electronics

