

CERTIFICATION by CUSTODIAN of RECORDS

Tucson Police Department Crime Laboratory 1306 W. Miracle Mile, Tucson, AZ 85705

SC

I hereby certify that these documents are a true and correct copy of the periodic maintenance* for the Intoxilyzer 8000 Serial Number 80-000334 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.	1/18/1)
A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.	
Date: 11-10-11 Cal 1 0.099 AC Cal 2 0.098 AC	
Between the dates of 10-17-11 and 1)-10-11 the instrument was operating properly and accurately.	-
The instrument was taken out of service for repair on N/A. The instrument was returned to service on N/A.	
Clare Conly	
Clare L. Conley	
11-10-11	
ADPS QAS Permit # 33314	
SUBSCRIBED AND SWORN to before me this, day of	
Notary Public Notary Public OFFICIAL SE DIANE QUII Notary Public Notary Public	ROZ RIZONA NTY

^{*}Includes documents specified in A.A.C. R13-10-104.A.5

EXHIBIT G-4

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

th?

<

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Mare Conley	AGENCY Tucson Police Department
DATE	TIME0859
INTOXILYZER SERIAL # 80-000334	
(/) 1. Display reads "PUSH BUTTON TO	START"
DIAGNOSTIC TESTS	
<pre>(/) 1. Clock time check. (/) 2. Date check.</pre>	
OPERATIONAL TESTS	
 1. Alcohol-free subject test res 2. Error recognition logic syste Not a Successfully Completed 3. Proper sample recognition sys Not a Successfully Completed Deficient sample printed. 4. Standard Calibration Check st Result: 0. 098 AC Instrument is operating properly and access	m functioning. Test Sequence printed tem. Test Sequence printed andard 0. 100 AC
COMMENTS None.	
SIGNATURE Clare Con ley	

DPS Form Exh G-4 (Rev. 05-01)

Location: TUCSON POLICE DEPT

Serial Number: 80-000334

11/10/2011

09:00:58

Standard Lot#: 20910100A1 Last Changed By: S. CIONE #50616 #34517 QAS: C. CONLEY

TUCSON PD

Operator: C. CONLEY #34517

TUCSON PD

Subject: TEST, SQAP DOB: 11/11/1111

Sex: M

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:01:43
Diagnostic Test	Pass	09:02:11
Air Blank	0.000	09:02:38
0.100 Cal Check	0.098	09:02:57
Air Blank	0.000	09:03:26
Subject Test	0.000	09:04:03
Air Blank	ABT*	09:04:10
Air Blank	0.000	09:04:38

^{*}Sequence Aborted

Not a Successfully Completed Test Sequence

Location: TUCSON POLICE DEPT

Serial Number: 80-000334

11/10/2011

09:04:44

Standard Lot#: 20910100A1
Last Changed By: S. CIONE #50616

QAS: C. CONLEY #34517
TUCSON PD

Operator: C. CONLEY #34517
TUCSON PD

Subject: TEST, SQAP DOB: 11/11/1111

Sex: M

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	09:05:15
Air Blank	RFI**	09:05:23

^{*}Improper Sample

Not a Successfully Completed Test Sequence

) / **S(**

^{**}RFI Detect

Location: TUCSON POLICE DEPT

Serial Number: 80-000334

11/10/2011

09:05:44

Standard Lot#: 20910100A1 Last Changed By: S. CIONE #50616 QAS: C. CONLEY #34517 TUCSON PD Operator: C. CONLEY #34517 TUCSON PD

Subject: TEST, SQAP

DOB: 11/11/1111

Sex: M

Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:06:37
Diagnostic Test	Pass	09:07:05
Air Blank	0.000	09:07:32
0.100 Cal Check	0.099	09:07:51
Air Blank	0.000	09:08:20
Subject Test	0.090*	09:11:29
Air Blank	0.000	09:12:00
Air Blank	0.000	09:12:28
Subject Test	XXX**	09:12:58
Air Blank	0.000	09:13:31

^{*}Deficient Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

^{**}Invalid Sample

Location: TUCSON POLICE DEPT

Serial Number: 80-000334

11/10/2011

09:13:48

Standard Lot#: 20910100A1

Last Changed By: S. CIONE #50616

QAS: C. CONLEY #34517

TUCSON PD

Test	g/210L	Time
Air Blank	0.000	09:14:15
Cal Check	0.098	09:14:34
Air Blank	0.000	09:15:04

Operator's Signature

SC

EXHIBIT G-2

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

6.67 Ch

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

PEL 8000

QA Specialist Mare Conley AGENCY Tucson Police Department DATE 1-10-11 Intoxilyzer Serial # 80-000334 Location ODW Dut HQ 1. Ensure that gas tank is attached to instrument and contains a

- 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0. AC.

 OR

 Pour a standard alcohol concentration solution AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C + 0.2°C
- (\checkmark) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (\checkmark) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (\checkmark) 4. Air blank completed.
- (\checkmark) 5. Calibration check completed. Test results 0. 099 AC.
- (/) 6. Air blank completed.
- (\(\subseternight) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Mare Conley

DPS Form Exh G-2 (Rev. 05-01)

Location: TUCSON POLICE DEPT

Serial Number: 80-000334 11/10/2011

08:59:32

Last Changed By: S. CIONE #50616

QAS: C. CONLEY #34517

TUCSON PD

Test	g/210L	Time
Air Blank	0.000	08:59:59
Cal Check	0.099	09:00:18
Air Blank	0.000	09:00:47

Operator's Signature

7667 en

