



CERTIFICATION by CUSTODIAN of RECORDS  
Tucson Police Department Crime Laboratory

By  
5-12-08

I hereby certify that these documents are a true and correct copy of the periodic maintenance\* for the Intoxilyzer 8000 Serial Number 80-001163 maintained by the Tucson Police Department Crime Laboratory pursuant to the requirements of the Arizona Department of Public Safety for a quality assurance program.

A quality assurance specialist (QAS) tested the instrument on the following dates. The acceptable accuracy limits, as defined in A.A.C. R13-10-104.A.2, for a 0.100AC solution are 0.090AC to 0.110AC.

Date: 4/21/08 Cal 1 N/A AC Cal 2 N/A AC

Between the dates of N/A and N/A the instrument was operating properly and accurately.

The instrument was taken out of service for repair on 4/21/08.  
The instrument was returned to service on -

Sarah Cione


Sarah L. Cione

5/12/08

ADPS QAS Permit # 33355

SUBSCRIBED AND SWORN to before me this 12, day of May, 2008.

Kathy A Kinney  
Notary Public



\*Includes documents specified in A.A.C. R13-10-104.A.5

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.



TUCSON POLICE DEPARTMENT  
CRIME LABORATORY

1 of 25c

INTOXILYZER SERVICE RECORD

TS4  
512-08

Date: 4/21/08

QAS: Sarah Cione

Intoxilyzer Serial Number: 80-001163

Reason For Service: Officer Ohara notified crime lab that instrument would not go to Ready Mode. Breath hose would heat up, but temperature is not registering. Breath hose gets VERY warm. Regulator has slow leak.

Date Out of Service: 4/21/08

Calibration Check completed prior to removal from service

Date: NO

SQAP completed prior to removal from service

Date: NO

Corrective Action:

Sent to CMI, Inc. (4/22/08)

COBRA data downloaded before sending to CMI



2 of 25c

# Return Material Authorization Form

This form **MUST** be completed and enclosed with item(s) being returned.  
Failure to complete and return this form may cause delays in crediting / repairing.

1. Item(s) Returned:  Instrument  Supplies  Other  
 (Note: please ship items in their original shipping container(s) or a similar protective box.)

2. Instrument Model: Intox 8000 Serial Number 80-001163

3. Bill To Address: Finance / Accounts Payable Ship To Address: (NO PO BOXES)  
P.O. Box 27450 Tucson Police Department  
TUCSON, AZ 85726-7450 Crime Lab  
270 S. Stone Ave  
TUCSON, AZ 85701-1917

Customer # 105075  
 Name Terry Gallegos  
 Phone: (520) 791-4494 Fax: (520) 791-4642  
 Email: Terry.Gallegos@tucsonaz.gov

4. Reason for Return: Note - If you are sending an item for repair, please give a detailed description of the problem. Please list any special instructions that you may have concerning this return. Note - Returns for credit may be subject to a 15% restocking charge.  
Breath Hose heats up, but temperature is not registering. Instrument never goes to "Ready Mode". Breath hose gets VERY warm. Also regulator has a slow leak (~200 psi/day).

5. Purchase Order Number (attach a copy of P.O. if applicable) 20669

6. I Authorize All Repairs:  Yes  No

7. I Authorize Repairs Up To: \$ \_\_\_\_\_

8. I Require An Estimate Regardless of Cost  Yes  No

**Note: An evaluation fee will apply to estimates that are not repaired.**

9. Authorized By:  
Terry Gallegos  
 Name (Please Print)  
[Signature]  
 Signature Date 4-22-08

Pursuant to Section 28.1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.