

189
CS

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD QUALITY ASSURANCE PROCEDURE

TBY
6.9.2009

QA SPECIALIST Cheryl STRAUB AGENCY Tucson Police Department
DATE 06/04/09 TIME 08:08
INTOXILYZER SERIAL # 80-002261 LOCATION DUE HA

() 1. Display reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- () 1. Clock time check.
- () 2. Date check.

OPERATIONAL TESTS

- () 1. Alcohol-free subject test result 0. 000 AC.
- () 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
- () 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
- () 4. Standard Calibration Check standard 0. 100 AC
Result: 0. 100 AC

Instrument is operating properly and accurately. YES NO

COMMENTS RESIDUAL mouth alcohol when performing
Alcohol free subject Test. Test repeated

SIGNATURE Cheryl Straub

289
CS

BY

INTOXILYZER 8000
Location: TUCSON PD
Serial Number : 80-002261
06/04/2009

08:11:08

Standard Lot#: 808101G
Last Changed By: S. CIONE #50616

QAS: C. STRAUB #50772
TUCSON PD

Operator: C. STRAUB #50772
TUCSON PD

Subject: TEST,TEST
DOB: 11/11/1111
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:12:05
Diagnostic Test	Pass	08:12:32
Air Blank	0.000	08:13:00
0.100 Cal Check	0.100	08:13:19
Air Blank	0.000	08:13:48
Subject Test	0.014	08:14:05
Air Blank	0.000	08:14:35
Five Minute Wait	ABT*	
Air Blank	0.000	08:17:00

*Sequence Aborted

- Residual mouth Alcohol -
eg. 50772

Not a Successfully Completed Test Sequence

298
CS
339
CS

B4

INTOXILYZER. 8000
Location: TUCSON PD
Serial Number : 80-002261
06/04/2009

08:18:11

Standard Lot#: 808101G
Last Changed By: S. CIONE #50616

QAS: C. STRAUB #50772
TUCSON PD

Operator: C. STRAUB #50772
TUCSON PD

Subject: TEST,TEST
DOB: 11/11/1111
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:19:10
Diagnostic Test	Pass	08:19:37
Air Blank	0.000	08:20:05
0.100 Cal Check	0.100	08:20:24
Air Blank	0.000	08:20:53
Subject Test	0.000	08:21:23
Air Blank	0.000	08:21:50
Five Minute Wait	ABT*	
Air Blank	0.000	08:22:40

*Sequence Aborted

Not a Successfully Completed Test Sequence

388
C8
489
C8

154

INTOXILYZER 8000
Location: TUCSON PD
Serial Number : 80-002261
06/04/2009 08:22:46

Standard Lot#: 808101G
Last Changed By: S. CIONE #50616

QAS: C. STRAUB #50772
TUCSON PD

Operator: C. STRAUB #50772
TUCSON PD

Subject: TEST,TEST
DOB: 11/11/1111
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	08:23:36
Air Blank	IPS**	08:23:47

*RFI Detect
**Improper Sample

Not a Successfully Completed Test Sequence

478
C8
579
C8

TS4

INTOXILYZER 8000
Location: TUCSON PD
Serial Number : 80-002261
06/04/2009 08:24:42

Standard Lot#: 808101G
Last Changed By: S. CIONE #50616

QAS: C. STRAUB #50772
TUCSON PD

Operator: C. STRAUB #50772
TUCSON PD

Subject: TEST,TEST
DOB: 11/11/1111
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:25:38
Diagnostic Test	Pass	08:26:05
Air Blank	0.000	08:26:32
0.100 Cal Check	0.100	08:26:51
Air Blank	0.000	08:27:20
Subject Test	XXX*	08:28:54
Air Blank	0.000	08:29:27

*Invalid Sample
Begin new deprivation period

Not a Successfully Completed Test Sequence

588
C8
689
C8

B4

INTOXILYZER 8000
Location: TUCSON PD
Serial Number : 80-002261
06/04/2009

08:31:22

Standard Lot#: 808101G
Last Changed By: S. CIONE #50616

QAS: C. STRAUB #50772
TUCSON PD

Operator: C. STRAUB #50772
TUCSON PD

Subject: TEST,TEST
DOB: 11/11/1111
Sex: F Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:32:25
Diagnostic Test	Pass	08:32:52
Air Blank	0.000	08:33:20
0.100 Cal Check	0.100	08:33:38
Air Blank	0.000	08:34:08
Subject Test	0.014*	08:37:16
Air Blank	0.000	08:37:47
Air Blank	0.000	08:38:14
Subject Test	DEF**	08:41:22
Air Blank	0.000	08:41:51

*Deficient Sample
**Deficient Sample

Not a Successfully Completed Test Sequence

698
CS
739
CS

INTOXILYZER 8000
Location: TUCSON PD
Serial Number: 80-002261
06/04/2009 08:42:55

TS4

Standard Lot#: 808101G
Last Changed By: S. CIONE #50616

QAS: C. STRAUB #50772
TUCSON PD

Test	g/210L	Time
Air Blank	0.000	08:43:22
Cal Check	0.100	08:43:41
Air Blank	0.000	08:44:10



Operator's Signature

EXHIBIT G-2

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY
A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

TSY
6-9-2009

798
CS
879
CS

QA Specialist Cheryl STRAUB AGENCY Tucson Police Department
DATE 06/04/09 TIME 08:08
Intoxilyzer Serial # 80-002261 Location DUF HQ

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0. 100 AC.
OR
Pour a standard alcohol concentration solution AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C + 0.2°C
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 100 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Cheryl Straub

DPS Form Exh G-2 (Rev. 05-01)

Pursuant to section 28-1327, Arizona Revised Statutes, this document is a certified duplicate of the information contained in computer storage devices of the Tucson Police Department.

~~898~~
CS
989
CS

INTOXILYZER 8000
Location: TUCSON PD
Serial Number: 80-002261
06/04/2009 08:08:41

134

Standard Lot#: 808101G
Last Changed By: S. CIONE #50616

QAS: C. STRAUB #50772
TUCSON PD

Test	g/210L	Time
Air Blank	0.000	08:09:09
Cal Check	0.100	08:09:28
Air Blank	0.000	08:09:57

Cheryl Straub

Operator's Signature